




GENEMEDI

Innovative Solution
for Therapeutics & Diagnostics Industry

www.genemedi.net

Antibodies And Antigens Veterinary Diagnostics Catalog

Ruminant
Companion Animals
Swine
Horse
Avian
Fish
Multiple animal



About GeneMedi

GeneMedi specializes in creating superior antibody, protein, and vector-based bioproducts, revolutionizing diagnostics and therapeutic solutions.



At GeneMedi, innovation, product integrity, and scalable solutions form the cornerstone of our mission to advance the field of diagnostics and therapeutics. Our portfolio of antibodies, proteins, and vector-derived products is built on a foundation of unparalleled expertise in the following areas:

Innovative Antigen Design and Robust Assay Development

Our strategic focus on **biomarkers** and **targets analysis** enables the creation of highly specific antigens and the development of robust assays, ensuring our products achieve superior performance in clinical and research settings.



Streamlined Molecular Discovery with Emphasis on Stability

Rapid Protein & Antibody Identification: Our proprietary platforms, TAURUS for accelerated antibody discovery and LIBRA for AI-driven protein evolution, are designed to identify and optimize molecules with optimal stability and functionality.


Cutting-Edge AAV & GCT Discoveries: The G-NEXT platform is our answer to the industry's need for innovative AAV vectors, offering improved stability, efficiency, and safety for groundbreaking gene therapy approaches.



Scalable Production and Uncompromising Quality

High-Volume Protein & Antibody Manufacturing: Our facilities are equipped to handle large-scale production up to 1000L per batch, ensuring high levels of purity and stability through stringent quality controls.

Advanced Vector Manufacturing Capabilities: With a focus on **AAV**, **Lentivirus**, and VLP production up to 200L per batch, we employ sophisticated purification techniques to guarantee vector efficacy and integrity.



Comprehensive Solutions for Diverse Application

Diagnostics: Our diagnostic solutions leverage **CLIA**, **LFA**, and **ELISA** platforms for comprehensive assay validation and clinical sample consistency, setting new standards in diagnostic accuracy.

Therapeutics: We specialize in analytical development for **ADCs**, **AAV therapies**, and **CAR-X technologies**, ensuring our therapeutic products are characterized by their specificity, potency, and safety.



GeneMedi (GM) is dedicated to delivering innovative and scalable biotechnological solutions, driving forward the fields of diagnostics and therapeutics with confidence and expertise. Contact with GM to reach your reliable industrial partner.

Table Of Contents

Part 1 Product list



Ruminants ----- 3

- Bovines
- Ovines
- Goat
- Deer



Avian ----- 7

- Avian
- Bird
- Poultry



Companion animals ----- 22

- Canine
- Feline
- Rabbit



Fish ----- 82

- Fish



Swine ----- 50

- Swine



Multiple animals ----- 82

- Cat/Feline, Dog/Canine, Rabbit, Bovines/Cattle, Equine/Horse, Ovines/Sheep, Caprine/Goat, Fish,



Horse ----- 56

- Equine/Horse

Part 2 Product list-Pseudotype

Pseudotype | Animal Host ----- 95

Part 2 Validation

Detection of SAA protein in cats, dogs and horses ----- 96

Ruminants

Cat.No	Products Name	Products Type
GMP-VT-P003-Tg002-Ag	Chlamydomphila spp. antigen	protein
GMP-VT-P003-Tg002-Ab	Anti-Chlamydomphila spp. antibody	antibody(Fc/mFc)
GMP-VT-P040-Tg001-Ag	Borrelia bovis antigen	protein
GMP-VT-P040-Tg001-Ab	Anti-Borrelia bovis antibody	antibody(Fc/mFc)
GMP-VT-P041-Tg001-Ag	BoHV-4 antigen	protein
GMP-VT-P041-Tg001-Ab	Anti-BoHV-4 antibody	antibody(Fc/mFc)
GMP-VT-P042-Tg001-Ag	Coxiella Com1 Protein antigen	protein
GMP-VT-P042-Tg001-Ab	Anti-Coxiella Com1 Protein antibody	antibody(Fc/mFc)
GMP-VT-P043-Tg001-Ag	Brucella ovis BP26 antigen	protein
GMP-VT-P043-Tg001-Ab	Anti-Brucella ovis BP26 antibody	antibody(Fc/mFc)
GMP-VT-P044-Tg001-Ag	RVA antigen	protein
GMP-VT-P044-Tg001-Ab	Anti-RVA antibody	antibody(Fc/mFc)
GMP-VT-P284-Tg001-Ag	E.coli F5 antigen	protein
GMP-VT-P284-Tg001-Ab	Anti-E.coli F5 antibody	antibody(Fc/mFc)
GMP-VT-P046-Tg001-Ag	Chlamydiaceae antigen	protein
GMP-VT-P046-Tg001-Ab	Anti-Chlamydiaceae antibody	antibody(Fc/mFc)
GMP-VT-P047-Tg001-Ag	Chlamydomphila abortus Outer Membrane Protein 91B antigen	protein
GMP-VT-P047-Tg001-Ab	Anti-Chlamydomphila abortus Outer Membrane Protein 91B antibody	antibody(Fc/mFc)
GMP-VT-P047-Tg002-Ag	Chlamydomphila abortus Processed Outer Membrane Protein 90 antigen	protein
GMP-VT-P047-Tg002-Ab	Anti-Chlamydomphila abortus Processed Outer Membrane Protein 90 antibody	antibody(Fc/mFc)
GMP-VT-P047-Tg003-Ag	Chlamydomphila abortus Major Outer Membrane Protein antigen	protein
GMP-VT-P047-Tg003-Ab	Anti-Chlamydomphila abortus Major Outer Membrane Protein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P048-Tg001-Ag	Mycobacterium avium subsp. paratuberculosis antigen	protein
GMP-VT-P048-Tg001-Ab	Anti-Mycobacterium avium subsp. paratuberculosis antibody	antibody(Fc/mFc)
GMP-VT-P049-Tg001-Ag	Mycobacterium antigen	protein
GMP-VT-P049-Tg001-Ab	Anti-Mycobacterium antibody	antibody(Fc/mFc)
GMP-VT-P050-Tg001-Ag	Mycobacterium avium ssp. paratuberculosis antigen	protein
GMP-VT-P050-Tg001-Ab	Anti-Mycobacterium avium ssp. paratuberculosis antibody	antibody(Fc/mFc)
GMP-VT-P051-Tg001-Ag	Besnoitia Besnoiti Protein Disulfide Isomerase antigen	protein
GMP-VT-P051-Tg001-Ab	Anti-Besnoitia Besnoiti Protein Disulfide Isomerase antibody	antibody(Fc/mFc)
GMP-VT-P052-Tg001-Ag	Rotavirus Viral Protein 6 antigen	protein
GMP-VT-P052-Tg001-Ab	Anti-Rotavirus Viral Protein 6 antibody	antibody(Fc/mFc)
GMP-VT-P053-Tg001-Ag	BVDV Glycoprotein 55 antigen	protein
GMP-VT-P053-Tg001-Ab	Anti-BVDV Glycoprotein 55 antibody	antibody(Fc/mFc)
GMP-VT-P053-Tg002-Ag	BVDV Glycoprotein 53 antigen	protein
GMP-VT-P053-Tg002-Ab	Anti-BVDV Glycoprotein 53 antibody	antibody(Fc/mFc)
GMP-VT-P053-Tg003-Ag	BVDV p80 antigen	protein
GMP-VT-P053-Tg003-Ab	Anti-BVDV p80 antibody	antibody(Fc/mFc)
GMP-VT-P054-Tg001-Ag	Mycoplasma pneumoniae antigen	protein
GMP-VT-P054-Tg001-Ab	Anti-Mycoplasma pneumoniae antibody	antibody(Fc/mFc)
GMP-VT-P287-Tg001-Ag	E.coli F45 antigen	protein
GMP-VT-P287-Tg001-Ab	Anti-E.coli F45 antibody	antibody(Fc/mFc)
GMP-VT-P055-Tg001-Ag	Tritrichomoas foetus antigen	protein
GMP-VT-P055-Tg001-Ab	Anti-Tritrichomoas foetus antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P056-Tg001-Ag	Pentatrichomonas hominis antigen	protein
GMP-VT-P056-Tg001-Ab	Anti-Pentatrichomonas hominis antibody	antibody(Fc/mFc)
GMP-VT-P057-Tg001-Ag	Salmonella Abortusovis antigen	protein
GMP-VT-P057-Tg001-Ab	Anti-Salmonella Abortusovis antibody	antibody(Fc/mFc)
GMP-VT-P058-Tg001-Ag	CAEV Glycoprotein 135 antigen	protein
GMP-VT-P058-Tg001-Ab	Anti-CAEV Glycoprotein 135 antibody	antibody(Fc/mFc)
GMP-VT-P058-Tg002-Ag	CAEV Protein 25 antigen	protein
GMP-VT-P058-Tg002-Ab	Anti-CAEV Protein 25 antibody	antibody(Fc/mFc)
GMP-VT-P058-Tg003-Ag	CAEV Protein 28 antigen	protein
GMP-VT-P058-Tg003-Ab	Anti-CAEV Protein 28 antibody	antibody(Fc/mFc)
GMP-VT-P058-Tg004-Ag	CAEV Capsid Protein antigen	protein
GMP-VT-P058-Tg004-Ab	Anti-CAEV Capsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P058-Tg005-Ag	CAEV Nucleocapsid Protein antigen	protein
GMP-VT-P058-Tg005-Ab	Anti-CAEV Nucleocapsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P058-Tg006-Ag	CAEV Matrix Proteins antigen	protein
GMP-VT-P058-Tg006-Ab	Anti-CAEV Matrix Proteins antibody	antibody(Fc/mFc)
GMP-VT-P059-Tg001-Ag	SBV Schmallenberg Virus Nucleoprotein antigen	protein
GMP-VT-P059-Tg001-Ab	Anti-SBV Schmallenberg Virus Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P060-Tg001-Ag	PPRV Nucleocapsid Protein antigen	protein
GMP-VT-P060-Tg001-Ab	Anti-PPRV Nucleocapsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P061-Tg001-Ag	BCV Nucleocapsid Protein antigen	protein
GMP-VT-P061-Tg001-Ab	Anti-BCV Nucleocapsid Protein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P062-Tg001-Ag	BTV Viral Protein 2 antigen	protein
GMP-VT-P062-Tg001-Ab	Anti-BTV Viral Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P062-Tg002-Ag	BTV Viral Protein 3 antigen	protein
GMP-VT-P062-Tg002-Ab	Anti-BTV Viral Protein 3 antibody	antibody(Fc/mFc)
GMP-VT-P062-Tg003-Ag	BTV Viral Protein 7 antigen	protein
GMP-VT-P062-Tg003-Ab	Anti-BTV Viral Protein 7 antibody	antibody(Fc/mFc)
GMP-VT-P062-Tg004-Ag	BTV Non-structural Protein 1 antigen	protein
GMP-VT-P062-Tg004-Ab	Anti-BTV Non-structural Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P062-Tg005-Ag	BTV Non-structural Protein 2 antigen	protein
GMP-VT-P062-Tg005-Ab	Anti-BTV Non-structural Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P062-Tg006-Ag	BTV Non-structural Protein 3 antigen	protein
GMP-VT-P062-Tg006-Ab	Anti-BTV Non-structural Protein 3 antibody	antibody(Fc/mFc)
GMP-VT-P063-Tg001-Ag	EBLV/BLV Capsid Protein p24 antigen	protein
GMP-VT-P063-Tg001-Ab	Anti-EBLV/BLV Capsid Protein p24 antibody	antibody(Fc/mFc)
GMP-VT-P064-Tg001-Ag	BoHV-1 Envelope Glycoprotein B antigen	protein
GMP-VT-P064-Tg001-Ab	Anti-BoHV-1 Envelope Glycoprotein B antibody	antibody(Fc/mFc)
GMP-VT-P064-Tg002-Ag	BoHV-1 Glycoprotein E antigen	protein
GMP-VT-P064-Tg002-Ab	Anti-BoHV-1 Glycoprotein E antibody	antibody(Fc/mFc)
GMP-VT-P065-Tg001-Ag	BoHV-2 Envelope Glycoprotein E antigen	protein
GMP-VT-P065-Tg001-Ab	Anti-BoHV-2 Envelope Glycoprotein E antibody	antibody(Fc/mFc)
GMP-VT-P066-Tg001-Ag	BRSV Nucleocapsid Protein antigen	protein
GMP-VT-P066-Tg001-Ab	Anti-BRSV Nucleocapsid Protein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P067-Tg001-Ag	IBRV Glycoprotein III antigen	protein
GMP-VT-P067-Tg001-Ab	Anti-IBRV Glycoprotein III antibody	antibody(Fc/mFc)
GMP-VT-P068-Tg001-Ag	BPI3 Nucleoprotein antigen	protein
GMP-VT-P068-Tg001-Ab	Anti-BPI3 Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P069-Tg001-Ag	BAdV-3 Hexon Protein antigen	protein
GMP-VT-P069-Tg001-Ab	Anti-BAdV-3 Hexon Protein antibody	antibody(Fc/mFc)
GMP-VT-P070-Tg001-Ag	Mycoplasma mycoides subsp. mycoides Lipoprotein LppQ antigen	protein
GMP-VT-P070-Tg001-Ab	Anti-Mycoplasma mycoides subsp. mycoides Lipoprotein LppQ antibody	antibody(Fc/mFc)
GMP-VT-P071-Tg001-Ag	tiny tapeworms of the genus Echinococcus Antigen B antigen	protein
GMP-VT-P071-Tg001-Ab	Anti-tiny tapeworms of the genus Echinococcus Antigen B antibody	antibody(Fc/mFc)
GMP-VT-P072-Tg001-Ag	BRV Viral Protein 6 antigen	protein
GMP-VT-P072-Tg001-Ab	Anti-BRV Viral Protein 6 antibody	antibody(Fc/mFc)
GMP-VT-P073-Tg001-Ag	Mycoplasma bovis Variable Surface Protein A antigen	protein
GMP-VT-P073-Tg001-Ab	Anti-Mycoplasma bovis Variable Surface Protein A antibody	antibody(Fc/mFc)
GMP-VT-P073-Tg002-Ag	Mycoplasma bovis Mycoplasma Immunogenic Lipase A antigen	protein
GMP-VT-P073-Tg002-Ab	Anti-Mycoplasma bovis Mycoplasma Immunogenic Lipase A antibody	antibody(Fc/mFc)
GMP-VT-P073-Tg003-Ag	Mycoplasma bovis Pyruvate Dehydrogenase E1 Component Beta Subunit antigen	protein
GMP-VT-P073-Tg003-Ab	Anti-Mycoplasma bovis Pyruvate Dehydrogenase E1 Component Beta Subunit antibody	antibody(Fc/mFc)
GMP-VT-P074-Tg001-Ag	Fasciola hepatica Cathepsin L1 antigen	protein
GMP-VT-P074-Tg001-Ab	Anti-Fasciola hepatica Cathepsin L1 antibody	antibody(Fc/mFc)
GMP-VT-P075-Tg001-Ag	AKAV Capsid Protein antigen	protein
GMP-VT-P075-Tg001-Ab	Anti-AKAV Capsid Protein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P076-Tg001-Ag	BDV Protein 40 antigen	protein
GMP-VT-P076-Tg001-Ab	Anti-BDV Protein 40 antibody	antibody(Fc/mFc)
GMP-VT-P077-Tg001-Ag	BRVA Viral Protein 7 antigen	protein
GMP-VT-P077-Tg001-Ab	Anti-BRVA Viral Protein 7 antibody	antibody(Fc/mFc)
GMP-VT-P078-Tg001-Ag	Campylobacter fetus Surface Layer Protein antigen	protein
GMP-VT-P078-Tg001-Ab	Anti-Campylobacter fetus Surface Layer Protein antibody	antibody(Fc/mFc)
GMP-VT-P288-Tg001-Ag	E.coli K99 Fimbrial Protein FanC antigen	protein
GMP-VT-P288-Tg001-Ab	Anti-E.coli K99 Fimbrial Protein FanC antibody	antibody(Fc/mFc)
GMP-VT-P079-Tg001-Ag	Histophilus somni Lipoprotein p31 antigen	protein
GMP-VT-P079-Tg001-Ab	Anti-Histophilus somni Lipoprotein p31 antibody	antibody(Fc/mFc)
GMP-VT-P079-Tg002-Ag	Histophilus somni Lipoprotein p40 antigen	protein
GMP-VT-P079-Tg002-Ab	Anti-Histophilus somni Lipoprotein p40 antibody	antibody(Fc/mFc)
GMP-VT-P080-Tg001-Ag	Hypoderma bovis and H. lineatum Hypodermins A antigen	protein
GMP-VT-P080-Tg001-Ab	Anti-Hypoderma bovis and H. lineatum Hypodermins A antibody	antibody(Fc/mFc)
GMP-VT-P080-Tg002-Ag	Hypoderma bovis and H. lineatum Hypodermins B antigen	protein
GMP-VT-P080-Tg002-Ab	Anti-Hypoderma bovis and H. lineatum Hypodermins B antibody	antibody(Fc/mFc)
GMP-VT-P080-Tg003-Ag	Hypoderma bovis and H. lineatum Hypodermins C antigen	protein
GMP-VT-P080-Tg003-Ab	Anti-Hypoderma bovis and H. lineatum Hypodermins C antibody	antibody(Fc/mFc)
GMP-VT-P081-Tg001-Ag	Mycoplasma agalactiae Protein 48 antigen	protein
GMP-VT-P081-Tg001-Ab	Anti-Mycoplasma agalactiae Protein 48 antibody	antibody(Fc/mFc)
GMP-VT-P082-Tg001-Ag	MVV Protein 25 antigen	protein
GMP-VT-P082-Tg001-Ab	Anti-MVV Protein 25 antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P083-Tg001-Ag	Mycoplasma capricolum subsp. capripneumoniae Whole Cell Antigen antigen	protein
GMP-VT-P083-Tg001-Ab	Anti-Mycoplasma capricolum subsp. capripneumoniae Whole Cell Antigen antibody	antibody(Fc/mFc)
GMP-VT-P084-Tg001-Ag	Mannheimia haemolytica Leukotoxin antigen	protein
GMP-VT-P084-Tg001-Ab	Anti-Mannheimia haemolytica Leukotoxin antibody	antibody(Fc/mFc)
GMP-VT-P085-Tg001-Ag	BDV Envelope Protein E2 antigen	protein
GMP-VT-P085-Tg001-Ab	Anti-BDV Envelope Protein E2 antibody	antibody(Fc/mFc)
GMP-VT-P086-Tg001-Ag	Taenia solium Ts8B2 antigen	protein
GMP-VT-P086-Tg001-Ab	Anti-Taenia solium Ts8B2 antibody	antibody(Fc/mFc)
GMP-VT-P086-Tg002-Ag	Taenia solium T24H antigen	protein
GMP-VT-P086-Tg002-Ab	Anti-Taenia solium T24H antibody	antibody(Fc/mFc)
GMP-VT-P289-Tg001-Ag	E.coli CS31A ClpG antigen	protein
GMP-VT-P289-Tg001-Ab	Anti-E.coli CS31A ClpG antibody	antibody(Fc/mFc)
GMP-VT-P087-Tg001-Ag	Leptospira sp Lipoprotein L32 antigen	protein
GMP-VT-P087-Tg001-Ab	Anti-Leptospira sp Lipoprotein L32 antibody	antibody(Fc/mFc)
GMP-VT-P088-Tg001-Ag	Mycoplasma ovipneumoniae Protein 71 antigen	protein
GMP-VT-P088-Tg001-Ab	Anti-Mycoplasma ovipneumoniae Protein 71 antibody	antibody(Fc/mFc)
GMP-VT-P089-Tg001-Ag	CPV Major Envelope Protein antigen	protein
GMP-VT-P089-Tg001-Ab	Anti-CPV Major Envelope Protein antibody	antibody(Fc/mFc)
GMP-VT-P089-Tg002-Ag	CPV Protein 32 antigen	protein
GMP-VT-P089-Tg002-Ab	Anti-CPV Protein 32 antibody	antibody(Fc/mFc)
GMP-VT-P090-Tg001-Ag	RPV Nucleoprotein antigen	protein
GMP-VT-P090-Tg001-Ab	Anti-RPV Nucleoprotein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P090-Tg002-Ag	RPV Hemagglutinin Protein antigen	protein
GMP-VT-P090-Tg002-Ab	Anti-RPV Hemagglutinin Protein antibody	antibody(Fc/mFc)
GMP-VT-P091-Tg001-Ag	BIV Transmembrane Protein antigen	protein
GMP-VT-P091-Tg001-Ab	Anti-BIV Transmembrane Protein antibody	antibody(Fc/mFc)
GMP-VT-P091-Tg002-Ag	BIV Gag Protein antigen	protein
GMP-VT-P091-Tg002-Ab	Anti-BIV Gag Protein antibody	antibody(Fc/mFc)
GMP-VT-P092-Tg001-Ag	Palyam serogroup Viral Protein 7 antigen	protein
GMP-VT-P092-Tg001-Ab	Anti-Palyam serogroup Viral Protein 7 antibody	antibody(Fc/mFc)
GMP-VT-P093-Tg001-Ag	Bovine noroviruses Capsid Proteins antigen	protein
GMP-VT-P093-Tg001-Ab	Anti-Bovine noroviruses Capsid Proteins antibody	antibody(Fc/mFc)
GMP-VT-P281-Tg001-Ag	prion Prion Protein antigen	protein
GMP-VT-P281-Tg001-Ab	Anti-prion Prion Protein antibody	antibody(Fc/mFc)
GMP-VT-P281-Tg002-Ag	prion Prion-Related Protein (Scrapie) antigen	protein
GMP-VT-P281-Tg002-Ab	Anti-prion Prion-Related Protein (Scrapie) antibody	antibody(Fc/mFc)
GMP-VT-P281-Tg003-Ag	prion Prion Protein Resistant antigen	protein
GMP-VT-P281-Tg003-Ab	Anti-prion Prion Protein Resistant antibody	antibody(Fc/mFc)
GMP-BOV-AGP-Ag	Bovine AGP antigen	protein
GMP-BOV-AGP-Ab	Anti-Bovine AGP antibody	antibody(Fc/mFc)
GMP-BOV-lactoferrin-Ag	Bovine lactoferrin antigen	protein
GMP-BOV-lactoferrin-Ab	Anti-Bovine lactoferrin antibody	antibody(Fc/mFc)
GMP-BOV-A1M-Ag	Bovine Alpha-1-microglobulin (A1M) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-A1M-Ab	Anti-Bovine Alpha-1-microglobulin (A1M) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-ACE2/hFc-Ag	Bovine Angiotensin Converting Enzyme 2 (ACE2/hFc) antigen	protein
GMP-BOV-ACE2/hFc-Ab	Anti-Bovine Angiotensin Converting Enzyme 2 (ACE2/hFc) antibody	antibody(Fc/mFc)
GMP-BOV-ACE2/his-Ag	Bovine Angiotensin Converting Enzyme 2 (His) (ACE2/his) antigen	protein
GMP-BOV-ACE2/his-Ab	Anti-Bovine Angiotensin Converting Enzyme 2 (His) (ACE2/his) antibody	antibody(Fc/mFc)
GMP-BOV-ACRV1-Ag	Bovine acrosome protein (ACRV1) antigen	protein
GMP-BOV-ACRV1-Ab	Anti-Bovine acrosome protein (ACRV1) antibody	antibody(Fc/mFc)
GMP-BOV-ACTH-Ag	Bovine adrenocorticotrophic hormone (ACTH) antigen	protein
GMP-BOV-ACTH-Ab	Anti-Bovine adrenocorticotrophic hormone (ACTH) antibody	antibody(Fc/mFc)
GMP-BOV-AD7c-NTP-Ag	Bovine AD-associated neuronal thread protein (AD7c-NTP) antigen	protein
GMP-BOV-AD7c-NTP-Ab	Anti-Bovine AD-associated neuronal thread protein (AD7c-NTP) antibody	antibody(Fc/mFc)
GMP-BOV-ADPN-Ag	Bovine Adiponectin (ADPN) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-ADPN-Ab	Anti-Bovine Adiponectin (ADPN) antibody	antibody(Fc/mFc)
GMP-BOV-AFP-Ag	Bovine Alpha-fetoprotein (AFP) antigen	protein
GMP-BOV-AFP-Ab	Anti-Bovine Alpha-fetoprotein (AFP) antibody	antibody(Fc/mFc)
GMP-BOV-AI-Ag	Bovine Angiotensin (AI) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-AI-Ab	Anti-Bovine Angiotensin (AI) antibody	antibody(Fc/mFc)
GMP-BOV-albumin-Ag	Bovine albumin protein (albumin) antigen	protein
GMP-BOV-albumin-Ab	Anti-Bovine albumin protein (albumin) antibody	antibody(Fc/mFc)
GMP-BOV-AMH-Ag	Bovine Anti-Müllerian hormone (AMH) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-AMH-Ab	Anti-Bovine Anti-Müllerian hormone (AMH) antibody	antibody(Fc/mFc)
GMP-BOV-ANP-Ag	Bovine atrial natriuretic peptide (ANP) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-ANP-Ab	Anti-Bovine atrial natriuretic peptide (ANP) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-ApoA-Ag	Bovine Apolipoprotein A (ApoA) antigen	protein
GMP-BOV-ApoA-Ab	Anti-Bovine Apolipoprotein A (ApoA) antibody	antibody(Fc/mFc)
GMP-BOV-ApoE-Ag	Bovine Apolipoprotein E (ApoE) antigen	protein
GMP-BOV-ApoE-Ab	Anti-Bovine Apolipoprotein E (ApoE) antibody	antibody(Fc/mFc)
GMP-BOV-AQP4-Ag	Bovine aquaporin-4 (AQP4) antigen	protein
GMP-BOV-AQP4-Ab	Anti-Bovine aquaporin-4 (AQP4) antibody	antibody(Fc/mFc)
GMP-BOV-AB38-Ag	Bovine beta-amyloid 38 (A β 38) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-AB38-Ab	Anti-Bovine beta-amyloid 38 (A β 38) antibody	antibody(Fc/mFc)
GMP-BOV-AB40-Ag	Bovine beta-amyloid 40 (A β 40) antigen	protein
GMP-BOV-AB40-Ab	Anti-Bovine beta-amyloid 40 (A β 40) antibody	antibody(Fc/mFc)
GMP-BOV-AB42-Ag	Bovine beta-amyloid 42 (A β 42) antigen	protein
GMP-BOV-AB42-Ab	Anti-Bovine beta-amyloid 42 (A β 42) antibody	antibody(Fc/mFc)
GMP-BOV-B2M-Ag	Bovine Beta-2 microglobulin (B2M) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-B2M-Ab	Anti-Bovine Beta-2 microglobulin (B2M) antibody	antibody(Fc/mFc)
GMP-BOV-BNP-Ag	Bovine B-type natriuretic peptide (BNP) antigen	protein
GMP-BOV-BNP-Ab	Anti-Bovine B-type natriuretic peptide (BNP) antibody	antibody(Fc/mFc)
GMP-BOV-CA-125-Ag	Bovine cancer antigen 125 (CA-125) antigen	protein
GMP-BOV-CA-125-Ab	Anti-Bovine cancer antigen 125 (CA-125) antibody	antibody(Fc/mFc)
GMP-BOV-CA15-3-Ag	Bovine Cancer antigen 15-3 (CA15-3) antigen	protein
GMP-BOV-CA15-3-Ab	Anti-Bovine Cancer antigen 15-3 (CA15-3) antibody	antibody(Fc/mFc)
GMP-BOV-CA199-Ag	Bovine Carbohydrate antigen 19-9 (CA199) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-CA199-Ab	Anti-Bovine Carbohydrate antigen 19-9 (CA199) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-CA50-Ag	Bovine Carbohydrate Antigen 50 (CA50) antigen	protein
GMP-BOV-CA50-Ab	Anti-Bovine Carbohydrate Antigen 50 (CA50) antibody	antibody(Fc/mFc)
GMP-BOV-CA-724-Ag	Bovine Cancer Antigen 72-4 (CA-724) antigen	protein
GMP-BOV-CA-724-Ab	Anti-Bovine Cancer Antigen 72-4 (CA-724) antibody	antibody(Fc/mFc)
GMP-BOV-Calcitonin-Ag	Bovine Calcitonin (Calcitonin) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-Calcitonin-Ab	Anti-Bovine Calcitonin (Calcitonin) antibody	antibody(Fc/mFc)
GMP-BOV-Calmodulin-Ag	Bovine Calmodulin antigen	protein
GMP-BOV-Calmodulin-Ab	Anti-Bovine Calmodulin antibody	antibody(Fc/mFc)
GMP-BOV-CBL-Ag	Bovine cystathionine- β -synthase enzyme (CBL) antigen	protein
GMP-BOV-CBL-Ab	Anti-Bovine cystathionine- β -synthase enzyme (CBL) antibody	antibody(Fc/mFc)
GMP-BOV-CBS-Ag	Bovine cystathionine- β -synthase enzyme (CBS) antigen	protein
GMP-BOV-CBS-Ab	Anti-Bovine cystathionine- β -synthase enzyme (CBS) antibody	antibody(Fc/mFc)
GMP-BOV-CEA-Ag	Bovine Carcinoembryonic antigen (CEA) antigen	protein
GMP-BOV-CEA-Ab	Anti-Bovine Carcinoembryonic antigen (CEA) antibody	antibody(Fc/mFc)
GMP-BOV-cIV-Ag	Bovine collagen type IV (cIV) antigen	protein
GMP-BOV-cIV-Ab	Anti-Bovine collagen type IV (cIV) antibody	antibody(Fc/mFc)
GMP-BOV-CK-MB-Ag	Bovine Creatine Kinase MB (CK-MB) antigen	protein
GMP-BOV-CK-MB-Ab	Anti-Bovine Creatine Kinase MB (CK-MB) antibody	antibody(Fc/mFc)
GMP-BOV-CNP-Ag	Bovine C-type natriuretic peptide (CNP) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-CNP-Ab	Anti-Bovine C-type natriuretic peptide (CNP) antibody	antibody(Fc/mFc)
GMP-BOV-Copeptin-Ag	Bovine Copeptin (Copeptin) antigen	protein
GMP-BOV-Copeptin-Ab	Anti-Bovine Copeptin (Copeptin) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-C-Peptide-Ag	Bovine C-peptide (C-Peptide) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-C-Peptide-Ab	Anti-Bovine C-peptide (C-Peptide) antibody	antibody(Fc/mFc)
GMP-BOV-CRP-Ag	Bovine CRP (C-Reactive Protein) antigen	protein
GMP-BOV-CRP-Ab	Anti-Bovine CRP (C-Reactive Protein) antibody	antibody(Fc/mFc)
GMP-BOV-cTnI-Ag	Bovine Cardiac troponin I (cTnI) antigen	protein
GMP-BOV-cTnI-Ab	Anti-Bovine Cardiac troponin I (cTnI) antibody	antibody(Fc/mFc)
GMP-BOV-cTnT-Ag	Bovine cardiac troponin T (cTnT) antigen	protein
GMP-BOV-cTnT-Ab	Anti-Bovine cardiac troponin T (cTnT) antibody	antibody(Fc/mFc)
GMP-BOV-CYFRA21-1-Ag	Bovine soluble fragment of cytokeratin-19 (CYFRA21-1) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-CYFRA21-1-Ab	Anti-Bovine soluble fragment of cytokeratin-19 (CYFRA21-1) antibody	antibody(Fc/mFc)
GMP-BOV-CysC-Ag	Bovine CysC (cystatin C) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-CysC-Ab	Anti-Bovine CysC (cystatin C) antibody	antibody(Fc/mFc)
GMP-BOV-D-Dimer-Ag	Bovine D dimer (D-Dimer) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-D-Dimer-Ab	Anti-Bovine D dimer (D-Dimer) antibody	antibody(Fc/mFc)
GMP-BOV-Dig-Ag	Bovine Digoxigenin (Dig) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-Dig-Ab	Anti-Bovine Digoxigenin (Dig) antibody	antibody(Fc/mFc)
GMP-BOV-FER-Ag	Bovine Ferritin (FER) antigen	protein
GMP-BOV-FER-Ab	Anti-Bovine Ferritin (FER) antibody	antibody(Fc/mFc)
GMP-BOV-fFN-Ag	Bovine Fetal fibronectin (fFN) antigen	protein
GMP-BOV-fFN-Ab	Anti-Bovine Fetal fibronectin (fFN) antibody	antibody(Fc/mFc)
GMP-BOV-FIB-Ag	Bovine Fibrinogen (FIB) antigen	protein
GMP-BOV-FIB-Ab	Anti-Bovine Fibrinogen (FIB) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-FSH-Ag	Bovine follicle stimulating hormone (FSH) antigen	protein
GMP-BOV-FSH-Ab	Anti-Bovine follicle stimulating hormone (FSH) antibody	antibody(Fc/mFc)
GMP-BOV-G100-Ag	Bovine Gastrin (G100) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-G100-Ab	Anti-Bovine Gastrin (G100) antibody	antibody(Fc/mFc)
GMP-BOV-GDF-15-Ag	Bovine Growth differentiation factor-15 (GDF-15) antigen	protein
GMP-BOV-GDF-15-Ab	Anti-Bovine Growth differentiation factor-15 (GDF-15) antibody	antibody(Fc/mFc)
GMP-BOV-GFAP-Ag	Bovine Glial fibrillary acid protein (GFAP) antigen	protein
GMP-BOV-GFAP-Ab	Anti-Bovine Glial fibrillary acid protein (GFAP) antibody	antibody(Fc/mFc)
GMP-BOV-HABP-Ag	Bovine Hyaluronic acid binding protein (HABP) antigen	protein
GMP-BOV-HABP-Ab	Anti-Bovine Hyaluronic acid binding protein (HABP) antibody	antibody(Fc/mFc)
GMP-BOV-Haptoglobin-Ag	Bovine Haptoglobin antigen	protein
GMP-BOV-Haptoglobin-Ab	Anti-Bovine Haptoglobin antibody	antibody(Fc/mFc)
GMP-BOV-HBP-Ag	Bovine Heparin-binding protein (HBP) antigen	protein
GMP-BOV-HBP-Ab	Anti-Bovine Heparin-binding protein (HBP) antibody	antibody(Fc/mFc)
GMP-BOV-HE4-Ag	Bovine epididymis protein 4 (HE4) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-HE4-Ab	Anti-Bovine epididymis protein 4 (HE4) antibody	antibody(Fc/mFc)
GMP-BOV-hEGF-Ag	Bovine epidermal growth factor (hEGF) antigen	protein
GMP-BOV-hEGF-Ab	Anti-Bovine epidermal growth factor (hEGF) antibody	antibody(Fc/mFc)
GMP-BOV-hemoglobin-Ag	Bovine hemoglobin antigen	protein
GMP-BOV-hemoglobin-Ab	Anti-Bovine hemoglobin antibody	antibody(Fc/mFc)
GMP-BOV-HER2-Ag	Bovine epidermal growth factor receptor 2 (HER2) antigen	protein
GMP-BOV-HER2-Ab	Anti-Bovine epidermal growth factor receptor 2 (HER2) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-H-FABP-Ag	Bovine Heart-type Fatty Acid-Binding Protein (H-FABP) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-H-FABP-Ab	Anti-Bovine Heart-type Fatty Acid-Binding Protein (H-FABP) antibody	antibody(Fc/mFc)
GMP-BOV-HSP90a-Ag	Bovine Heat shock proteins 90-alpha (HSP90α) antigen	protein
GMP-BOV-HSP90a-Ab	Anti-Bovine Heat shock proteins 90-alpha (HSP90α) antibody	antibody(Fc/mFc)
GMP-BOV-IFN-γ-Ag	Bovine IFN-γ (Interferon gamma) antigen	protein
GMP-BOV-IFN-γ-Ab	Anti-Bovine IFN-γ (Interferon gamma) antibody	antibody(Fc/mFc)
GMP-BOV-IgG-Ag	Bovine IgG antigen	protein
GMP-BOV-IgG-Ab	Anti-Bovine IgG antibody	antibody(Fc/mFc)
GMP-BOV-IL-6-Ag	Bovine Interleukin 6 (IL-6) antigen	protein
GMP-BOV-IL-6-Ab	Anti-Bovine Interleukin 6 (IL-6) antibody	antibody(Fc/mFc)
GMP-BOV-Ins-Ag	Bovine Insulin (Ins) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-Ins-Ab	Anti-Bovine Insulin (Ins) antibody	antibody(Fc/mFc)
GMP-BOV-LN-Ag	Bovine laminins (LN) antigen	protein
GMP-BOV-LN-Ab	Anti-Bovine laminins (LN) antibody	antibody(Fc/mFc)
GMP-BOV-MCP-1-Ag	Bovine Monocyte chemoattractant protein-1 (MCP-1) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-MCP-1-Ab	Anti-Bovine Monocyte chemoattractant protein-1 (MCP-1) antibody	antibody(Fc/mFc)
GMP-BOV-MMP3-Ag	Bovine matrix metalloproteinase-3 (MMP3) antigen	protein
GMP-BOV-MMP3-Ab	Anti-Bovine matrix metalloproteinase-3 (MMP3) antibody	antibody(Fc/mFc)
GMP-BOV-MOG-Ag	Bovine myelin oligodendrocyte glycoprotein (MOG) antigen	protein
GMP-BOV-MOG-Ab	Anti-Bovine myelin oligodendrocyte glycoprotein (MOG) antibody	antibody(Fc/mFc)
GMP-BOV-MPO-Ag	Bovine Myeloperoxidase (MPO) antigen	protein
GMP-BOV-MPO-Ab	Anti-Bovine Myeloperoxidase (MPO) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-Myo-Ag	Bovine Myoglobin (Myo) antigen	protein
GMP-BOV-Myo-Ab	Anti-Bovine Myoglobin (Myo) antibody	antibody(Fc/mFc)
GMP-BOV-neurogranin-Ag	Bovine neurogranin protein (neurogranin) antigen	protein
GMP-BOV-neurogranin-Ab	Anti-Bovine neurogranin protein (neurogranin) antibody	antibody(Fc/mFc)
GMP-BOV-NFL-Ag	Bovine neurofilament light chain (NFL) antigen	protein
GMP-BOV-NFL-Ab	Anti-Bovine neurofilament light chain (NFL) antibody	antibody(Fc/mFc)
GMP-BOV-NGAL-Ag	Bovine Neutrophil gelatinase-associated lipocalin (NGAL) antigen	protein
GMP-BOV-NGAL-Ab	Anti-Bovine Neutrophil gelatinase-associated lipocalin (NGAL) antibody	antibody(Fc/mFc)
GMP-BOV-N-MID-OC-Ag	Bovine N-terminal midfragment of Osteocalcin (N-MID OC) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-N-MID-OC-Ab	Anti-Bovine N-terminal midfragment of Osteocalcin (N-MID OC) antibody	antibody(Fc/mFc)
GMP-BOV-NSE-Ag	Bovine Neuron specific enolase (NSE) antigen	protein
GMP-BOV-NSE-Ab	Anti-Bovine Neuron specific enolase (NSE) antibody	antibody(Fc/mFc)
GMP-BOV-NT-proBNP-Ag	Bovine NT-proBNP (N-terminal prohormone of brain natriuretic peptide) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-NT-proBNP-Ab	Anti-Bovine NT-proBNP (N-terminal prohormone of brain natriuretic peptide) antibody	antibody(Fc/mFc)
GMP-BOV-PA-Ag	Bovine Prealbumin (PA) antigen	protein
GMP-BOV-PA-Ab	Anti-Bovine Prealbumin (PA) antibody	antibody(Fc/mFc)
GMP-BOV-Pancreatic-elastase-Ag	Bovine PE-1 (Pancreatic elastase) antigen	protein
GMP-BOV-Pancreatic-elastase-Ab	Anti-Bovine PE-1 (Pancreatic elastase) antibody	antibody(Fc/mFc)
GMP-BOV-pancreatic-lipase-Ag	Bovine pancreatic lipase antigen	protein
GMP-BOV-pancreatic-lipase-Ab	Anti-Bovine pancreatic lipase antibody	antibody(Fc/mFc)
GMP-BOV-PCT-Ag	Bovine procalcitonin (PCT) antigen	protein
GMP-BOV-PCT-Ab	Anti-Bovine procalcitonin (PCT) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-PG-I-Ag	Bovine pepsinogen I (PG I) antigen	protein
GMP-BOV-PG-I-Ab	Anti-Bovine pepsinogen I (PG I) antibody	antibody(Fc/mFc)
GMP-BOV-PG-II-Ag	Bovine Pepsinogen II (PG II) antigen	protein
GMP-BOV-PG-II-Ab	Anti-Bovine Pepsinogen II (PG II) antibody	antibody(Fc/mFc)
GMP-BOV-PIIINP-Ag	Bovine Procollagen III N-terminal peptide (PIIINP) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-PIIINP-Ab	Anti-Bovine Procollagen III N-terminal peptide (PIIINP) antibody	antibody(Fc/mFc)
GMP-BOV-PINP-Ag	Bovine amino-terminal propeptide (PINP) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-PINP-Ab	Anti-Bovine amino-terminal propeptide (PINP) antibody	antibody(Fc/mFc)
GMP-BOV-PIVKA-II-Ag	Bovine Protein induced by vitamin K absence-II (PIVKA-II) antigen	protein
GMP-BOV-PIVKA-II-Ab	Anti-Bovine Protein induced by vitamin K absence-II (PIVKA-II) antibody	antibody(Fc/mFc)
GMP-BOV-PLGF-1-Ag	Bovine placental growth factor-1 (PLGF-1) antigen	protein
GMP-BOV-PLGF-1-Ab	Anti-Bovine placental growth factor-1 (PLGF-1) antibody	antibody(Fc/mFc)
GMP-BOV-Pre-Spec-Pro-B-Ag	Bovine Pregnancy Specific Protein B antigen	protein
GMP-BOV-Pre-Spec-Pro-B-Ab	Anti-Bovine Pregnancy Specific Protein B antibody	antibody(Fc/mFc)
GMP-BOV-Preg-ass-gly-1-Ag	Bovine Pregnancy-associated glycoprotein 1 antigen	protein
GMP-BOV-Preg-ass-gly-1-Ab	Anti-Bovine Pregnancy-associated glycoprotein 1 antibody	antibody(Fc/mFc)
GMP-BOV-proBNP-Ag	Bovine pro B-type natriuretic peptide (proBNP) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-proBNP-Ab	Anti-Bovine pro B-type natriuretic peptide (proBNP) antibody	antibody(Fc/mFc)
GMP-BOV-Prog-Ag	Bovine Prog (Progesterone) antigen	protein
GMP-BOV-Prog-Ab	Anti-Bovine Prog (Progesterone) antibody	antibody(Fc/mFc)
GMP-BOV-ProGRP-Ag	Bovine Gastrin-releasing peptide (ProGRP) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-ProGRP-Ab	Anti-Bovine Gastrin-releasing peptide (ProGRP) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-Proinsulin-Ag	Bovine Proinsulin (Proinsulin) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-Proinsulin-Ab	Anti-Bovine Proinsulin (Proinsulin) antibody	antibody(Fc/mFc)
GMP-BOV-PSA-Ag	Bovine Prostate-specific antigen (PSA) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-PSA-Ab	Anti-Bovine Prostate-specific antigen (PSA) antibody	antibody(Fc/mFc)
GMP-BOV-p-tau181-Ag	Bovine Phospho-tau181 (p-tau181) antigen	protein
GMP-BOV-p-tau181-Ab	Anti-Bovine Phospho-tau181 (p-tau181) antibody	antibody(Fc/mFc)
GMP-BOV-p-tau202-205-Ag	Bovine Phospho-tau202 and 205 (p-tau202 and 205) antigen	protein
GMP-BOV-p-tau202-205-Ab	Anti-Bovine Phospho-tau202 and 205 (p-tau202 and 205) antibody	antibody(Fc/mFc)
GMP-BOV-p-tau212-214-Ag	Bovine Phospho-tau212 and 214 (p-tau212 and 214) antigen	protein
GMP-BOV-p-tau212-214-Ab	Anti-Bovine Phospho-tau212 and 214 (p-tau212 and 214) antibody	antibody(Fc/mFc)
GMP-BOV-p-tau217-Ag	Bovine Phospho-tau217 (p-tau217) antigen	protein
GMP-BOV-p-tau217-Ab	Anti-Bovine Phospho-tau217 (p-tau217) antibody	antibody(Fc/mFc)
GMP-BOV-p-tau231-Ag	Bovine Phospho-tau231 (p-tau231) antigen	protein
GMP-BOV-p-tau231-Ab	Anti-Bovine Phospho-tau231 (p-tau231) antibody	antibody(Fc/mFc)
GMP-BOV-p-tau413-Ag	Bovine Phospho-tau413 (p-tau413) antigen	protein
GMP-BOV-p-tau413-Ab	Anti-Bovine Phospho-tau413 (p-tau413) antibody	antibody(Fc/mFc)
GMP-BOV-p-tau422-Ag	Bovine Phospho-tau422 (p-tau422) antigen	protein
GMP-BOV-p-tau422-Ab	Anti-Bovine Phospho-tau422 (p-tau422) antibody	antibody(Fc/mFc)
GMP-BOV-PTH-Ag	Bovine parathyroid hormone (PTH) antigen	protein
GMP-BOV-PTH-Ab	Anti-Bovine parathyroid hormone (PTH) antibody	antibody(Fc/mFc)
GMP-BOV-relaxin-Ag	Bovine RLN (Relaxin) antigen	protein
GMP-BOV-relaxin-Ab	Anti-Bovine RLN (Relaxin) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-S100-A8-Ag	Bovine Calprotectin (S100-A8) antigen	protein
GMP-BOV-S100-A8-Ab	Anti-Bovine Calprotectin (S100-A8) antibody	antibody(Fc/mFc)
GMP-BOV-S100-A9-Ag	Bovine Calprotectin heavy chain (S100-A9) antigen	protein
GMP-BOV-S100-A9-Ab	Anti-Bovine Calprotectin heavy chain (S100-A9) antibody	antibody(Fc/mFc)
GMP-BOV-S100B-Ag	Bovine S100 proteins (S100B) antigen	protein
GMP-BOV-S100B-Ab	Anti-Bovine S100 proteins (S100B) antibody	antibody(Fc/mFc)
GMP-BOV-SAHH-Ag	Bovine S-adenosylhomocysteine (AdoHcy) hydrolase (SAHH) antigen	protein
GMP-BOV-SAHH-Ab	Anti-Bovine S-adenosylhomocysteine (AdoHcy) hydrolase (SAHH) antibody	antibody(Fc/mFc)
GMP-BOV-SCC-Ag-Ag	Bovine squamous cell carcinoma antigen (SCC Ag) antigen	protein
GMP-BOV-SCC-Ag-Ab	Anti-Bovine squamous cell carcinoma antigen (SCC Ag) antibody	antibody(Fc/mFc)
GMP-BOV-Serum-amyloid-A-Ag	Bovine SAA (Serum amyloid A) antigen	protein
GMP-BOV-Serum-amyloid-A-Ab	Anti-Bovine SAA (Serum amyloid A) antibody	antibody(Fc/mFc)
GMP-BOV-sFlt-1-Ag	Bovine soluble vascular endothelial growth factor receptor-2 (sFlt-1) antigen	protein
GMP-BOV-sFlt-1-Ab	Anti-Bovine soluble vascular endothelial growth factor receptor-2 (sFlt-1) antibody	antibody(Fc/mFc)
GMP-BOV-SPA-Ag	Bovine Staphylococcus protein A (SPA) antigen	protein
GMP-BOV-SPA-Ab	Anti-Bovine Staphylococcus protein A (SPA) antibody	antibody(Fc/mFc)
GMP-BOV-Spirochaeta-Ag	Bovine Spirochaeta antigen	protein
GMP-BOV-Spirochaeta-Ab	Anti-Bovine Spirochaeta antibody	antibody(Fc/mFc)
GMP-BOV-sST2-Ag	Bovine soluble ST2 (sST2) antigen	protein
GMP-BOV-sST2-Ab	Anti-Bovine soluble ST2 (sST2) antibody	antibody(Fc/mFc)
GMP-BOV-sTREM2-Ag	Bovine soluble fragment of triggering receptor expressed on myeloid cells 2 (sTREM2) antigen	protein
GMP-BOV-sTREM2-Ab	Anti-Bovine soluble fragment of triggering receptor expressed on myeloid cells 2 (sTREM2) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-Tau-Ag	Bovine tau proteins (Tau) antigen	protein
GMP-BOV-Tau-Ab	Anti-Bovine tau proteins (Tau) antibody	antibody(Fc/mFc)
GMP-BOV-TG-Ag	Bovine Thyroglobulin (TG) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-TG-Ab	Anti-Bovine Thyroglobulin (TG) antibody	antibody(Fc/mFc)
GMP-BOV-TIMP-1-Ag	Bovine Matrix Metalloproteinase Inhibitor-1 (TIMP-1) antigen	protein
GMP-BOV-TIMP-1-Ab	Anti-Bovine Matrix Metalloproteinase Inhibitor-1 (TIMP-1) antibody	antibody(Fc/mFc)
GMP-BOV-TK1-Ag	Bovine Thymidine Kinase 1 (TK1) antigen	protein
GMP-BOV-TK1-Ab	Anti-Bovine Thymidine Kinase 1 (TK1) antibody	antibody(Fc/mFc)
GMP-BOV-TNFa-Ag	Bovine tumor necrosis factor alpha (TNF α) antigen	protein
GMP-BOV-TNFa-Ab	Anti-Bovine tumor necrosis factor alpha (TNF α) antibody	antibody(Fc/mFc)
GMP-BOV-TPO-Ag	Bovine Thyroid peroxidase (TPO) antigen	protein
GMP-BOV-TPO-Ab	Anti-Bovine Thyroid peroxidase (TPO) antibody	antibody(Fc/mFc)
GMP-BOV-TrxR-Ag	Bovine thioredoxin reductase (TrxR) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-TrxR-Ab	Anti-Bovine thioredoxin reductase (TrxR) antibody	antibody(Fc/mFc)
GMP-BOV-TSH-Ag	Bovine Thyroid-stimulating hormone (TSH) antigen	protein
GMP-BOV-TSH-Ab	Anti-Bovine Thyroid-stimulating hormone (TSH) antibody	antibody(Fc/mFc)
GMP-BOV-VDBP-Ag	Bovine vitamin D binding protein (VDBP) antigen	protein
GMP-BOV-VDBP-Ab	Anti-Bovine vitamin D binding protein (VDBP) antibody	antibody(Fc/mFc)
GMP-BOV-VEGF-Ag	Bovine Vascular Endothelial Growth Factor (VEGF) antigen	protein
GMP-BOV-VEGF-Ab	Anti-Bovine Vascular Endothelial Growth Factor (VEGF) antibody	antibody(Fc/mFc)
GMP-BOV-VLP-1-Ag	Bovine Visinin-like protein 1 (VLP-1) antigen	protein
GMP-BOV-VLP-1-Ab	Anti-Bovine Visinin-like protein 1 (VLP-1) antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-BOV-YKL-40-Ag	Bovine Chitinase-3-like protein 1 (YKL-40) antigen	protein
GMP-BOV-YKL-40-Ab	Anti-Bovine Chitinase-3-like protein 1 (YKL-40) antibody	antibody(Fc/mFc)
GMP-BOV-a-defensin-1-Ag	Bovine alpha defensin 1 (α -defensin 1) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-a-defensin-1-Ab	Anti-Bovine alpha defensin 1 (α -defensin 1) antibody	antibody(Fc/mFc)
GMP-BOV-a-defensin-5-Ag	Bovine alpha defensin 5 (α -defensin 5) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-a-defensin-5-Ab	Anti-Bovine alpha defensin 5 (α -defensin 5) antibody	antibody(Fc/mFc)
GMP-BOV-a-synuclein-Ag	Bovine Alpha-synuclein (α -synuclein) antigen	protein
GMP-BOV-a-synuclein-Ab	Anti-Bovine Alpha-synuclein (α -synuclein) antibody	antibody(Fc/mFc)
GMP-BOV-B-CTx-Ag	Bovine β -isomerized C-terminal telopeptides (β -CTx) antigen	peptide-(naked/BSA/OVA)
GMP-BOV-B-CTx-Ab	Anti-Bovine β -isomerized C-terminal telopeptides (β -CTx) antibody	antibody(Fc/mFc)

Companion Animals

Cat.No	Products Name	Products Type
GMP-VT-P001-Tg001-Ag	CHV antigen	protein
GMP-VT-P001-Tg001-Ab	Anti-CHV antibody	antibody(Fc/mFc)
GMP-VT-P002-Tg001-Ag	Anaplasma phagocytophilum Major Surface Protein 44 antigen	protein
GMP-VT-P002-Tg001-Ab	Anti-Anaplasma phagocytophilum Major Surface Protein 44 antibody	antibody(Fc/mFc)
GMP-VT-P002-Tg002-Ag	Anaplasma phagocytophilum Major Surface Protein 44 antigen	protein
GMP-VT-P002-Tg002-Ab	Anti-Anaplasma phagocytophilum Major Surface Protein 44 antibody	antibody(Fc/mFc)
GMP-VT-P002-Tg003-Ag	Anaplasma phagocytophilum Major Surface Protein 4 antigen	protein
GMP-VT-P002-Tg003-Ab	Anti-Anaplasma phagocytophilum Major Surface Protein 4 antibody	antibody(Fc/mFc)
GMP-VT-P003-Tg001-Ag	Chlamydomphila spp. Outer Membrane Protein A antigen	protein

Cat.No	Products Name	Products Type
GMP-VT-P003-Tg001-Ab	Anti-Chlamydomphila spp. Outer Membrane Protein A antibody	antibody(Fc/mFc)
GMP-VT-P004-Tg001-Ag	RHDV Viral Protein 60 antigen	protein
GMP-VT-P004-Tg001-Ab	Anti-RHDV Viral Protein 60 antibody	antibody(Fc/mFc)
GMP-VT-P005-Tg001-Ag	CCV antigen	protein
GMP-VT-P005-Tg001-Ab	Anti-CCV antibody	antibody(Fc/mFc)
GMP-VT-P006-Tg001-Ag	Angiostrongylus vasorum antigen	protein
GMP-VT-P006-Tg001-Ab	Anti-Angiostrongylus vasorum antibody	antibody(Fc/mFc)
GMP-VT-P007-Tg001-Ag	CAV-2 antigen	protein
GMP-VT-P007-Tg001-Ab	Anti-CAV-2 antibody	antibody(Fc/mFc)
GMP-VT-P008-Tg001-Ag	CAV-1 antigen	protein
GMP-VT-P008-Tg001-Ab	Anti-CAV-1 antibody	antibody(Fc/mFc)
GMP-VT-P009-Tg001-Ag	Cowbone Ridge virus antigen	protein
GMP-VT-P009-Tg001-Ab	Anti-Cowbone Ridge virus antibody	antibody(Fc/mFc)
GMP-VT-P010-Tg001-Ag	Malassezia antigen	protein
GMP-VT-P010-Tg001-Ab	Anti-Malassezia antibody	antibody(Fc/mFc)
GMP-VT-P011-Tg001-Ag	CPV antigen	protein
GMP-VT-P011-Tg001-Ab	Anti-CPV antibody	antibody(Fc/mFc)
GMP-VT-P012-Tg001-Ag	CDV Nucleocapsid Protein antigen	protein
GMP-VT-P012-Tg001-Ab	Anti-CDV Nucleocapsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P013-Tg001-Ag	CAV Hexon Protein antigen	protein
GMP-VT-P013-Tg001-Ab	Anti-CAV Hexon Protein antibody	antibody(Fc/mFc)
GMP-VT-P014-Tg001-Ag	FPV Viral Protein 2 antigen	protein

Cat.No	Products Name	Products Type
GMP-VT-P014-Tg001-Ab	Anti-FPV Viral Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P015-Tg001-Ag	FeLV P27 antigen	protein
GMP-VT-P015-Tg001-Ab	Anti-FeLV P27 antibody	antibody(Fc/mFc)
GMP-VT-P016-Tg001-Ag	CPIV-2 Nucleocapsid Protein antigen	protein
GMP-VT-P016-Tg001-Ab	Anti-CPIV-2 Nucleocapsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P017-Tg001-Ag	Ehrlichia Major Antigenic Protein 2 antigen	protein
GMP-VT-P017-Tg001-Ab	Anti-Ehrlichia Major Antigenic Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P018-Tg001-Ag	Rabies virus Glycoprotein antigen	protein
GMP-VT-P018-Tg001-Ab	Anti-Rabies virus Glycoprotein antibody	antibody(Fc/mFc)
GMP-VT-P019-Tg001-Ag	FHV Glycoprotein B antigen	protein
GMP-VT-P019-Tg001-Ab	Anti-FHV Glycoprotein B antibody	antibody(Fc/mFc)
GMP-VT-P020-Tg001-Ag	FPV Viral Protein 2 antigen	protein
GMP-VT-P020-Tg001-Ab	Anti-FPV Viral Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P021-Tg001-Ag	FCV Viral Protein 1 antigen	protein
GMP-VT-P021-Tg001-Ab	Anti-FCV Viral Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P022-Tg001-Ag	FIV Capsid Protein p24 antigen	protein
GMP-VT-P022-Tg001-Ab	Anti-FIV Capsid Protein p24 antibody	antibody(Fc/mFc)
GMP-VT-P023-Tg001-Ag	Leishmania infantum Paraflagellar Rod Proteins antigen	protein
GMP-VT-P023-Tg001-Ab	Anti-Leishmania infantum Paraflagellar Rod Proteins antibody	antibody(Fc/mFc)
GMP-VT-P023-Tg002-Ag	Leishmania infantum Paraflagellar Rod Proteins antigen	protein
GMP-VT-P023-Tg002-Ab	Anti-Leishmania infantum Paraflagellar Rod Proteins antibody	antibody(Fc/mFc)
GMP-VT-P023-Tg003-Ag	Leishmania infantum Kinesin-Related Protein K39 antigen	protein

Cat.No	Products Name	Products Type
GMP-VT-P023-Tg003-Ab	Anti-Leishmania infantum Kinesin-Related Protein K39 antibody	antibody(Fc/mFc)
GMP-VT-P024-Tg001-Ag	Dirofilaria immitis Dirofilaria Glycoprotein K antigen	protein
GMP-VT-P024-Tg001-Ab	Anti-Dirofilaria immitis Dirofilaria Glycoprotein K antibody	antibody(Fc/mFc)
GMP-VT-P024-Tg002-Ag	Dirofilaria immitis Dirofilaria Immunogenic Protein 33 antigen	protein
GMP-VT-P024-Tg002-Ab	Anti-Dirofilaria immitis Dirofilaria Immunogenic Protein 33 antibody	antibody(Fc/mFc)
GMP-VT-P024-Tg003-Ag	Dirofilaria immitis Dirofilaria Immunogenic Protein 33 antigen	protein
GMP-VT-P024-Tg003-Ab	Anti-Dirofilaria immitis Dirofilaria Immunogenic Protein 33 antibody	antibody(Fc/mFc)
GMP-VT-P025-Tg001-Ag	FCoV/FIPV Spike Protein antigen	protein
GMP-VT-P025-Tg001-Ab	Anti-FCoV/FIPV Spike Protein antibody	antibody(Fc/mFc)
GMP-VT-P026-Tg002-Ag	Leptospira Lipoprotein L21 antigen	protein
GMP-VT-P026-Tg002-Ab	Anti-Leptospira Lipoprotein L21 antibody	antibody(Fc/mFc)
GMP-VT-P027-Tg001-Ag	Brucella canis SlyX Protein antigen	protein
GMP-VT-P027-Tg001-Ab	Anti-Brucella canis SlyX Protein antibody	antibody(Fc/mFc)
GMP-VT-P028-Tg001-Ag	Canine leptospira Outer Membrane Protein L1 antigen	protein
GMP-VT-P028-Tg001-Ab	Anti-Canine leptospira Outer Membrane Protein L1 antibody	antibody(Fc/mFc)
GMP-VT-P029-Tg001-Ag	Ehrlichia canis P30 antigen	protein
GMP-VT-P029-Tg001-Ab	Anti-Ehrlichia canis P30 antibody	antibody(Fc/mFc)
GMP-VT-P030-Tg001-Ag	CPV/CPV-2 Viral Protein 2 antigen	protein
GMP-VT-P030-Tg001-Ab	Anti-CPV/CPV-2 Viral Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P031-Tg001-Ag	Anaplasma platys Outer Membrane Protein 1X antigen	protein
GMP-VT-P031-Tg001-Ab	Anti-Anaplasma platys Outer Membrane Protein 1X antibody	antibody(Fc/mFc)
GMP-VT-P032-Tg001-Ag	Bartonella henselae Protein 26 antigen	protein

Cat.No	Products Name	Products Type
GMP-VT-P032-Tg001-Ab	Anti-Bartonella henselae Protein 26 antibody	antibody(Fc/mFc)
GMP-VT-P033-Tg001-Ag	Chlamydomphila felis CF0218 antigen	protein
GMP-VT-P033-Tg001-Ab	Anti-Chlamydomphila felis CF0218 antibody	antibody(Fc/mFc)
GMP-VT-P034-Tg001-Ag	Cryptococcus neoformans Heat Shock Protein 70 antigen	protein
GMP-VT-P034-Tg001-Ab	Anti-Cryptococcus neoformans Heat Shock Protein 70 antibody	antibody(Fc/mFc)
GMP-VT-P035-Tg001-Ag	Francisella tularensis Outer Membrane Protein A antigen	protein
GMP-VT-P035-Tg001-Ab	Anti-Francisella tularensis Outer Membrane Protein A antibody	antibody(Fc/mFc)
GMP-VT-P036-Tg001-Ag	Rickettsia conorii Outer Membrane Proteins A and B antigen	protein
GMP-VT-P036-Tg001-Ab	Anti-Rickettsia conorii Outer Membrane Proteins A and B antibody	antibody(Fc/mFc)
GMP-VT-P037-Tg001-Ag	Rickettsia massiliae Translocated Icm/Dot-Like Substrates antigen	protein
GMP-VT-P037-Tg001-Ab	Anti-Rickettsia massiliae Translocated Icm/Dot-Like Substrates antibody	antibody(Fc/mFc)
GMP-VT-P038-Tg001-Ag	Rickettsia raoultii antigen	protein
GMP-VT-P038-Tg001-Ab	Anti-Rickettsia raoultii antibody	antibody(Fc/mFc)
GMP-VT-P039-Tg001-Ag	Rickettsia rickettsii Outer Membrane Protein antigen	protein
GMP-VT-P039-Tg001-Ab	Anti-Rickettsia rickettsii Outer Membrane Protein antibody	antibody(Fc/mFc)
GMP-FEL-AGP-Ag	Feline AGP antigen	protein
GMP-FEL-AGP-Ab	Anti-Feline AGP antibody	antibody(Fc/mFc)
GMP-CAN-AGP-Ag	Canine AGP antigen	protein
GMP-CAN-AGP-Ab	Anti-Canine AGP antibody	antibody(Fc/mFc)
GMP-FEL-lactoferrin-Ag	Feline lactoferrin antigen	protein
GMP-FEL-lactoferrin-Ab	Anti-Feline lactoferrin antibody	antibody(Fc/mFc)
GMP-CAN-lactoferrin-Ag	Canine lactoferrin antigen	protein

Cat.No	Products Name	Products Type
GMP-CAN-lactoferrin-Ab	Anti-Canine lactoferrin antibody	antibody(Fc/mFc)
GMP-FEL-A1M-Ag	Feline Alpha-1-microglobulin (A1M) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-A1M-Ab	Anti-Feline Alpha-1-microglobulin (A1M) antibody	antibody(Fc/mFc)
GMP-CAN-A1M-Ag	Canine Alpha-1-microglobulin (A1M) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-A1M-Ab	Anti-Canine Alpha-1-microglobulin (A1M) antibody	antibody(Fc/mFc)
GMP-FEL-ACE2/hFc-Ag	Feline Angiotensin Converting Enzyme 2 (ACE2/hFc) antigen	protein
GMP-FEL-ACE2/hFc-Ab	Anti-Feline Angiotensin Converting Enzyme 2 (ACE2/hFc) antibody	antibody(Fc/mFc)
GMP-CAN-ACE2/hFc-Ag	Canine Angiotensin Converting Enzyme 2 (ACE2/hFc) antigen	protein
GMP-CAN-ACE2/hFc-Ab	Anti-Canine Angiotensin Converting Enzyme 2 (ACE2/hFc) antibody	antibody(Fc/mFc)
GMP-FEL-ACE2/his-Ag	Feline Angiotensin Converting Enzyme 2 (His) (ACE2/his) antigen	protein
GMP-FEL-ACE2/his-Ab	Anti-Feline Angiotensin Converting Enzyme 2 (His) (ACE2/his) antibody	antibody(Fc/mFc)
GMP-CAN-ACE2/his-Ag	Canine Angiotensin Converting Enzyme 2 (His) (ACE2/his) antigen	protein
GMP-CAN-ACE2/his-Ab	Anti-Canine Angiotensin Converting Enzyme 2 (His) (ACE2/his) antibody	antibody(Fc/mFc)
GMP-FEL-ACRV1-Ag	Feline acrosome protein (ACRV1) antigen	protein
GMP-FEL-ACRV1-Ab	Anti-Feline acrosome protein (ACRV1) antibody	antibody(Fc/mFc)
GMP-CAN-ACRV1-Ag	Canine acrosome protein (ACRV1) antigen	protein
GMP-CAN-ACRV1-Ab	Anti-Canine acrosome protein (ACRV1) antibody	antibody(Fc/mFc)
GMP-FEL-ACTH-Ag	Feline adrenocorticotrophic hormone (ACTH) antigen	protein
GMP-FEL-ACTH-Ab	Anti-Feline adrenocorticotrophic hormone (ACTH) antibody	antibody(Fc/mFc)
GMP-CAN-ACTH-Ag	Canine adrenocorticotrophic hormone (ACTH) antigen	protein
GMP-CAN-ACTH-Ab	Anti-Canine adrenocorticotrophic hormone (ACTH) antibody	antibody(Fc/mFc)
GMP-FEL-AD7c-NTP-Ag	Feline AD-associated neuronal thread protein (AD7c-NTP) antigen	protein

Cat.No	Products Name	Products Type
GMP-FEL-AD7c-NTP-Ab	Anti-Feline AD-associated neuronal thread protein (AD7c-NTP) antibody	antibody(Fc/mFc)
GMP-CAN-AD7c-NTP-Ag	Canine AD-associated neuronal thread protein (AD7c-NTP) antigen	protein
GMP-CAN-AD7c-NTP-Ab	Anti-Canine AD-associated neuronal thread protein (AD7c-NTP) antibody	antibody(Fc/mFc)
GMP-FEL-ADPN-Ag	Feline Adiponectin (ADPN) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-ADPN-Ab	Anti-Feline Adiponectin (ADPN) antibody	antibody(Fc/mFc)
GMP-CAN-ADPN-Ag	Canine Adiponectin (ADPN) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-ADPN-Ab	Anti-Canine Adiponectin (ADPN) antibody	antibody(Fc/mFc)
GMP-FEL-AFP-Ag	Feline Alpha-fetoprotein (AFP) antigen	protein
GMP-FEL-AFP-Ab	Anti-Feline Alpha-fetoprotein (AFP) antibody	antibody(Fc/mFc)
GMP-CAN-AFP-Ag	Canine Alpha-fetoprotein (AFP) antigen	protein
GMP-CAN-AFP-Ab	Anti-Canine Alpha-fetoprotein (AFP) antibody	antibody(Fc/mFc)
GMP-FEL-AI-Ag	Feline Angiotensin (AI) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-AI-Ab	Anti-Feline Angiotensin (AI) antibody	antibody(Fc/mFc)
GMP-CAN-AI-Ag	Canine Angiotensin (AI) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-AI-Ab	Anti-Canine Angiotensin (AI) antibody	antibody(Fc/mFc)
GMP-FEL-albumin-Ag	Feline albumin protein (albumin) antigen	protein
GMP-FEL-albumin-Ab	Anti-Feline albumin protein (albumin) antibody	antibody(Fc/mFc)
GMP-CAN-albumin-Ag	Canine albumin protein (albumin) antigen	protein
GMP-CAN-albumin-Ab	Anti-Canine albumin protein (albumin) antibody	antibody(Fc/mFc)
GMP-FEL-AMH-Ag	Feline Anti-Müllerian hormone (AMH) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-AMH-Ab	Anti-Feline Anti-Müllerian hormone (AMH) antibody	antibody(Fc/mFc)
GMP-CAN-AMH-Ag	Canine Anti-Müllerian hormone (AMH) antigen	peptide-(naked/BSA/OVA)

Cat.No	Products Name	Products Type
GMP-CAN-AMH-Ab	Anti-Canine Anti-Mu"llerian hormone (AMH) antibody	antibody(Fc/mFc)
GMP-FEL-ANP-Ag	Feline atrial natriuretic peptide (ANP) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-ANP-Ab	Anti-Feline atrial natriuretic peptide (ANP) antibody	antibody(Fc/mFc)
GMP-CAN-ANP-Ag	Canine atrial natriuretic peptide (ANP) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-ANP-Ab	Anti-Canine atrial natriuretic peptide (ANP) antibody	antibody(Fc/mFc)
GMP-FEL-ApoA-Ag	Feline Apolipoprotein A (ApoA) antigen	protein
GMP-FEL-ApoA-Ab	Anti-Feline Apolipoprotein A (ApoA) antibody	antibody(Fc/mFc)
GMP-CAN-ApoA-Ag	Canine Apolipoprotein A (ApoA) antigen	protein
GMP-CAN-ApoA-Ab	Anti-Canine Apolipoprotein A (ApoA) antibody	antibody(Fc/mFc)
GMP-FEL-ApoE-Ag	Feline Apolipoprotein E (ApoE) antigen	protein
GMP-FEL-ApoE-Ab	Anti-Feline Apolipoprotein E (ApoE) antibody	antibody(Fc/mFc)
GMP-CAN-ApoE-Ag	Canine Apolipoprotein E (ApoE) antigen	protein
GMP-CAN-ApoE-Ab	Anti-Canine Apolipoprotein E (ApoE) antibody	antibody(Fc/mFc)
GMP-FEL-AQP4-Ag	Feline aquaporin-4 (AQP4) antigen	protein
GMP-FEL-AQP4-Ab	Anti-Feline aquaporin-4 (AQP4) antibody	antibody(Fc/mFc)
GMP-CAN-AQP4-Ag	Canine aquaporin-4 (AQP4) antigen	protein
GMP-CAN-AQP4-Ab	Anti-Canine aquaporin-4 (AQP4) antibody	antibody(Fc/mFc)
GMP-FEL-AB38-Ag	Feline beta-amyloid 38 (A β 38) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-AB38-Ab	Anti-Feline beta-amyloid 38 (A β 38) antibody	antibody(Fc/mFc)
GMP-CAN-AB38-Ag	Canine beta-amyloid 38 (A β 38) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-AB38-Ab	Anti-Canine beta-amyloid 38 (A β 38) antibody	antibody(Fc/mFc)
GMP-FEL-AB40-Ag	Feline beta-amyloid 40 (A β 40) antigen	protein

Cat.No	Products Name	Products Type
GMP-FEL-AB40-Ab	Anti-Feline beta-amyloid 40 (A β 40) antibody	antibody(Fc/mFc)
GMP-CAN-AB40-Ag	Canine beta-amyloid 40 (A β 40) antigen	protein
GMP-CAN-AB40-Ab	Anti-Canine beta-amyloid 40 (A β 40) antibody	antibody(Fc/mFc)
GMP-FEL-AB42-Ag	Feline beta-amyloid 42 (A β 42) antigen	protein
GMP-FEL-AB42-Ab	Anti-Feline beta-amyloid 42 (A β 42) antibody	antibody(Fc/mFc)
GMP-CAN-AB42-Ag	Canine beta-amyloid 42 (A β 42) antigen	protein
GMP-CAN-AB42-Ab	Anti-Canine beta-amyloid 42 (A β 42) antibody	antibody(Fc/mFc)
GMP-FEL-B2M-Ag	Feline Beta-2 microglobulin (B2M) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-B2M-Ab	Anti-Feline Beta-2 microglobulin (B2M) antibody	antibody(Fc/mFc)
GMP-CAN-B2M-Ag	Canine Beta-2 microglobulin (B2M) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-B2M-Ab	Anti-Canine Beta-2 microglobulin (B2M) antibody	antibody(Fc/mFc)
GMP-FEL-BNP-Ag	Feline B-type natriuretic peptide (BNP) antigen	protein
GMP-FEL-BNP-Ab	Anti-Feline B-type natriuretic peptide (BNP) antibody	antibody(Fc/mFc)
GMP-CAN-BNP-Ag	Canine B-type natriuretic peptide (BNP) antigen	protein
GMP-CAN-BNP-Ab	Anti-Canine B-type natriuretic peptide (BNP) antibody	antibody(Fc/mFc)
GMP-FEL-CA-125-Ag	Feline cancer antigen 125 (CA-125) antigen	protein
GMP-FEL-CA-125-Ab	Anti-Feline cancer antigen 125 (CA-125) antibody	antibody(Fc/mFc)
GMP-CAN-CA-125-Ag	Canine cancer antigen 125 (CA-125) antigen	protein
GMP-CAN-CA-125-Ab	Anti-Canine cancer antigen 125 (CA-125) antibody	antibody(Fc/mFc)
GMP-FEL-CA15-3-Ag	Feline Cancer antigen 15-3 (CA15-3) antigen	protein
GMP-FEL-CA15-3-Ab	Anti-Feline Cancer antigen 15-3 (CA15-3) antibody	antibody(Fc/mFc)
GMP-CAN-CA15-3-Ag	Canine Cancer antigen 15-3 (CA15-3) antigen	protein

Cat.No	Products Name	Products Type
GMP-CAN-CA15-3-Ab	Anti-Canine Cancer antigen 15-3 (CA15-3) antibody	antibody(Fc/mFc)
GMP-FEL-CA199-Ag	Feline Carbohydrate antigen 19-9 (CA199) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-CA199-Ab	Anti-Feline Carbohydrate antigen 19-9 (CA199) antibody	antibody(Fc/mFc)
GMP-CAN-CA199-Ag	Canine Carbohydrate antigen 19-9 (CA199) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-CA199-Ab	Anti-Canine Carbohydrate antigen 19-9 (CA199) antibody	antibody(Fc/mFc)
GMP-FEL-CA50-Ag	Feline Carbohydrate Antigen 50 (CA50) antigen	protein
GMP-FEL-CA50-Ab	Anti-Feline Carbohydrate Antigen 50 (CA50) antibody	antibody(Fc/mFc)
GMP-CAN-CA50-Ag	Canine Carbohydrate Antigen 50 (CA50) antigen	protein
GMP-CAN-CA50-Ab	Anti-Canine Carbohydrate Antigen 50 (CA50) antibody	antibody(Fc/mFc)
GMP-FEL-CA-724-Ag	Feline Cancer Antigen 72-4 (CA-724) antigen	protein
GMP-FEL-CA-724-Ab	Anti-Feline Cancer Antigen 72-4 (CA-724) antibody	antibody(Fc/mFc)
GMP-CAN-CA-724-Ag	Canine Cancer Antigen 72-4 (CA-724) antigen	protein
GMP-CAN-CA-724-Ab	Anti-Canine Cancer Antigen 72-4 (CA-724) antibody	antibody(Fc/mFc)
GMP-FEL-Calcitonin-Ag	Feline Calcitonin (Calcitonin) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-Calcitonin-Ab	Anti-Feline Calcitonin (Calcitonin) antibody	antibody(Fc/mFc)
GMP-CAN-Calcitonin-Ag	Canine Calcitonin (Calcitonin) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-Calcitonin-Ab	Anti-Canine Calcitonin (Calcitonin) antibody	antibody(Fc/mFc)
GMP-FEL-Calmodulin-Ag	Feline Calmodulin antigen	protein
GMP-FEL-Calmodulin-Ab	Anti-Feline Calmodulin antibody	antibody(Fc/mFc)
GMP-CAN-Calmodulin-Ag	Canine Calmodulin antigen	protein
GMP-CAN-Calmodulin-Ab	Anti-Canine Calmodulin antibody	antibody(Fc/mFc)
GMP-FEL-CBL-Ag	Feline cystathionine-β-synthase enzyme (CBL) antigen	protein

Cat.No	Products Name	Products Type
GMP-FEL-CBL-Ab	Anti-Feline cystathionine-β-synthase enzyme (CBL) antibody	antibody(Fc/mFc)
GMP-CAN-CBL-Ag	Canine cystathionine-β-synthase enzyme (CBL) antigen	protein
GMP-CAN-CBL-Ab	Anti-Canine cystathionine-β-synthase enzyme (CBL) antibody	antibody(Fc/mFc)
GMP-FEL-CBS-Ag	Feline cystathionine-β-synthase enzyme (CBS) antigen	protein
GMP-FEL-CBS-Ab	Anti-Feline cystathionine-β-synthase enzyme (CBS) antibody	antibody(Fc/mFc)
GMP-CAN-CBS-Ag	Canine cystathionine-β-synthase enzyme (CBS) antigen	protein
GMP-CAN-CBS-Ab	Anti-Canine cystathionine-β-synthase enzyme (CBS) antibody	antibody(Fc/mFc)
GMP-FEL-CEA-Ag	Feline Carcinoembryonic antigen (CEA) antigen	protein
GMP-FEL-CEA-Ab	Anti-Feline Carcinoembryonic antigen (CEA) antibody	antibody(Fc/mFc)
GMP-CAN-CEA-Ag	Canine Carcinoembryonic antigen (CEA) antigen	protein
GMP-CAN-CEA-Ab	Anti-Canine Carcinoembryonic antigen (CEA) antibody	antibody(Fc/mFc)
GMP-FEL-clV-Ag	Feline collagen type IV (clV) antigen	protein
GMP-FEL-clV-Ab	Anti-Feline collagen type IV (clV) antibody	antibody(Fc/mFc)
GMP-CAN-clV-Ag	Canine collagen type IV (clV) antigen	protein
GMP-CAN-clV-Ab	Anti-Canine collagen type IV (clV) antibody	antibody(Fc/mFc)
GMP-FEL-CK-MB-Ag	Feline Creatine Kinase MB (CK-MB) antigen	protein
GMP-FEL-CK-MB-Ab	Anti-Feline Creatine Kinase MB (CK-MB) antibody	antibody(Fc/mFc)
GMP-CAN-CK-MB-Ag	Canine Creatine Kinase MB (CK-MB) antigen	protein
GMP-CAN-CK-MB-Ab	Anti-Canine Creatine Kinase MB (CK-MB) antibody	antibody(Fc/mFc)
GMP-FEL-CNP-Ag	Feline C-type natriuretic peptide (CNP) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-CNP-Ab	Anti-Feline C-type natriuretic peptide (CNP) antibody	antibody(Fc/mFc)
GMP-CAN-CNP-Ag	Canine C-type natriuretic peptide (CNP) antigen	peptide-(naked/BSA/OVA)

Cat.No	Products Name	Products Type
GMP-CAN-CNP-Ab	Anti-Canine C-type natriuretic peptide (CNP) antibody	antibody(Fc/mFc)
GMP-FEL-Copeptin-Ag	Feline Copeptin (Copeptin) antigen	protein
GMP-FEL-Copeptin-Ab	Anti-Feline Copeptin (Copeptin) antibody	antibody(Fc/mFc)
GMP-CAN-Copeptin-Ag	Canine Copeptin (Copeptin) antigen	protein
GMP-CAN-Copeptin-Ab	Anti-Canine Copeptin (Copeptin) antibody	antibody(Fc/mFc)
GMP-FEL-C-Peptide-Ag	Feline C-peptide (C-Peptide) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-C-Peptide-Ab	Anti-Feline C-peptide (C-Peptide) antibody	antibody(Fc/mFc)
GMP-CAN-C-Peptide-Ag	Canine C-peptide (C-Peptide) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-C-Peptide-Ab	Anti-Canine C-peptide (C-Peptide) antibody	antibody(Fc/mFc)
GMP-CAN-CRP-Ag	Canine CRP (C-Reactive Protein) antigen	protein
GMP-CAN-CRP-Ab	Anti-Canine CRP (C-Reactive Protein) antibody	antibody(Fc/mFc)
GMP-FEL-CRP-Ag	Feline CRP (C-Reactive Protein) antigen	protein
GMP-FEL-CRP-Ab	Anti-Feline CRP (C-Reactive Protein) antibody	antibody(Fc/mFc)
GMP-FEL-cTnI-Ag	Feline Cardiac troponin I (cTnI) antigen	protein
GMP-FEL-cTnI-Ab	Anti-Feline Cardiac troponin I (cTnI) antibody	antibody(Fc/mFc)
GMP-CAN-cTnI-Ag	Canine Cardiac troponin I (cTnI) antigen	protein
GMP-CAN-cTnI-Ab	Anti-Canine Cardiac troponin I (cTnI) antibody	antibody(Fc/mFc)
GMP-FEL-cTnT-Ag	Feline cardiac troponin T (cTnT) antigen	protein
GMP-FEL-cTnT-Ab	Anti-Feline cardiac troponin T (cTnT) antibody	antibody(Fc/mFc)
GMP-CAN-cTnT-Ag	Canine cardiac troponin T (cTnT) antigen	protein
GMP-CAN-cTnT-Ab	Anti-Canine cardiac troponin T (cTnT) antibody	antibody(Fc/mFc)
GMP-FEL-CYFRA21-1-Ag	Feline soluble fragment of cytokeratin-19 (CYFRA21-1) antigen	peptide-(naked/BSA/OVA)

Cat.No	Products Name	Products Type
GMP-FEL-CYFRA21-1-Ab	Anti-Feline soluble fragment of cytokeratin-19 (CYFRA21-1) antibody	antibody(Fc/mFc)
GMP-CAN-CYFRA21-1-Ag	Canine soluble fragment of cytokeratin-19 (CYFRA21-1) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-CYFRA21-1-Ab	Anti-Canine soluble fragment of cytokeratin-19 (CYFRA21-1) antibody	antibody(Fc/mFc)
GMP-CAN-CysC-Ag	Canine CysC (cystatin C) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-CysC-Ab	Anti-Canine CysC (cystatin C) antibody	antibody(Fc/mFc)
GMP-FEL-CysC-Ag	Feline CysC (cystatin C) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-CysC-Ab	Anti-Feline CysC (cystatin C) antibody	antibody(Fc/mFc)
GMP-FEL-D-Dimer-Ag	Feline D dimer (D-Dimer) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-D-Dimer-Ab	Anti-Feline D dimer (D-Dimer) antibody	antibody(Fc/mFc)
GMP-CAN-D-Dimer-Ag	Canine D dimer (D-Dimer) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-D-Dimer-Ab	Anti-Canine D dimer (D-Dimer) antibody	antibody(Fc/mFc)
GMP-FEL-Dig-Ag	Feline Digoxigenin (Dig) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-Dig-Ab	Anti-Feline Digoxigenin (Dig) antibody	antibody(Fc/mFc)
GMP-CAN-Dig-Ag	Canine Digoxigenin (Dig) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-Dig-Ab	Anti-Canine Digoxigenin (Dig) antibody	antibody(Fc/mFc)
GMP-FEL-IgG-Ag	Feline IgG antigen	protein
GMP-FEL-IgG-Ab	Anti-Feline IgG antibody	antibody(Fc/mFc)
GMP-FEL-FER-Ag	Feline Ferritin (FER) antigen	protein
GMP-FEL-FER-Ab	Anti-Feline Ferritin (FER) antibody	antibody(Fc/mFc)
GMP-CAN-FER-Ag	Canine Ferritin (FER) antigen	protein
GMP-CAN-FER-Ab	Anti-Canine Ferritin (FER) antibody	antibody(Fc/mFc)
GMP-FEL-fFN-Ag	Feline Fetal fibronectin (fFN) antigen	protein

Cat.No	Products Name	Products Type
GMP-FEL-fFN-Ab	Anti-Feline Fetal fibronectin (fFN) antibody	antibody(Fc/mFc)
GMP-CAN-fFN-Ag	Canine Fetal fibronectin (fFN) antigen	protein
GMP-CAN-fFN-Ab	Anti-Canine Fetal fibronectin (fFN) antibody	antibody(Fc/mFc)
GMP-FEL-FIB-Ag	Feline Fibrinogen (FIB) antigen	protein
GMP-FEL-FIB-Ab	Anti-Feline Fibrinogen (FIB) antibody	antibody(Fc/mFc)
GMP-CAN-FIB-Ag	Canine Fibrinogen (FIB) antigen	protein
GMP-CAN-FIB-Ab	Anti-Canine Fibrinogen (FIB) antibody	antibody(Fc/mFc)
GMP-FEL-FSH-Ag	Feline follicle stimulating hormone (FSH) antigen	protein
GMP-FEL-FSH-Ab	Anti-Feline follicle stimulating hormone (FSH) antibody	antibody(Fc/mFc)
GMP-CAN-FSH-Ag	Canine follicle stimulating hormone (FSH) antigen	protein
GMP-CAN-FSH-Ab	Anti-Canine follicle stimulating hormone (FSH) antibody	antibody(Fc/mFc)
GMP-FEL-G100-Ag	Feline Gastrin (G100) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-G100-Ab	Anti-Feline Gastrin (G100) antibody	antibody(Fc/mFc)
GMP-CAN-G100-Ag	Canine Gastrin (G100) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-G100-Ab	Anti-Canine Gastrin (G100) antibody	antibody(Fc/mFc)
GMP-FEL-GDF-15-Ag	Feline Growth differentiation factor-15 (GDF-15) antigen	protein
GMP-FEL-GDF-15-Ab	Anti-Feline Growth differentiation factor-15 (GDF-15) antibody	antibody(Fc/mFc)
GMP-CAN-GDF-15-Ag	Canine Growth differentiation factor-15 (GDF-15) antigen	protein
GMP-CAN-GDF-15-Ab	Anti-Canine Growth differentiation factor-15 (GDF-15) antibody	antibody(Fc/mFc)
GMP-FEL-GFAP-Ag	Feline Glial fibrillary acid protein (GFAP) antigen	protein
GMP-FEL-GFAP-Ab	Anti-Feline Glial fibrillary acid protein (GFAP) antibody	antibody(Fc/mFc)
GMP-CAN-GFAP-Ag	Canine Glial fibrillary acid protein (GFAP) antigen	protein

Cat.No	Products Name	Products Type
GMP-CAN-GFAP-Ab	Anti-Canine Glial fibrillary acid protein (GFAP) antibody	antibody(Fc/mFc)
GMP-FEL-HABP-Ag	Feline Hyaluronic acid binding protein (HABP) antigen	protein
GMP-FEL-HABP-Ab	Anti-Feline Hyaluronic acid binding protein (HABP) antibody	antibody(Fc/mFc)
GMP-CAN-HABP-Ag	Canine Hyaluronic acid binding protein (HABP) antigen	protein
GMP-CAN-HABP-Ab	Anti-Canine Hyaluronic acid binding protein (HABP) antibody	antibody(Fc/mFc)
GMP-FEL-Haptoglobin-Ag	Feline Haptoglobin antigen	protein
GMP-FEL-Haptoglobin-Ab	Anti-Feline Haptoglobin antibody	antibody(Fc/mFc)
GMP-FEL-HBP-Ag	Feline Heparin-binding protein (HBP) antigen	protein
GMP-FEL-HBP-Ab	Anti-Feline Heparin-binding protein (HBP) antibody	antibody(Fc/mFc)
GMP-CAN-HBP-Ag	Canine Heparin-binding protein (HBP) antigen	protein
GMP-CAN-HBP-Ab	Anti-Canine Heparin-binding protein (HBP) antibody	antibody(Fc/mFc)
GMP-FEL-HE4-Ag	Feline epididymis protein 4 (HE4) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-HE4-Ab	Anti-Feline epididymis protein 4 (HE4) antibody	antibody(Fc/mFc)
GMP-CAN-HE4-Ag	Canine epididymis protein 4 (HE4) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-HE4-Ab	Anti-Canine epididymis protein 4 (HE4) antibody	antibody(Fc/mFc)
GMP-FEL-hEGF-Ag	Feline epidermal growth factor (hEGF) antigen	protein
GMP-FEL-hEGF-Ab	Anti-Feline epidermal growth factor (hEGF) antibody	antibody(Fc/mFc)
GMP-CAN-hEGF-Ag	Canine epidermal growth factor (hEGF) antigen	protein
GMP-CAN-hEGF-Ab	Anti-Canine epidermal growth factor (hEGF) antibody	antibody(Fc/mFc)
GMP-CAN-Hemoglobin-Ag	Canine hemoglobin antigen	protein
GMP-CAN-Hemoglobin-Ab	Anti-Canine hemoglobin antibody	antibody(Fc/mFc)
GMP-FEL-hemoglobin-Ag	Feline hemoglobin antigen	protein

Cat.No	Products Name	Products Type
GMP-FEL-hemoglobin-Ab	Anti-Feline hemoglobin antibody	antibody(Fc/mFc)
GMP-FEL-HER2-Ag	Feline epidermal growth factor receptor 2 (HER2) antigen	protein
GMP-FEL-HER2-Ab	Anti-Feline epidermal growth factor receptor 2 (HER2) antibody	antibody(Fc/mFc)
GMP-CAN-HER2-Ag	Canine epidermal growth factor receptor 2 (HER2) antigen	protein
GMP-CAN-HER2-Ab	Anti-Canine epidermal growth factor receptor 2 (HER2) antibody	antibody(Fc/mFc)
GMP-FEL-H-FABP-Ag	Feline Heart-type Fatty Acid-Binding Protein (H-FABP) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-H-FABP-Ab	Anti-Feline Heart-type Fatty Acid-Binding Protein (H-FABP) antibody	antibody(Fc/mFc)
GMP-CAN-H-FABP-Ag	Canine Heart-type Fatty Acid-Binding Protein (H-FABP) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-H-FABP-Ab	Anti-Canine Heart-type Fatty Acid-Binding Protein (H-FABP) antibody	antibody(Fc/mFc)
GMP-FEL-HSP90a-Ag	Feline Heat shock proteins 90-alpha (HSP90α) antigen	protein
GMP-FEL-HSP90a-Ab	Anti-Feline Heat shock proteins 90-alpha (HSP90α) antibody	antibody(Fc/mFc)
GMP-CAN-HSP90a-Ag	Canine Heat shock proteins 90-alpha (HSP90α) antigen	protein
GMP-CAN-HSP90a-Ab	Anti-Canine Heat shock proteins 90-alpha (HSP90α) antibody	antibody(Fc/mFc)
GMP-FEL-IFN-γ-Ag	Feline IFN-γ (Interferon gamma) antigen	protein
GMP-FEL-IFN-γ-Ab	Anti-Feline IFN-γ (Interferon gamma) antibody	antibody(Fc/mFc)
GMP-CAN-IFN-γ-Ag	Canine IFN-γ (Interferon gamma) antigen	protein
GMP-CAN-IFN-γ-Ab	Anti-Canine IFN-γ (Interferon gamma) antibody	antibody(Fc/mFc)
GMP-CAN-IgG-Ag	Canine IgG antigen	protein
GMP-CAN-IgG-Ab	Anti-Canine IgG antibody	antibody(Fc/mFc)
GMP-FEL-IL-6-Ag	Feline Interleukin 6 (IL-6) antigen	protein
GMP-FEL-IL-6-Ab	Anti-Feline Interleukin 6 (IL-6) antibody	antibody(Fc/mFc)
GMP-CAN-IL-6-Ag	Canine Interleukin 6 (IL-6) antigen	protein

Cat.No	Products Name	Products Type
GMP-CAN-IL-6-Ab	Anti-Canine Interleukin 6 (IL-6) antibody	antibody(Fc/mFc)
GMP-FEL-Ins-Ag	Feline Insulin (Ins) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-Ins-Ab	Anti-Feline Insulin (Ins) antibody	antibody(Fc/mFc)
GMP-CAN-Ins-Ag	Canine Insulin (Ins) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-Ins-Ab	Anti-Canine Insulin (Ins) antibody	antibody(Fc/mFc)
GMP-FEL-LN-Ag	Feline laminins (LN) antigen	protein
GMP-FEL-LN-Ab	Anti-Feline laminins (LN) antibody	antibody(Fc/mFc)
GMP-CAN-LN-Ag	Canine laminins (LN) antigen	protein
GMP-CAN-LN-Ab	Anti-Canine laminins (LN) antibody	antibody(Fc/mFc)
GMP-FEL-MCP-1-Ag	Feline Monocyte chemoattractant protein-1 (MCP-1) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-MCP-1-Ab	Anti-Feline Monocyte chemoattractant protein-1 (MCP-1) antibody	antibody(Fc/mFc)
GMP-CAN-MCP-1-Ag	Canine Monocyte chemoattractant protein-1 (MCP-1) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-MCP-1-Ab	Anti-Canine Monocyte chemoattractant protein-1 (MCP-1) antibody	antibody(Fc/mFc)
GMP-FEL-MMP3-Ag	Feline matrix metalloproteinase-3 (MMP3) antigen	protein
GMP-FEL-MMP3-Ab	Anti-Feline matrix metalloproteinase-3 (MMP3) antibody	antibody(Fc/mFc)
GMP-CAN-MMP3-Ag	Canine matrix metalloproteinase-3 (MMP3) antigen	protein
GMP-CAN-MMP3-Ab	Anti-Canine matrix metalloproteinase-3 (MMP3) antibody	antibody(Fc/mFc)
GMP-FEL-MOG-Ag	Feline myelin oligodendrocyte glycoprotein (MOG) antigen	protein
GMP-FEL-MOG-Ab	Anti-Feline myelin oligodendrocyte glycoprotein (MOG) antibody	antibody(Fc/mFc)
GMP-CAN-MOG-Ag	Canine myelin oligodendrocyte glycoprotein (MOG) antigen	protein
GMP-CAN-MOG-Ab	Anti-Canine myelin oligodendrocyte glycoprotein (MOG) antibody	antibody(Fc/mFc)
GMP-FEL-MPO-Ag	Feline Myeloperoxidase (MPO) antigen	protein

Cat.No	Products Name	Products Type
GMP-FEL-MPO-Ab	Anti-Feline Myeloperoxidase (MPO) antibody	antibody(Fc/mFc)
GMP-CAN-MPO-Ag	Canine Myeloperoxidase (MPO) antigen	protein
GMP-CAN-MPO-Ab	Anti-Canine Myeloperoxidase (MPO) antibody	antibody(Fc/mFc)
GMP-FEL-Myo-Ag	Feline Myoglobin (Myo) antigen	protein
GMP-FEL-Myo-Ab	Anti-Feline Myoglobin (Myo) antibody	antibody(Fc/mFc)
GMP-CAN-Myo-Ag	Canine Myoglobin (Myo) antigen	protein
GMP-CAN-Myo-Ab	Anti-Canine Myoglobin (Myo) antibody	antibody(Fc/mFc)
GMP-FEL-neurogranin-Ag	Feline neurogranin protein (neurogranin) antigen	protein
GMP-FEL-neurogranin-Ab	Anti-Feline neurogranin protein (neurogranin) antibody	antibody(Fc/mFc)
GMP-CAN-neurogranin-Ag	Canine neurogranin protein (neurogranin) antigen	protein
GMP-CAN-neurogranin-Ab	Anti-Canine neurogranin protein (neurogranin) antibody	antibody(Fc/mFc)
GMP-FEL-NFL-Ag	Feline neurofilament light chain (NFL) antigen	protein
GMP-FEL-NFL-Ab	Anti-Feline neurofilament light chain (NFL) antibody	antibody(Fc/mFc)
GMP-CAN-NFL-Ag	Canine neurofilament light chain (NFL) antigen	protein
GMP-CAN-NFL-Ab	Anti-Canine neurofilament light chain (NFL) antibody	antibody(Fc/mFc)
GMP-FEL-NGAL-Ag	Feline Neutrophil gelatinase-associated lipocalin (NGAL) antigen	protein
GMP-FEL-NGAL-Ab	Anti-Feline Neutrophil gelatinase-associated lipocalin (NGAL) antibody	antibody(Fc/mFc)
GMP-CAN-NGAL-Ag	Canine Neutrophil gelatinase-associated lipocalin (NGAL) antigen	protein
GMP-CAN-NGAL-Ab	Anti-Canine Neutrophil gelatinase-associated lipocalin (NGAL) antibody	antibody(Fc/mFc)
GMP-FEL-N-MID-OC-Ag	Feline N-terminal midfragment of Osteocalcin (N-MID OC) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-N-MID-OC-Ab	Anti-Feline N-terminal midfragment of Osteocalcin (N-MID OC) antibody	antibody(Fc/mFc)
GMP-CAN-N-MID-OC-Ag	Canine N-terminal midfragment of Osteocalcin (N-MID OC) antigen	peptide-(naked/BSA/OVA)

Cat.No	Products Name	Products Type
GMP-CAN-N-MID-OC-Ab	Anti-Canine N-terminal midfragment of Osteocalcin (N-MID OC) antibody	antibody(Fc/mFc)
GMP-FEL-NSE-Ag	Feline Neuron specific enolase (NSE) antigen	protein
GMP-FEL-NSE-Ab	Anti-Feline Neuron specific enolase (NSE) antibody	antibody(Fc/mFc)
GMP-CAN-NSE-Ag	Canine Neuron specific enolase (NSE) antigen	protein
GMP-CAN-NSE-Ab	Anti-Canine Neuron specific enolase (NSE) antibody	antibody(Fc/mFc)
GMP-FEL-NT-proBNP-Ag	Feline NT-proBNP (N-terminal prohormone of brain natriuretic peptide) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-NT-proBNP-Ab	Anti-Feline NT-proBNP (N-terminal prohormone of brain natriuretic peptide) antibody	antibody(Fc/mFc)
GMP-CAN-NT-proBNP-Ag	Canine NT-proBNP (N-terminal prohormone of brain natriuretic peptide) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-NT-proBNP-Ab	Anti-Canine NT-proBNP (N-terminal prohormone of brain natriuretic peptide) antibody	antibody(Fc/mFc)
GMP-FEL-PA-Ag	Feline Prealbumin (PA) antigen	protein
GMP-FEL-PA-Ab	Anti-Feline Prealbumin (PA) antibody	antibody(Fc/mFc)
GMP-CAN-PA-Ag	Canine Prealbumin (PA) antigen	protein
GMP-CAN-PA-Ab	Anti-Canine Prealbumin (PA) antibody	antibody(Fc/mFc)
GMP-FEL-Pancreatic-elastase-Ag	Feline PE-1 (Pancreatic elastase) antigen	protein
GMP-FEL-Pancreatic-elastase-Ab	Anti-Feline PE-1 (Pancreatic elastase) antibody	antibody(Fc/mFc)
GMP-CAN-Pancreatic-elastase-Ag	Canine PE-1 (Pancreatic elastase) antigen	protein
GMP-CAN-Pancreatic-elastase-Ab	Anti-Canine PE-1 (Pancreatic elastase) antibody	antibody(Fc/mFc)
GMP-CAN-Pancreatic-lipase-Ag	Canine pancreatic lipase antigen	protein
GMP-CAN-Pancreatic-lipase-Ab	Anti-Canine pancreatic lipase antibody	antibody(Fc/mFc)
GMP-FEL-Pancreatic-lipase-Ag	Feline pancreatic lipase antigen	protein
GMP-FEL-Pancreatic-lipase-Ab	Anti-Feline pancreatic lipase antibody	antibody(Fc/mFc)
GMP-FEL-PCT-Ag	Feline procalcitonin (PCT) antigen	protein

Cat.No	Products Name	Products Type
GMP-FEL-PCT-Ab	Anti-Feline procalcitonin (PCT) antibody	antibody(Fc/mFc)
GMP-CAN-PCT-Ag	Canine procalcitonin (PCT) antigen	protein
GMP-CAN-PCT-Ab	Anti-Canine procalcitonin (PCT) antibody	antibody(Fc/mFc)
GMP-FEL-PG-I-Ag	Feline pepsinogen I (PG I) antigen	protein
GMP-FEL-PG-I-Ab	Anti-Feline pepsinogen I (PG I) antibody	antibody(Fc/mFc)
GMP-CAN-PG-I-Ag	Canine pepsinogen I (PG I) antigen	protein
GMP-CAN-PG-I-Ab	Anti-Canine pepsinogen I (PG I) antibody	antibody(Fc/mFc)
GMP-FEL-PG-II-Ag	Feline Pepsinogen II (PG II) antigen	protein
GMP-FEL-PG-II-Ab	Anti-Feline Pepsinogen II (PG II) antibody	antibody(Fc/mFc)
GMP-CAN-PG-II-Ag	Canine Pepsinogen II (PG II) antigen	protein
GMP-CAN-PG-II-Ab	Anti-Canine Pepsinogen II (PG II) antibody	antibody(Fc/mFc)
GMP-FEL-PIIINP-Ag	Feline Procollagen III N-terminal peptide (PIIINP) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-PIIINP-Ab	Anti-Feline Procollagen III N-terminal peptide (PIIINP) antibody	antibody(Fc/mFc)
GMP-CAN-PIIINP-Ag	Canine Procollagen III N-terminal peptide (PIIINP) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-PIIINP-Ab	Anti-Canine Procollagen III N-terminal peptide (PIIINP) antibody	antibody(Fc/mFc)
GMP-FEL-PINP-Ag	Feline amino-terminal propeptide (PINP) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-PINP-Ab	Anti-Feline amino-terminal propeptide (PINP) antibody	antibody(Fc/mFc)
GMP-CAN-PINP-Ag	Canine amino-terminal propeptide (PINP) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-PINP-Ab	Anti-Canine amino-terminal propeptide (PINP) antibody	antibody(Fc/mFc)
GMP-FEL-PIVKA- II -Ag	Feline Protein induced by vitamin K absence-II (PIVKA- II) antigen	protein
GMP-FEL-PIVKA- II -Ab	Anti-Feline Protein induced by vitamin K absence-II (PIVKA- II) antibody	antibody(Fc/mFc)
GMP-CAN-PIVKA- II -Ag	Canine Protein induced by vitamin K absence-II (PIVKA- II) antigen	protein

Cat.No	Products Name	Products Type
GMP-CAN-PIVKA-II-Ab	Anti-Canine Protein induced by vitamin K absence-II (PIVKA-II) antibody	antibody(Fc/mFc)
GMP-FEL-PLGF-1-Ag	Feline placental growth factor-1 (PLGF-1) antigen	protein
GMP-FEL-PLGF-1-Ab	Anti-Feline placental growth factor-1 (PLGF-1) antibody	antibody(Fc/mFc)
GMP-CAN-PLGF-1-Ag	Canine placental growth factor-1 (PLGF-1) antigen	protein
GMP-CAN-PLGF-1-Ab	Anti-Canine placental growth factor-1 (PLGF-1) antibody	antibody(Fc/mFc)
GMP-FEL-Pre-Spec-Pro-B-Ag	Feline Pregnancy Specific Protein B antigen	protein
GMP-FEL-Pre-Spec-Pro-B-Ab	Anti-Feline Pregnancy Specific Protein B antibody	antibody(Fc/mFc)
GMP-CAN-Pre-Spec-Pro-B-Ag	Canine Pregnancy Specific Protein B antigen	protein
GMP-CAN-Pre-Spec-Pro-B-Ab	Anti-Canine Pregnancy Specific Protein B antibody	antibody(Fc/mFc)
GMP-FEL-Preg-ass-gly-1-Ag	Feline Pregnancy-associated glycoprotein 1 antigen	protein
GMP-FEL-Preg-ass-gly-1-Ab	Anti-Feline Pregnancy-associated glycoprotein 1 antibody	antibody(Fc/mFc)
GMP-CAN-Preg-ass-gly-1-Ag	Canine Pregnancy-associated glycoprotein 1 antigen	protein
GMP-CAN-Preg-ass-gly-1-Ab	Anti-Canine Pregnancy-associated glycoprotein 1 antibody	antibody(Fc/mFc)
GMP-FEL-proBNP-Ag	Feline pro B-type natriuretic peptide (proBNP) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-proBNP-Ab	Anti-Feline pro B-type natriuretic peptide (proBNP) antibody	antibody(Fc/mFc)
GMP-CAN-proBNP-Ag	Canine pro B-type natriuretic peptide (proBNP) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-proBNP-Ab	Anti-Canine pro B-type natriuretic peptide (proBNP) antibody	antibody(Fc/mFc)
GMP-FEL-Prog-Ag	Feline Prog (Progesterone) antigen	protein
GMP-FEL-Prog-Ab	Anti-Feline Prog (Progesterone) antibody	antibody(Fc/mFc)
GMP-CAN-Prog-Ag	Canine Prog (Progesterone) antigen	protein
GMP-CAN-Prog-Ab	Anti-Canine Prog (Progesterone) antibody	antibody(Fc/mFc)
GMP-FEL-ProGRP-Ag	Feline Gastrin-releasing peptide (ProGRP) antigen	peptide-(naked/BSA/OVA)

Cat.No	Products Name	Products Type
GMP-FEL-ProGRP-Ab	Anti-Feline Gastrin-releasing peptide (ProGRP) antibody	antibody(Fc/mFc)
GMP-CAN-ProGRP-Ag	Canine Gastrin-releasing peptide (ProGRP) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-ProGRP-Ab	Anti-Canine Gastrin-releasing peptide (ProGRP) antibody	antibody(Fc/mFc)
GMP-FEL-Proinsulin-Ag	Feline Proinsulin (Proinsulin) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-Proinsulin-Ab	Anti-Feline Proinsulin (Proinsulin) antibody	antibody(Fc/mFc)
GMP-CAN-Proinsulin-Ag	Canine Proinsulin (Proinsulin) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-Proinsulin-Ab	Anti-Canine Proinsulin (Proinsulin) antibody	antibody(Fc/mFc)
GMP-FEL-PSA-Ag	Feline Prostate-specific antigen (PSA) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-PSA-Ab	Anti-Feline Prostate-specific antigen (PSA) antibody	antibody(Fc/mFc)
GMP-CAN-PSA-Ag	Canine Prostate-specific antigen (PSA) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-PSA-Ab	Anti-Canine Prostate-specific antigen (PSA) antibody	antibody(Fc/mFc)
GMP-FEL-p-tau181-Ag	Feline Phospho-tau181 (p-tau181) antigen	protein
GMP-FEL-p-tau181-Ab	Anti-Feline Phospho-tau181 (p-tau181) antibody	antibody(Fc/mFc)
GMP-CAN-p-tau181-Ag	Canine Phospho-tau181 (p-tau181) antigen	protein
GMP-CAN-p-tau181-Ab	Anti-Canine Phospho-tau181 (p-tau181) antibody	antibody(Fc/mFc)
GMP-FEL-p-tau202-205-Ag	Feline Phospho-tau202 and 205 (p-tau202 and 205) antigen	protein
GMP-FEL-p-tau202-205-Ab	Anti-Feline Phospho-tau202 and 205 (p-tau202 and 205) antibody	antibody(Fc/mFc)
GMP-CAN-p-tau202-205-Ag	Canine Phospho-tau202 and 205 (p-tau202 and 205) antigen	protein
GMP-CAN-p-tau202-205-Ab	Anti-Canine Phospho-tau202 and 205 (p-tau202 and 205) antibody	antibody(Fc/mFc)
GMP-FEL-p-tau212-214-Ag	Feline Phospho-tau212 and 214 (p-tau212 and 214) antigen	protein
GMP-FEL-p-tau212-214-Ab	Anti-Feline Phospho-tau212 and 214 (p-tau212 and 214) antibody	antibody(Fc/mFc)
GMP-CAN-p-tau212-214-Ag	Canine Phospho-tau212 and 214 (p-tau212 and 214) antigen	protein

Cat.No	Products Name	Products Type
GMP-CAN-p-tau212-214-Ab	Anti-Canine Phospho-tau212 and 214 (p-tau212 and 214) antibody	antibody(Fc/mFc)
GMP-FEL-p-tau217-Ag	Feline Phospho-tau217 (p-tau217) antigen	protein
GMP-FEL-p-tau217-Ab	Anti-Feline Phospho-tau217 (p-tau217) antibody	antibody(Fc/mFc)
GMP-CAN-p-tau217-Ag	Canine Phospho-tau217 (p-tau217) antigen	protein
GMP-CAN-p-tau217-Ab	Anti-Canine Phospho-tau217 (p-tau217) antibody	antibody(Fc/mFc)
GMP-FEL-p-tau231-Ag	Feline Phospho-tau231 (p-tau231) antigen	protein
GMP-FEL-p-tau231-Ab	Anti-Feline Phospho-tau231 (p-tau231) antibody	antibody(Fc/mFc)
GMP-CAN-p-tau231-Ag	Canine Phospho-tau231 (p-tau231) antigen	protein
GMP-CAN-p-tau231-Ab	Anti-Canine Phospho-tau231 (p-tau231) antibody	antibody(Fc/mFc)
GMP-FEL-p-tau413-Ag	Feline Phospho-tau413 (p-tau413) antigen	protein
GMP-FEL-p-tau413-Ab	Anti-Feline Phospho-tau413 (p-tau413) antibody	antibody(Fc/mFc)
GMP-CAN-p-tau413-Ag	Canine Phospho-tau413 (p-tau413) antigen	protein
GMP-CAN-p-tau413-Ab	Anti-Canine Phospho-tau413 (p-tau413) antibody	antibody(Fc/mFc)
GMP-FEL-p-tau422-Ag	Feline Phospho-tau422 (p-tau422) antigen	protein
GMP-FEL-p-tau422-Ab	Anti-Feline Phospho-tau422 (p-tau422) antibody	antibody(Fc/mFc)
GMP-CAN-p-tau422-Ag	Canine Phospho-tau422 (p-tau422) antigen	protein
GMP-CAN-p-tau422-Ab	Anti-Canine Phospho-tau422 (p-tau422) antibody	antibody(Fc/mFc)
GMP-FEL-PTH-Ag	Feline parathyroid hormone (PTH) antigen	protein
GMP-FEL-PTH-Ab	Anti-Feline parathyroid hormone (PTH) antibody	antibody(Fc/mFc)
GMP-CAN-PTH-Ag	Canine parathyroid hormone (PTH) antigen	protein
GMP-CAN-PTH-Ab	Anti-Canine parathyroid hormone (PTH) antibody	antibody(Fc/mFc)
GMP-CAN-Relaxin-Ag	Canine RLN (Relaxin) antigen	protein

Cat.No	Products Name	Products Type
GMP-CAN-Relaxin-Ab	Anti-Canine RLN (Relaxin) antibody	antibody(Fc/mFc)
GMP-FEL-relaxin-Ag	Feline RLN (Relaxin) antigen	protein
GMP-FEL-relaxin-Ab	Anti-Feline RLN (Relaxin) antibody	antibody(Fc/mFc)
GMP-CAN-RLN-Ag	Canine RLN (Relaxin) antigen	protein
GMP-CAN-RLN-Ab	Anti-Canine RLN (Relaxin) antibody	antibody(Fc/mFc)
GMP-FEL-S100-A8-Ag	Feline Calprotectin (S100-A8) antigen	protein
GMP-FEL-S100-A8-Ab	Anti-Feline Calprotectin (S100-A8) antibody	antibody(Fc/mFc)
GMP-CAN-S100-A8-Ag	Canine Calprotectin (S100-A8) antigen	protein
GMP-CAN-S100-A8-Ab	Anti-Canine Calprotectin (S100-A8) antibody	antibody(Fc/mFc)
GMP-FEL-S100-A9-Ag	Feline Calprotectin heavy chain (S100-A9) antigen	protein
GMP-FEL-S100-A9-Ab	Anti-Feline Calprotectin heavy chain (S100-A9) antibody	antibody(Fc/mFc)
GMP-CAN-S100-A9-Ag	Canine Calprotectin heavy chain (S100-A9) antigen	protein
GMP-CAN-S100-A9-Ab	Anti-Canine Calprotectin heavy chain (S100-A9) antibody	antibody(Fc/mFc)
GMP-FEL-S100B-Ag	Feline S100 proteins (S100B) antigen	protein
GMP-FEL-S100B-Ab	Anti-Feline S100 proteins (S100B) antibody	antibody(Fc/mFc)
GMP-CAN-S100B-Ag	Canine S100 proteins (S100B) antigen	protein
GMP-CAN-S100B-Ab	Anti-Canine S100 proteins (S100B) antibody	antibody(Fc/mFc)
GMP-FEL-SAHH-Ag	Feline S-adenosylhomocysteine (AdoHcy) hydrolase (SAHH) antigen	protein
GMP-FEL-SAHH-Ab	Anti-Feline S-adenosylhomocysteine (AdoHcy) hydrolase (SAHH) antibody	antibody(Fc/mFc)
GMP-CAN-SAHH-Ag	Canine S-adenosylhomocysteine (AdoHcy) hydrolase (SAHH) antigen	protein
GMP-CAN-SAHH-Ab	Anti-Canine S-adenosylhomocysteine (AdoHcy) hydrolase (SAHH) antibody	antibody(Fc/mFc)
GMP-FEL-SCC-Ag-Ag	Feline squamous cell carcinoma antigen (SCC Ag) antigen	protein

Cat.No	Products Name	Products Type
GMP-FEL-SCC-Ag-Ab	Anti-Feline squamous cell carcinoma antigen (SCC Ag) antibody	antibody(Fc/mFc)
GMP-CAN-SCC-Ag-Ag	Canine squamous cell carcinoma antigen (SCC Ag) antigen	protein
GMP-CAN-SCC-Ag-Ab	Anti-Canine squamous cell carcinoma antigen (SCC Ag) antibody	antibody(Fc/mFc)
GMP-FEL-Serum-amyloid-A-Ag	Feline SAA (Serum amyloid A) antigen	protein
GMP-FEL-Serum-amyloid-A-Ab	Anti-Feline SAA (Serum amyloid A) antibody	antibody(Fc/mFc)
GMP-CAN-Serum-amyloid-A-Ag	Canine SAA (Serum amyloid A) antigen	protein
GMP-CAN-Serum-amyloid-A-Ab	Anti-Canine SAA (Serum amyloid A) antibody	antibody(Fc/mFc)
GMP-FEL-sFlt-1-Ag	Feline soluble vascular endothelial growth factor receptor-2 (sFlt-1) antigen	protein
GMP-FEL-sFlt-1-Ab	Anti-Feline soluble vascular endothelial growth factor receptor-2 (sFlt-1) antibody	antibody(Fc/mFc)
GMP-CAN-sFlt-1-Ag	Canine soluble vascular endothelial growth factor receptor-2 (sFlt-1) antigen	protein
GMP-CAN-sFlt-1-Ab	Anti-Canine soluble vascular endothelial growth factor receptor-2 (sFlt-1) antibody	antibody(Fc/mFc)
GMP-FEL-SPA-Ag	Feline Staphylococcus protein A (SPA) antigen	protein
GMP-FEL-SPA-Ab	Anti-Feline Staphylococcus protein A (SPA) antibody	antibody(Fc/mFc)
GMP-CAN-SPA-Ag	Canine Staphylococcus protein A (SPA) antigen	protein
GMP-CAN-SPA-Ab	Anti-Canine Staphylococcus protein A (SPA) antibody	antibody(Fc/mFc)
GMP-FEL-Spirochaeta-Ag	Feline Spirochaeta antigen	protein
GMP-FEL-Spirochaeta-Ab	Anti-Feline Spirochaeta antibody	antibody(Fc/mFc)
GMP-CAN-Spirochaeta-Ag	Canine Spirochaeta antigen	protein
GMP-CAN-Spirochaeta-Ab	Anti-Canine Spirochaeta antibody	antibody(Fc/mFc)
GMP-FEL-sST2-Ag	Feline soluble ST2 (sST2) antigen	protein
GMP-FEL-sST2-Ab	Anti-Feline soluble ST2 (sST2) antibody	antibody(Fc/mFc)
GMP-CAN-sST2-Ag	Canine soluble ST2 (sST2) antigen	protein

Cat.No	Products Name	Products Type
GMP-CAN-sST2-Ab	Anti-Canine soluble ST2 (sST2) antibody	antibody(Fc/mFc)
GMP-FEL-sTREM2-Ag	Feline soluble fragment of triggering receptor expressed on myeloid cells 2 (sTREM2) antigen	protein
GMP-FEL-sTREM2-Ab	Anti-Feline soluble fragment of triggering receptor expressed on myeloid cells 2 (sTREM2) antibody	antibody(Fc/mFc)
GMP-CAN-sTREM2-Ag	Canine soluble fragment of triggering receptor expressed on myeloid cells 2 (sTREM2) antigen	protein
GMP-CAN-sTREM2-Ab	Anti-Canine soluble fragment of triggering receptor expressed on myeloid cells 2 (sTREM2) antibody	antibody(Fc/mFc)
GMP-FEL-Tau-Ag	Feline tau proteins (Tau) antigen	protein
GMP-FEL-Tau-Ab	Anti-Feline tau proteins (Tau) antibody	antibody(Fc/mFc)
GMP-CAN-Tau-Ag	Canine tau proteins (Tau) antigen	protein
GMP-CAN-Tau-Ab	Anti-Canine tau proteins (Tau) antibody	antibody(Fc/mFc)
GMP-FEL-TG-Ag	Feline Thyroglobulin (TG) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-TG-Ab	Anti-Feline Thyroglobulin (TG) antibody	antibody(Fc/mFc)
GMP-CAN-TG-Ag	Canine Thyroglobulin (TG) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-TG-Ab	Anti-Canine Thyroglobulin (TG) antibody	antibody(Fc/mFc)
GMP-FEL-TIMP-1-Ag	Feline Matrix Metalloproteinase Inhibitor-1 (TIMP-1) antigen	protein
GMP-FEL-TIMP-1-Ab	Anti-Feline Matrix Metalloproteinase Inhibitor-1 (TIMP-1) antibody	antibody(Fc/mFc)
GMP-CAN-TIMP-1-Ag	Canine Matrix Metalloproteinase Inhibitor-1 (TIMP-1) antigen	protein
GMP-CAN-TIMP-1-Ab	Anti-Canine Matrix Metalloproteinase Inhibitor-1 (TIMP-1) antibody	antibody(Fc/mFc)
GMP-FEL-TK1-Ag	Feline Thymidine Kinase 1 (TK1) antigen	protein
GMP-FEL-TK1-Ab	Anti-Feline Thymidine Kinase 1 (TK1) antibody	antibody(Fc/mFc)
GMP-CAN-TK1-Ag	Canine Thymidine Kinase 1 (TK1) antigen	protein
GMP-CAN-TK1-Ab	Anti-Canine Thymidine Kinase 1 (TK1) antibody	antibody(Fc/mFc)
GMP-FEL-TNFa-Ag	Feline tumor necrosis factor alpha (TNFα) antigen	protein

Cat.No	Products Name	Products Type
GMP-FEL-TNFa-Ab	Anti-Feline tumor necrosis factor alpha (TNFα) antibody	antibody(Fc/mFc)
GMP-CAN-TNFa-Ag	Canine tumor necrosis factor alpha (TNFα) antigen	protein
GMP-CAN-TNFa-Ab	Anti-Canine tumor necrosis factor alpha (TNFα) antibody	antibody(Fc/mFc)
GMP-FEL-TPO-Ag	Feline Thyroid peroxidase (TPO) antigen	protein
GMP-FEL-TPO-Ab	Anti-Feline Thyroid peroxidase (TPO) antibody	antibody(Fc/mFc)
GMP-CAN-TPO-Ag	Canine Thyroid peroxidase (TPO) antigen	protein
GMP-CAN-TPO-Ab	Anti-Canine Thyroid peroxidase (TPO) antibody	antibody(Fc/mFc)
GMP-FEL-TrxR-Ag	Feline thioredoxin reductase (TrxR) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-TrxR-Ab	Anti-Feline thioredoxin reductase (TrxR) antibody	antibody(Fc/mFc)
GMP-CAN-TrxR-Ag	Canine thioredoxin reductase (TrxR) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-TrxR-Ab	Anti-Canine thioredoxin reductase (TrxR) antibody	antibody(Fc/mFc)
GMP-FEL-TSH-Ag	Feline Thyroid-stimulating hormone (TSH) antigen	protein
GMP-FEL-TSH-Ab	Anti-Feline Thyroid-stimulating hormone (TSH) antibody	antibody(Fc/mFc)
GMP-CAN-TSH-Ag	Canine Thyroid-stimulating hormone (TSH) antigen	protein
GMP-CAN-TSH-Ab	Anti-Canine Thyroid-stimulating hormone (TSH) antibody	antibody(Fc/mFc)
GMP-FEL-VDBP-Ag	Feline vitamin D binding protein (VDBP) antigen	protein
GMP-FEL-VDBP-Ab	Anti-Feline vitamin D binding protein (VDBP) antibody	antibody(Fc/mFc)
GMP-CAN-VDBP-Ag	Canine vitamin D binding protein (VDBP) antigen	protein
GMP-CAN-VDBP-Ab	Anti-Canine vitamin D binding protein (VDBP) antibody	antibody(Fc/mFc)
GMP-FEL-VEGF-Ag	Feline Vascular Endothelial Growth Factor (VEGF) antigen	protein
GMP-FEL-VEGF-Ab	Anti-Feline Vascular Endothelial Growth Factor (VEGF) antibody	antibody(Fc/mFc)
GMP-CAN-VEGF-Ag	Canine Vascular Endothelial Growth Factor (VEGF) antigen	protein

Cat.No	Products Name	Products Type
GMP-CAN-VEGF-Ab	Anti-Canine Vascular Endothelial Growth Factor (VEGF) antibody	antibody(Fc/mFc)
GMP-FEL-VLP-1-Ag	Feline Visinin-like protein 1 (VLP-1) antigen	protein
GMP-FEL-VLP-1-Ab	Anti-Feline Visinin-like protein 1 (VLP-1) antibody	antibody(Fc/mFc)
GMP-CAN-VLP-1-Ag	Canine Visinin-like protein 1 (VLP-1) antigen	protein
GMP-CAN-VLP-1-Ab	Anti-Canine Visinin-like protein 1 (VLP-1) antibody	antibody(Fc/mFc)
GMP-FEL-YKL-40-Ag	Feline Chitinase-3-like protein 1 (YKL-40) antigen	protein
GMP-FEL-YKL-40-Ab	Anti-Feline Chitinase-3-like protein 1 (YKL-40) antibody	antibody(Fc/mFc)
GMP-CAN-YKL-40-Ag	Canine Chitinase-3-like protein 1 (YKL-40) antigen	protein
GMP-CAN-YKL-40-Ab	Anti-Canine Chitinase-3-like protein 1 (YKL-40) antibody	antibody(Fc/mFc)
GMP-FEL-a-defensin-1-Ag	Feline alpha defensin 1 (α -defensin 1) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-a-defensin-1-Ab	Anti-Feline alpha defensin 1 (α -defensin 1) antibody	antibody(Fc/mFc)
GMP-CAN-a-defensin-1-Ag	Canine alpha defensin 1 (α -defensin 1) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-a-defensin-1-Ab	Anti-Canine alpha defensin 1 (α -defensin 1) antibody	antibody(Fc/mFc)
GMP-FEL-a-defensin-5-Ag	Feline alpha defensin 5 (α -defensin 5) antigen	peptide-(naked/BSA/OVA)
GMP-FEL-a-defensin-5-Ab	Anti-Feline alpha defensin 5 (α -defensin 5) antibody	antibody(Fc/mFc)
GMP-CAN-a-defensin-5-Ag	Canine alpha defensin 5 (α -defensin 5) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-a-defensin-5-Ab	Anti-Canine alpha defensin 5 (α -defensin 5) antibody	antibody(Fc/mFc)
GMP-FEL-a-synuclein-Ag	Feline Alpha-synuclein (α -synuclein) antigen	protein
GMP-FEL-a-synuclein-Ab	Anti-Feline Alpha-synuclein (α -synuclein) antibody	antibody(Fc/mFc)
GMP-CAN-a-synuclein-Ag	Canine Alpha-synuclein (α -synuclein) antigen	protein
GMP-CAN-a-synuclein-Ab	Anti-Canine Alpha-synuclein (α -synuclein) antibody	antibody(Fc/mFc)
GMP-FEL-B-CTx-Ag	Feline β -isomerized C-terminal telopeptides (β -CTx) antigen	peptide-(naked/BSA/OVA)

Cat.No	Products Name	Products Type
GMP-FEL-B-CTx-Ab	Anti-Feline β -isomerized C-terminal telopeptides (β -CTx) antibody	antibody(Fc/mFc)
GMP-CAN-B-CTx-Ag	Canine β -isomerized C-terminal telopeptides (β -CTx) antigen	peptide-(naked/BSA/OVA)
GMP-CAN-B-CTx-Ab	Anti-Canine β -isomerized C-terminal telopeptides (β -CTx) antibody	antibody(Fc/mFc)

Swine

Cat.No	Products Name	Products Type
GMP-VT-P026-Tg003-Ag	Leptospira Recombinant Envelope-Related Protein Y antigen	protein
GMP-VT-P026-Tg003-Ab	Anti-Leptospira Recombinant Envelope-Related Protein Y antibody	antibody(Fc/mFc)
GMP-VT-P046-Tg002-Ag	Chlamydiaceae Major Outer Membrane Protein antigen	protein
GMP-VT-P046-Tg002-Ab	Anti-Chlamydiaceae Major Outer Membrane Protein antibody	antibody(Fc/mFc)
GMP-VT-P054-Tg002-Ag	Mycoplasma pneumoniae antigen	protein
GMP-VT-P054-Tg002-Ab	Anti-Mycoplasma pneumoniae antibody	antibody(Fc/mFc)
GMP-VT-P094-Tg001-Ag	Mycoplasma hyosynoviae antigen	protein
GMP-VT-P094-Tg001-Ab	Anti-Mycoplasma hyosynoviae antibody	antibody(Fc/mFc)
GMP-VT-P095-Tg001-Ag	SVV Viral Protein 2 antigen	protein
GMP-VT-P095-Tg001-Ab	Anti-SVV Viral Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P095-Tg002-Ag	SVV Viral Protein 1 antigen	protein
GMP-VT-P095-Tg002-Ab	Anti-SVV Viral Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P096-Tg001-Ag	Spirochetes Bacterial heme lipoprotein 29.7 antigen	protein
GMP-VT-P096-Tg001-Ab	Anti-Spirochetes Bacterial heme lipoprotein 29.7 antibody	antibody(Fc/mFc)
GMP-VT-P097-Tg001-Ag	Blue ear virus Open Reading Frame 7 Protein antigen	protein
GMP-VT-P097-Tg001-Ab	Anti-Blue ear virus Open Reading Frame 7 Protein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P098-Tg001-Ag	Swine enteric coronavirus Spike Protein antigen	protein
GMP-VT-P098-Tg001-Ab	Anti-Swine enteric coronavirus Spike Protein antibody	antibody(Fc/mFc)
GMP-VT-P098-Tg002-Ag	Swine enteric coronavirus Nucleocapsid Protein antigen	protein
GMP-VT-P098-Tg002-Ab	Anti-Swine enteric coronavirus Nucleocapsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P099-Tg001-Ag	Brucella suis antigen	protein
GMP-VT-P099-Tg001-Ab	Anti-Brucella suis antibody	antibody(Fc/mFc)
GMP-VT-P101-Tg001-Ag	Pasteurella multocida Pasteurella Multocida Toxin antigen	protein
GMP-VT-P101-Tg001-Ab	Anti-Pasteurella multocida Pasteurella Multocida Toxin antibody	antibody(Fc/mFc)
GMP-VT-P102-Tg001-Ag	Brachyspira hyodysenteriae Bacterial heme lipoprotein 29.7 antigen	protein
GMP-VT-P102-Tg001-Ab	Anti-Brachyspira hyodysenteriae Bacterial heme lipoprotein 29.7 antibody	antibody(Fc/mFc)
GMP-VT-P103-Tg001-Ag	Ecephalitis antigen	protein
GMP-VT-P103-Tg001-Ab	Anti-Ecephalitis antibody	antibody(Fc/mFc)
GMP-VT-P104-Tg001-Ag	PDCoV Membrane Protein antigen	protein
GMP-VT-P104-Tg001-Ab	Anti-PDCoV Membrane Protein antibody	antibody(Fc/mFc)
GMP-VT-P104-Tg002-Ag	PDCoV Nucleocapsid Protein antigen	protein
GMP-VT-P104-Tg002-Ab	Anti-PDCoV Nucleocapsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P291-Tg001-Ag	JEV Envelope Protein antigen	protein
GMP-VT-P291-Tg001-Ab	Anti-JEV Envelope Protein antibody	antibody(Fc/mFc)
GMP-VT-P291-Tg002-Ag	JEV Non-structural Protein 1 antigen	protein
GMP-VT-P291-Tg002-Ab	Anti-JEV Non-structural Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P105-Tg001-Ag	PRV Glycoprotein E antigen	protein
GMP-VT-P105-Tg001-Ab	Anti-PRV Glycoprotein E antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P105-Tg002-Ag	PRV Glycoprotein III antigen	protein
GMP-VT-P105-Tg002-Ab	Anti-PRV Glycoprotein III antibody	antibody(Fc/mFc)
GMP-VT-P105-Tg003-Ag	PRV Glycoprotein 50 antigen	protein
GMP-VT-P105-Tg003-Ab	Anti-PRV Glycoprotein 50 antibody	antibody(Fc/mFc)
GMP-VT-P292-Tg001-Ag	PPV Viral Protein 2 antigen	protein
GMP-VT-P292-Tg001-Ab	Anti-PPV Viral Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P292-Tg002-Ag	PPV Non-structural Protein 1 antigen	protein
GMP-VT-P292-Tg002-Ab	Anti-PPV Non-structural Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P292-Tg003-Ag	PPV antigen	protein
GMP-VT-P292-Tg003-Ab	Anti-PPV antibody	antibody(Fc/mFc)
GMP-VT-P106-Tg001-Ag	CSFV Glycoprotein E2 antigen	protein
GMP-VT-P106-Tg001-Ab	Anti-CSFV Glycoprotein E2 antibody	antibody(Fc/mFc)
GMP-VT-P107-Tg001-Ag	ASFV pp62 antigen	protein
GMP-VT-P107-Tg001-Ab	Anti-ASFV pp62 antibody	antibody(Fc/mFc)
GMP-VT-P107-Tg002-Ag	ASFV p54 antigen	protein
GMP-VT-P107-Tg002-Ab	Anti-ASFV p54 antibody	antibody(Fc/mFc)
GMP-VT-P107-Tg003-Ag	ASFV p30 antigen	protein
GMP-VT-P107-Tg003-Ab	Anti-ASFV p30 antibody	antibody(Fc/mFc)
GMP-VT-P107-Tg004-Ag	ASFV p72 antigen	protein
GMP-VT-P107-Tg004-Ab	Anti-ASFV p72 antibody	antibody(Fc/mFc)
GMP-VT-P107-Tg005-Ag	ASFV p32 antigen	protein
GMP-VT-P107-Tg005-Ab	Anti-ASFV p32 antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P107-Tg006-Ag	ASFV p62 antigen	protein
GMP-VT-P107-Tg006-Ab	Anti-ASFV p62 antibody	antibody(Fc/mFc)
GMP-VT-P285-Tg001-Ag	E.coli F18 F18 fimbrial adhesin FedF antigen	protein
GMP-VT-P285-Tg001-Ab	Anti-E.coli F18 F18 fimbrial adhesin FedF antibody	antibody(Fc/mFc)
GMP-VT-P285-Tg002-Ag	E.coli F18 F18 fimbrial adhesin FaeG antigen	protein
GMP-VT-P285-Tg002-Ab	Anti-E.coli F18 F18 fimbrial adhesin FaeG antibody	antibody(Fc/mFc)
GMP-VT-P285-Tg003-Ag	E.coli F18 F18 fimbrial adhesin FanC antigen	protein
GMP-VT-P285-Tg003-Ab	Anti-E.coli F18 F18 fimbrial adhesin FanC antibody	antibody(Fc/mFc)
GMP-VT-P285-Tg004-Ag	E.coli F18 F18 fimbrial adhesin FasA antigen	protein
GMP-VT-P285-Tg004-Ab	Anti-E.coli F18 F18 fimbrial adhesin FasA antibody	antibody(Fc/mFc)
GMP-VT-P285-Tg005-Ag	E.coli F18 Fimbrial protein Fim41a antigen	protein
GMP-VT-P285-Tg005-Ab	Anti-E.coli F18 Fimbrial protein Fim41a antibody	antibody(Fc/mFc)
GMP-VT-P108-Tg001-Ag	PEDV Nucleoprotein antigen	protein
GMP-VT-P108-Tg001-Ab	Anti-PEDV Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P109-Tg001-Ag	TGEV Spike Glycoprotein antigen	protein
GMP-VT-P109-Tg001-Ab	Anti-TGEV Spike Glycoprotein antibody	antibody(Fc/mFc)
GMP-VT-P110-Tg001-Ag	RVC Viral Protein 6 antigen	protein
GMP-VT-P110-Tg001-Ab	Anti-RVC Viral Protein 6 antibody	antibody(Fc/mFc)
GMP-VT-P044-Tg002-Ag	RVA Viral Protein 6 antigen	protein
GMP-VT-P044-Tg002-Ab	Anti-RVA Viral Protein 6 antibody	antibody(Fc/mFc)
GMP-VT-P111-Tg001-Ag	PRRSV Nucleocapsid Protein antigen	protein
GMP-VT-P111-Tg001-Ab	Anti-PRRSV Nucleocapsid Protein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P112-Tg001-Ag	Streptococcus suis Surface Antigen One M antigen	protein
GMP-VT-P112-Tg001-Ab	Anti-Streptococcus suis Surface Antigen One M antibody	antibody(Fc/mFc)
GMP-VT-P113-Tg001-Ag	PCV Capsid Protein antigen	protein
GMP-VT-P113-Tg001-Ab	Anti-PCV Capsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P114-Tg001-Ag	PCV2 Capsid Protein antigen	protein
GMP-VT-P114-Tg001-Ab	Anti-PCV2 Capsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P115-Tg001-Ag	Mycoplasma hyopneumoniae Mhp366 Protein antigen	protein
GMP-VT-P115-Tg001-Ab	Anti-Mycoplasma hyopneumoniae Mhp366 Protein antibody	antibody(Fc/mFc)
GMP-VT-P115-Tg002-Ag	Mycoplasma hyopneumoniae P65 antigen	protein
GMP-VT-P115-Tg002-Ab	Anti-Mycoplasma hyopneumoniae P65 antibody	antibody(Fc/mFc)
GMP-VT-P115-Tg003-Ag	Mycoplasma hyopneumoniae P102 antigen	protein
GMP-VT-P115-Tg003-Ab	Anti-Mycoplasma hyopneumoniae P102 antibody	antibody(Fc/mFc)
GMP-VT-P115-Tg004-Ag	Mycoplasma hyopneumoniae P46 antigen	protein
GMP-VT-P115-Tg004-Ab	Anti-Mycoplasma hyopneumoniae P46 antibody	antibody(Fc/mFc)
GMP-VT-P116-Tg001-Ag	Mycoplasma hyorhinis Lawsonia surface antigen A antigen	protein
GMP-VT-P116-Tg001-Ab	Anti-Mycoplasma hyorhinis Lawsonia surface antigen A antibody	antibody(Fc/mFc)
GMP-VT-P117-Tg001-Ag	APP ApxII Toxin antigen	protein
GMP-VT-P117-Tg001-Ab	Anti-APP ApxII Toxin antibody	antibody(Fc/mFc)
GMP-VT-P118-Tg001-Ag	SIV/SIV-A Nucleoprotein antigen	protein
GMP-VT-P118-Tg001-Ab	Anti-SIV/SIV-A Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P119-Tg001-Ag	Haemophilus parasuis Transferrin Binding Protein B antigen	protein
GMP-VT-P119-Tg001-Ab	Anti-Haemophilus parasuis Transferrin Binding Protein B antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P119-Tg002-Ag	Haemophilus parasuis Apd Protein antigen	protein
GMP-VT-P119-Tg002-Ab	Anti-Haemophilus parasuis Apd Protein antibody	antibody(Fc/mFc)
GMP-VT-P120-Tg001-Ag	ADV Glycoprotein E antigen	protein
GMP-VT-P120-Tg001-Ab	Anti-ADV Glycoprotein E antibody	antibody(Fc/mFc)
GMP-VT-P121-Tg001-Ag	APP serotype 10 ApxII Toxin antigen	protein
GMP-VT-P121-Tg001-Ab	Anti-APP serotype 10 ApxII Toxin antibody	antibody(Fc/mFc)
GMP-VT-P122-Tg001-Ag	APP serotype 12 ApxII Toxin antigen	protein
GMP-VT-P122-Tg001-Ab	Anti-APP serotype 12 ApxII Toxin antibody	antibody(Fc/mFc)
GMP-VT-P123-Tg001-Ag	APP serotype 2 Outer Membrane Protein antigen	protein
GMP-VT-P123-Tg001-Ab	Anti-APP serotype 2 Outer Membrane Protein antibody	antibody(Fc/mFc)
GMP-VT-P124-Tg001-Ag	APP serotype 5 ApxI Toxin antigen	protein
GMP-VT-P124-Tg001-Ab	Anti-APP serotype 5 ApxI Toxin antibody	antibody(Fc/mFc)
GMP-VT-P125-Tg001-Ag	APP serotype 1to12 ApxIV Toxin antigen	protein
GMP-VT-P125-Tg001-Ab	Anti-APP serotype 1to12 ApxIV Toxin antibody	antibody(Fc/mFc)
GMP-VT-P126-Tg001-Ag	APP serotype1, 9, 11 ApxI Toxin antigen	protein
GMP-VT-P126-Tg001-Ab	Anti-APP serotype1, 9, 11 ApxI Toxin antibody	antibody(Fc/mFc)
GMP-VT-P127-Tg001-Ag	APP serotype 3 Cysteine Synthase G antigen	protein
GMP-VT-P127-Tg001-Ab	Anti-APP serotype 3 Cysteine Synthase G antibody	antibody(Fc/mFc)
GMP-VT-P128-Tg001-Ag	APP serotype 4 ApxIII Toxin antigen	protein
GMP-VT-P128-Tg001-Ab	Anti-APP serotype 4 ApxIII Toxin antibody	antibody(Fc/mFc)
GMP-VT-P129-Tg001-Ag	Brachyspira pilosicoli ClpX Protease antigen	protein
GMP-VT-P129-Tg001-Ab	Anti-Brachyspira pilosicoli ClpX Protease antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P283-Tg001-Ag	E.coli F4 Fimbrial Adhesin FaeG antigen	protein
GMP-VT-P283-Tg001-Ab	Anti-E.coli F4 Fimbrial Adhesin FaeG antibody	antibody(Fc/mFc)
GMP-VT-P286-Tg001-Ag	E.coli F41 F41 Fimbrial Adhesin antigen	protein
GMP-VT-P286-Tg001-Ab	Anti-E.coli F41 F41 Fimbrial Adhesin antibody	antibody(Fc/mFc)
GMP-VT-P284-Tg002-Ag	E.coli F5 Fimbrial Adhesin K99 antigen	protein
GMP-VT-P284-Tg002-Ab	Anti-E.coli F5 Fimbrial Adhesin K99 antibody	antibody(Fc/mFc)
GMP-VT-P130-Tg001-Ag	Erysipelothrix rhusiopathiae Surface Protective Antigen A antigen	protein
GMP-VT-P130-Tg001-Ab	Anti-Erysipelothrix rhusiopathiae Surface Protective Antigen A antibody	antibody(Fc/mFc)
GMP-VT-P131-Tg001-Ag	Mycoplasma hyorhinis Variable Lipoprotein A-G antigen	protein
GMP-VT-P131-Tg001-Ab	Anti-Mycoplasma hyorhinis Variable Lipoprotein A-G antibody	antibody(Fc/mFc)
GMP-VT-P132-Tg001-Ag	PRCV Spike Protein antigen	protein
GMP-VT-P132-Tg001-Ab	Anti-PRCV Spike Protein antibody	antibody(Fc/mFc)
GMP-VT-P133-Tg001-Ag	SVDV Virus-Like Particles antigen	protein
GMP-VT-P133-Tg001-Ab	Anti-SVDV Virus-Like Particles antibody	antibody(Fc/mFc)
GMP-VT-P134-Tg001-Ag	NiV Nucleoprotein, Fusion Protein, Attachment Protein, Matrix Protein antigen	protein
GMP-VT-P134-Tg001-Ab	Anti-NiV Nucleoprotein, Fusion Protein, Attachment Protein, Matrix Protein antibody	antibody(Fc/mFc)

Horse

Cat.No	Products Name	Products Type
GMP-VT-P135-Tg001-Ag	Taylorella equigenitalis antigen	protein
GMP-VT-P135-Tg001-Ab	Anti-Taylorella equigenitalis antibody	antibody(Fc/mFc)
GMP-VT-P136-Tg001-Ag	EHV Glycoprotein D antigen	protein

Cat.No	Products Name	Products Type
GMP-VT-P136-Tg001-Ab	Anti-EHV Glycoprotein D antibody	antibody(Fc/mFc)
GMP-VT-P136-Tg002-Ag	EHV Glycoprotein E antigen	protein
GMP-VT-P136-Tg002-Ab	Anti-EHV Glycoprotein E antibody	antibody(Fc/mFc)
GMP-VT-P136-Tg003-Ag	EHV Glycoprotein B antigen	protein
GMP-VT-P136-Tg003-Ab	Anti-EHV Glycoprotein B antibody	antibody(Fc/mFc)
GMP-VT-P136-Tg004-Ag	EHV Glycoprotein G antigen	protein
GMP-VT-P136-Tg004-Ab	Anti-EHV Glycoprotein G antibody	antibody(Fc/mFc)
GMP-VT-P137-Tg001-Ag	EHV-1 Glycoprotein D antigen	protein
GMP-VT-P137-Tg001-Ab	Anti-EHV-1 Glycoprotein D antibody	antibody(Fc/mFc)
GMP-VT-P137-Tg002-Ag	EHV-1 Glycoprotein E antigen	protein
GMP-VT-P137-Tg002-Ab	Anti-EHV-1 Glycoprotein E antibody	antibody(Fc/mFc)
GMP-VT-P137-Tg003-Ag	EHV-1 Glycoprotein B antigen	protein
GMP-VT-P137-Tg003-Ab	Anti-EHV-1 Glycoprotein B antibody	antibody(Fc/mFc)
GMP-VT-P138-Tg001-Ag	EHV-4 Glycoprotein D antigen	protein
GMP-VT-P138-Tg001-Ab	Anti-EHV-4 Glycoprotein D antibody	antibody(Fc/mFc)
GMP-VT-P138-Tg002-Ag	EHV-4 Glycoprotein E antigen	protein
GMP-VT-P138-Tg002-Ab	Anti-EHV-4 Glycoprotein E antibody	antibody(Fc/mFc)
GMP-VT-P138-Tg003-Ag	EHV-4 Glycoprotein B antigen	protein
GMP-VT-P138-Tg003-Ab	Anti-EHV-4 Glycoprotein B antibody	antibody(Fc/mFc)
GMP-VT-P139-Tg001-Ag	EIAV Gag Precursor Protein antigen	protein
GMP-VT-P139-Tg001-Ab	Anti-EIAV Gag Precursor Protein antibody	antibody(Fc/mFc)
GMP-VT-P139-Tg002-Ag	EIAV P26 Protein antigen	protein

Cat.No	Products Name	Products Type
GMP-VT-P139-Tg002-Ab	Anti-EIAV P26 Protein antibody	antibody(Fc/mFc)
GMP-VT-P140-Tg001-Ag	AHSV African Horse Sickness Virus Viral Protein 7 antigen	protein
GMP-VT-P140-Tg001-Ab	Anti-AHSV African Horse Sickness Virus Viral Protein 7 antibody	antibody(Fc/mFc)
GMP-VT-P140-Tg002-Ag	AHSV Non-Structural Protein 3 antigen	protein
GMP-VT-P140-Tg002-Ab	Anti-AHSV Non-Structural Protein 3 antibody	antibody(Fc/mFc)
GMP-VT-P141-Tg001-Ag	EIV-A Nucleoprotein antigen	protein
GMP-VT-P141-Tg001-Ab	Anti-EIV-A Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P142-Tg001-Ag	EAV Major Envelope Protein antigen	protein
GMP-VT-P142-Tg001-Ab	Anti-EAV Major Envelope Protein antibody	antibody(Fc/mFc)
GMP-VT-P142-Tg002-Ag	EAV Glycoprotein L antigen	protein
GMP-VT-P142-Tg002-Ab	Anti-EAV Glycoprotein L antibody	antibody(Fc/mFc)
GMP-VT-P142-Tg003-Ag	EAV Multiple M Proteins antigen	protein
GMP-VT-P142-Tg003-Ab	Anti-EAV Multiple M Proteins antibody	antibody(Fc/mFc)
GMP-VT-P142-Tg004-Ag	EAV Nucleocapsid Proteins antigen	protein
GMP-VT-P142-Tg004-Ab	Anti-EAV Nucleocapsid Proteins antibody	antibody(Fc/mFc)
GMP-VT-P143-Tg001-Ag	Babesia caballi Rhostry-Associated Protein 1 antigen	protein
GMP-VT-P143-Tg001-Ab	Anti-Babesia caballi Rhostry-Associated Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P144-Tg001-Ag	Theileria equi Merozoite Surface Antigen-1 antigen	protein
GMP-VT-P144-Tg001-Ab	Anti-Theileria equi Merozoite Surface Antigen-1 antibody	antibody(Fc/mFc)
GMP-VT-P144-Tg002-Ag	Theileria equi Merozoite Surface Antigen-2 antigen	protein
GMP-VT-P144-Tg002-Ab	Anti-Theileria equi Merozoite Surface Antigen-2 antibody	antibody(Fc/mFc)
GMP-VT-P145-Tg001-Ag	Klebsiella pneumoniae & Pseudomonas aeruginosa Metallo-beta-lactamase type 2 antigen	protein

Cat.No	Products Name	Products Type
GMP-VT-P145-Tg001-Ab	Anti-Klebsiella pneumoniae & Pseudomonas aeruginosa Metallo-beta-lactamase type 2 antibody	antibody(Fc/mFc)
GMP-VT-P145-Tg002-Ag	Klebsiella pneumoniae & Pseudomonas aeruginosa Elastase antigen	protein
GMP-VT-P145-Tg002-Ab	Anti-Klebsiella pneumoniae & Pseudomonas aeruginosa Elastase antibody	antibody(Fc/mFc)
GMP-VT-P146-Tg001-Ag	Streptococcus equi Streptococcus equi M protein antigen	protein
GMP-VT-P146-Tg001-Ab	Anti-Streptococcus equi Streptococcus equi M protein antibody	antibody(Fc/mFc)
GMP-VT-P147-Tg001-Ag	EEEV, WEEV Glycoprotein E1 antigen	protein
GMP-VT-P147-Tg001-Ab	Anti-EEEV, WEEV Glycoprotein E1 antibody	antibody(Fc/mFc)
GMP-VT-P147-Tg002-Ag	EEEV, WEEV Glycoprotein E2 antigen	protein
GMP-VT-P147-Tg002-Ab	Anti-EEEV, WEEV Glycoprotein E2 antibody	antibody(Fc/mFc)
GMP-VT-P148-Tg001-Ag	Equine/Horse morbillivirus Matrix and Fusion Proteins antigen	protein
GMP-VT-P148-Tg001-Ab	Anti-Equine/Horse morbillivirus Matrix and Fusion Proteins antibody	antibody(Fc/mFc)
GMP-EQU-AGP-Ag	Equine AGP antigen	protein
GMP-EQU-AGP-Ab	Anti-Equine AGP antibody	antibody(Fc/mFc)
GMP-EQU-lactoferrin-Ag	Equine lactoferrin antigen	protein
GMP-EQU-lactoferrin-Ab	Anti-Equine lactoferrin antibody	antibody(Fc/mFc)
GMP-EQU-A1M-Ag	Equine Alpha-1-microglobulin (A1M) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-A1M-Ab	Anti-Equine Alpha-1-microglobulin (A1M) antibody	antibody(Fc/mFc)
GMP-EQU-ACE2/hFc-Ag	Equine Angiotensin Converting Enzyme 2 (ACE2/hFc) antigen	protein
GMP-EQU-ACE2/hFc-Ab	Anti-Equine Angiotensin Converting Enzyme 2 (ACE2/hFc) antibody	antibody(Fc/mFc)
GMP-EQU-ACE2/his-Ag	Equine Angiotensin Converting Enzyme 2 (His) (ACE2/his) antigen	protein
GMP-EQU-ACE2/his-Ab	Anti-Equine Angiotensin Converting Enzyme 2 (His) (ACE2/his) antibody	antibody(Fc/mFc)
GMP-EQU-ACRV1-Ag	Equine acrosome protein (ACRV1) antigen	protein

Cat.No	Products Name	Products Type
GMP-EQU-ACRV1-Ab	Anti-Equine acrosome protein (ACRV1) antibody	antibody(Fc/mFc)
GMP-EQU-ACTH-Ag	Equine adrenocorticotrophic hormone (ACTH) antigen	protein
GMP-EQU-ACTH-Ab	Anti-Equine adrenocorticotrophic hormone (ACTH) antibody	antibody(Fc/mFc)
GMP-EQU-AD7c-NTP-Ag	Equine AD-associated neuronal thread protein (AD7c-NTP) antigen	protein
GMP-EQU-AD7c-NTP-Ab	Anti-Equine AD-associated neuronal thread protein (AD7c-NTP) antibody	antibody(Fc/mFc)
GMP-EQU-ADPN-Ag	Equine Adiponectin (ADPN) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-ADPN-Ab	Anti-Equine Adiponectin (ADPN) antibody	antibody(Fc/mFc)
GMP-EQU-AFP-Ag	Equine Alpha-fetoprotein (AFP) antigen	protein
GMP-EQU-AFP-Ab	Anti-Equine Alpha-fetoprotein (AFP) antibody	antibody(Fc/mFc)
GMP-EQU-AI-Ag	Equine Angiotensin (AI) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-AI-Ab	Anti-Equine Angiotensin (AI) antibody	antibody(Fc/mFc)
GMP-EQU-albumin-Ag	Equine albumin protein (albumin) antigen	protein
GMP-EQU-albumin-Ab	Anti-Equine albumin protein (albumin) antibody	antibody(Fc/mFc)
GMP-EQU-AMH-Ag	Equine Anti-Müllerian hormone (AMH) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-AMH-Ab	Anti-Equine Anti-Müllerian hormone (AMH) antibody	antibody(Fc/mFc)
GMP-EQU-ANP-Ag	Equine atrial natriuretic peptide (ANP) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-ANP-Ab	Anti-Equine atrial natriuretic peptide (ANP) antibody	antibody(Fc/mFc)
GMP-EQU-ApoA-Ag	Equine Apolipoprotein A (ApoA) antigen	protein
GMP-EQU-ApoA-Ab	Anti-Equine Apolipoprotein A (ApoA) antibody	antibody(Fc/mFc)
GMP-EQU-ApoE-Ag	Equine Apolipoprotein E (ApoE) antigen	protein
GMP-EQU-ApoE-Ab	Anti-Equine Apolipoprotein E (ApoE) antibody	antibody(Fc/mFc)
GMP-EQU-AQP4-Ag	Equine aquaporin-4 (AQP4) antigen	protein

Cat.No	Products Name	Products Type
GMP-EQU-AQP4-Ab	Anti-Equine aquaporin-4 (AQP4) antibody	antibody(Fc/mFc)
GMP-EQU-AB38-Ag	Equine beta-amyloid 38 (A β 38) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-AB38-Ab	Anti-Equine beta-amyloid 38 (A β 38) antibody	antibody(Fc/mFc)
GMP-EQU-AB40-Ag	Equine beta-amyloid 40 (A β 40) antigen	protein
GMP-EQU-AB40-Ab	Anti-Equine beta-amyloid 40 (A β 40) antibody	antibody(Fc/mFc)
GMP-EQU-AB42-Ag	Equine beta-amyloid 42 (A β 42) antigen	protein
GMP-EQU-AB42-Ab	Anti-Equine beta-amyloid 42 (A β 42) antibody	antibody(Fc/mFc)
GMP-EQU-B2M-Ag	Equine Beta-2 microglobulin (B2M) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-B2M-Ab	Anti-Equine Beta-2 microglobulin (B2M) antibody	antibody(Fc/mFc)
GMP-EQU-BNP-Ag	Equine B-type natriuretic peptide (BNP) antigen	protein
GMP-EQU-BNP-Ab	Anti-Equine B-type natriuretic peptide (BNP) antibody	antibody(Fc/mFc)
GMP-EQU-CA-125-Ag	Equine cancer antigen 125 (CA-125) antigen	protein
GMP-EQU-CA-125-Ab	Anti-Equine cancer antigen 125 (CA-125) antibody	antibody(Fc/mFc)
GMP-EQU-CA15-3-Ag	Equine Cancer antigen 15-3 (CA15-3) antigen	protein
GMP-EQU-CA15-3-Ab	Anti-Equine Cancer antigen 15-3 (CA15-3) antibody	antibody(Fc/mFc)
GMP-EQU-CA199-Ag	Equine Carbohydrate antigen 19-9 (CA199) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-CA199-Ab	Anti-Equine Carbohydrate antigen 19-9 (CA199) antibody	antibody(Fc/mFc)
GMP-EQU-CA50-Ag	Equine Carbohydrate Antigen 50 (CA50) antigen	protein
GMP-EQU-CA50-Ab	Anti-Equine Carbohydrate Antigen 50 (CA50) antibody	antibody(Fc/mFc)
GMP-EQU-CA-724-Ag	Equine Cancer Antigen 72-4 (CA-724) antigen	protein
GMP-EQU-CA-724-Ab	Anti-Equine Cancer Antigen 72-4 (CA-724) antibody	antibody(Fc/mFc)
GMP-EQU-Calcitonin-Ag	Equine Calcitonin (Calcitonin) antigen	peptide-(naked/BSA/OVA)

Cat.No	Products Name	Products Type
GMP-EQU-Calcitonin-Ab	Anti-Equine Calcitonin (Calcitonin) antibody	antibody(Fc/mFc)
GMP-EQU-Calmodulin-Ag	Equine Calmodulin antigen	protein
GMP-EQU-Calmodulin-Ab	Anti-Equine Calmodulin antibody	antibody(Fc/mFc)
GMP-EQU-CBL-Ag	Equine cystathionine- β -synthase enzyme (CBL) antigen	protein
GMP-EQU-CBL-Ab	Anti-Equine cystathionine- β -synthase enzyme (CBL) antibody	antibody(Fc/mFc)
GMP-EQU-CBS-Ag	Equine cystathionine- β -synthase enzyme (CBS) antigen	protein
GMP-EQU-CBS-Ab	Anti-Equine cystathionine- β -synthase enzyme (CBS) antibody	antibody(Fc/mFc)
GMP-EQU-CEA-Ag	Equine Carcinoembryonic antigen (CEA) antigen	protein
GMP-EQU-CEA-Ab	Anti-Equine Carcinoembryonic antigen (CEA) antibody	antibody(Fc/mFc)
GMP-EQU-cIV-Ag	Equine collagen type IV (cIV) antigen	protein
GMP-EQU-cIV-Ab	Anti-Equine collagen type IV (cIV) antibody	antibody(Fc/mFc)
GMP-EQU-CK-MB-Ag	Equine Creatine Kinase MB (CK-MB) antigen	protein
GMP-EQU-CK-MB-Ab	Anti-Equine Creatine Kinase MB (CK-MB) antibody	antibody(Fc/mFc)
GMP-EQU-CNP-Ag	Equine C-type natriuretic peptide (CNP) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-CNP-Ab	Anti-Equine C-type natriuretic peptide (CNP) antibody	antibody(Fc/mFc)
GMP-EQU-Copeptin-Ag	Equine Copeptin (Copeptin) antigen	protein
GMP-EQU-Copeptin-Ab	Anti-Equine Copeptin (Copeptin) antibody	antibody(Fc/mFc)
GMP-EQU-C-Peptide-Ag	Equine C-peptide (C-Peptide) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-C-Peptide-Ab	Anti-Equine C-peptide (C-Peptide) antibody	antibody(Fc/mFc)
GMP-EQU-CRP-Ag	Equine CRP (C-Reactive Protein) antigen	protein
GMP-EQU-CRP-Ab	Anti-Equine CRP (C-Reactive Protein) antibody	antibody(Fc/mFc)
GMP-EQU-cTnI-Ag	Equine Cardiac troponin I (cTnI) antigen	protein

Cat.No	Products Name	Products Type
GMP-EQU-cTnI-Ab	Anti-Equine Cardiac troponin I (cTnI) antibody	antibody(Fc/mFc)
GMP-EQU-cTnT-Ag	Equine cardiac troponin T (cTnT) antigen	protein
GMP-EQU-cTnT-Ab	Anti-Equine cardiac troponin T (cTnT) antibody	antibody(Fc/mFc)
GMP-EQU-CYFRA21-1-Ag	Equine soluble fragment of cytokeratin-19 (CYFRA21-1) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-CYFRA21-1-Ab	Anti-Equine soluble fragment of cytokeratin-19 (CYFRA21-1) antibody	antibody(Fc/mFc)
GMP-EQU-CysC-Ag	Equine CysC (cystatin C) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-CysC-Ab	Anti-Equine CysC (cystatin C) antibody	antibody(Fc/mFc)
GMP-EQU-D-Dimer-Ag	Equine D dimer (D-Dimer) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-D-Dimer-Ab	Anti-Equine D dimer (D-Dimer) antibody	antibody(Fc/mFc)
GMP-EQU-Dig-Ag	Equine Digoxigenin (Dig) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-Dig-Ab	Anti-Equine Digoxigenin (Dig) antibody	antibody(Fc/mFc)
GMP-EQU-FER-Ag	Equine Ferritin (FER) antigen	protein
GMP-EQU-FER-Ab	Anti-Equine Ferritin (FER) antibody	antibody(Fc/mFc)
GMP-EQU-fFN-Ag	Equine Fetal fibronectin (fFN) antigen	protein
GMP-EQU-fFN-Ab	Anti-Equine Fetal fibronectin (fFN) antibody	antibody(Fc/mFc)
GMP-EQU-FIB-Ag	Equine Fibrinogen (FIB) antigen	protein
GMP-EQU-FIB-Ab	Anti-Equine Fibrinogen (FIB) antibody	antibody(Fc/mFc)
GMP-EQU-FSH-Ag	Equine follicle stimulating hormone (FSH) antigen	protein
GMP-EQU-FSH-Ab	Anti-Equine follicle stimulating hormone (FSH) antibody	antibody(Fc/mFc)
GMP-EQU-G100-Ag	Equine Gastrin (G100) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-G100-Ab	Anti-Equine Gastrin (G100) antibody	antibody(Fc/mFc)
GMP-EQU-GDF-15-Ag	Equine Growth differentiation factor-15 (GDF-15) antigen	protein

Cat.No	Products Name	Products Type
GMP-EQU-GDF-15-Ab	Anti-Equine Growth differentiation factor-15 (GDF-15) antibody	antibody(Fc/mFc)
GMP-EQU-GFAP-Ag	Equine Glial fibrillary acid protein (GFAP) antigen	protein
GMP-EQU-GFAP-Ab	Anti-Equine Glial fibrillary acid protein (GFAP) antibody	antibody(Fc/mFc)
GMP-EQU-HABP-Ag	Equine Hyaluronic acid binding protein (HABP) antigen	protein
GMP-EQU-HABP-Ab	Anti-Equine Hyaluronic acid binding protein (HABP) antibody	antibody(Fc/mFc)
GMP-EQU-Haptoglobin-Ag	Equine Haptoglobin antigen	protein
GMP-EQU-Haptoglobin-Ab	Anti-Equine Haptoglobin antibody	antibody(Fc/mFc)
GMP-EQU-HBP-Ag	Equine Heparin-binding protein (HBP) antigen	protein
GMP-EQU-HBP-Ab	Anti-Equine Heparin-binding protein (HBP) antibody	antibody(Fc/mFc)
GMP-EQU-HE4-Ag	Equine epididymis protein 4 (HE4) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-HE4-Ab	Anti-Equine epididymis protein 4 (HE4) antibody	antibody(Fc/mFc)
GMP-EQU-hEGF-Ag	Equine epidermal growth factor (hEGF) antigen	protein
GMP-EQU-hEGF-Ab	Anti-Equine epidermal growth factor (hEGF) antibody	antibody(Fc/mFc)
GMP-EQU-hemoglobin-Ag	Equine hemoglobin antigen	protein
GMP-EQU-hemoglobin-Ab	Anti-Equine hemoglobin antibody	antibody(Fc/mFc)
GMP-EQU-HER2-Ag	Equine epidermal growth factor receptor 2 (HER2) antigen	protein
GMP-EQU-HER2-Ab	Anti-Equine epidermal growth factor receptor 2 (HER2) antibody	antibody(Fc/mFc)
GMP-EQU-H-FABP-Ag	Equine Heart-type Fatty Acid-Binding Protein (H-FABP) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-H-FABP-Ab	Anti-Equine Heart-type Fatty Acid-Binding Protein (H-FABP) antibody	antibody(Fc/mFc)
GMP-EQU-HSP90a-Ag	Equine Heat shock proteins 90-alpha (HSP90α) antigen	protein
GMP-EQU-HSP90a-Ab	Anti-Equine Heat shock proteins 90-alpha (HSP90α) antibody	antibody(Fc/mFc)
GMP-EQU-IFN-γ-Ag	Equine IFN-γ (Interferon gamma) antigen	protein

Cat.No	Products Name	Products Type
GMP-EQU-IFN-γ-Ab	Anti-Equine IFN-γ (Interferon gamma) antibody	antibody(Fc/mFc)
GMP-EQU-IgG-Ag	Equine IgG antigen	protein
GMP-EQU-IgG-Ab	Anti-Equine IgG antibody	antibody(Fc/mFc)
GMP-EQU-IL-6-Ag	Equine Interleukin 6 (IL-6) antigen	protein
GMP-EQU-IL-6-Ab	Anti-Equine Interleukin 6 (IL-6) antibody	antibody(Fc/mFc)
GMP-EQU-Ins-Ag	Equine Insulin (Ins) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-Ins-Ab	Anti-Equine Insulin (Ins) antibody	antibody(Fc/mFc)
GMP-EQU-LN-Ag	Equine laminins (LN) antigen	protein
GMP-EQU-LN-Ab	Anti-Equine laminins (LN) antibody	antibody(Fc/mFc)
GMP-EQU-MCP-1-Ag	Equine Monocyte chemoattractant protein-1 (MCP-1) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-MCP-1-Ab	Anti-Equine Monocyte chemoattractant protein-1 (MCP-1) antibody	antibody(Fc/mFc)
GMP-EQU-MMP3-Ag	Equine matrix metalloproteinase-3 (MMP3) antigen	protein
GMP-EQU-MMP3-Ab	Anti-Equine matrix metalloproteinase-3 (MMP3) antibody	antibody(Fc/mFc)
GMP-EQU-MOG-Ag	Equine myelin oligodendrocyte glycoprotein (MOG) antigen	protein
GMP-EQU-MOG-Ab	Anti-Equine myelin oligodendrocyte glycoprotein (MOG) antibody	antibody(Fc/mFc)
GMP-EQU-MPO-Ag	Equine Myeloperoxidase (MPO) antigen	protein
GMP-EQU-MPO-Ab	Anti-Equine Myeloperoxidase (MPO) antibody	antibody(Fc/mFc)
GMP-EQU-Myo-Ag	Equine Myoglobin (Myo) antigen	protein
GMP-EQU-Myo-Ab	Anti-Equine Myoglobin (Myo) antibody	antibody(Fc/mFc)
GMP-EQU-neurogranin-Ag	Equine neurogranin protein (neurogranin) antigen	protein
GMP-EQU-neurogranin-Ab	Anti-Equine neurogranin protein (neurogranin) antibody	antibody(Fc/mFc)
GMP-EQU-NFL-Ag	Equine neurofilament light chain (NFL) antigen	protein

Cat.No	Products Name	Products Type
GMP-EQU-NFL-Ab	Anti-Equine neurofilament light chain (NFL) antibody	antibody(Fc/mFc)
GMP-EQU-NGAL-Ag	Equine Neutrophil gelatinase-associated lipocalin (NGAL) antigen	protein
GMP-EQU-NGAL-Ab	Anti-Equine Neutrophil gelatinase-associated lipocalin (NGAL) antibody	antibody(Fc/mFc)
GMP-EQU-N-MID-OC-Ag	Equine N-terminal midfragment of Osteocalcin (N-MID OC) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-N-MID-OC-Ab	Anti-Equine N-terminal midfragment of Osteocalcin (N-MID OC) antibody	antibody(Fc/mFc)
GMP-EQU-NSE-Ag	Equine Neuron specific enolase (NSE) antigen	protein
GMP-EQU-NSE-Ab	Anti-Equine Neuron specific enolase (NSE) antibody	antibody(Fc/mFc)
GMP-EQU-NT-proBNP-Ag	Equine NT-proBNP (N-terminal prohormone of brain natriuretic peptide) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-NT-proBNP-Ab	Anti-Equine NT-proBNP (N-terminal prohormone of brain natriuretic peptide) antibody	antibody(Fc/mFc)
GMP-EQU-PA-Ag	Equine Prealbumin (PA) antigen	protein
GMP-EQU-PA-Ab	Anti-Equine Prealbumin (PA) antibody	antibody(Fc/mFc)
GMP-EQU-Pancreatic-elastase-Ag	Equine PE-1 (Pancreatic elastase) antigen	protein
GMP-EQU-Pancreatic-elastase-Ab	Anti-Equine PE-1 (Pancreatic elastase) antibody	antibody(Fc/mFc)
GMP-EQU-pancreatic-lipase-Ag	Equine pancreatic lipase antigen	protein
GMP-EQU-pancreatic-lipase-Ab	Anti-Equine pancreatic lipase antibody	antibody(Fc/mFc)
GMP-EQU-PCT-Ag	Equine procalcitonin (PCT) antigen	protein
GMP-EQU-PCT-Ab	Anti-Equine procalcitonin (PCT) antibody	antibody(Fc/mFc)
GMP-EQU-PG-I-Ag	Equine pepsinogen I (PG I) antigen	protein
GMP-EQU-PG-I-Ab	Anti-Equine pepsinogen I (PG I) antibody	antibody(Fc/mFc)
GMP-EQU-PG-II-Ag	Equine Pepsinogen II (PG II) antigen	protein
GMP-EQU-PG-II-Ab	Anti-Equine Pepsinogen II (PG II) antibody	antibody(Fc/mFc)
GMP-EQU-PIIINP-Ag	Equine Procollagen III N-terminal peptide (PIIINP) antigen	peptide-(naked/BSA/OVA)

Cat.No	Products Name	Products Type
GMP-EQU-PIIINP-Ab	Anti-Equine Procollagen III N-terminal peptide (PIIINP) antibody	antibody(Fc/mFc)
GMP-EQU-PINP-Ag	Equine amino-terminal propeptide (PINP) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-PINP-Ab	Anti-Equine amino-terminal propeptide (PINP) antibody	antibody(Fc/mFc)
GMP-EQU-PIVKA-II-Ag	Equine Protein induced by vitamin K absence-II (PIVKA-II) antigen	protein
GMP-EQU-PIVKA-II-Ab	Anti-Equine Protein induced by vitamin K absence-II (PIVKA-II) antibody	antibody(Fc/mFc)
GMP-EQU-PLGF-1-Ag	Equine placental growth factor-1 (PLGF-1) antigen	protein
GMP-EQU-PLGF-1-Ab	Anti-Equine placental growth factor-1 (PLGF-1) antibody	antibody(Fc/mFc)
GMP-EQU-Pre-Spec-Pro-B-Ag	Equine Pregnancy Specific Protein B antigen	protein
GMP-EQU-Pre-Spec-Pro-B-Ab	Anti-Equine Pregnancy Specific Protein B antibody	antibody(Fc/mFc)
GMP-EQU-Preg-ass-gly-1-Ag	Equine Pregnancy-associated glycoprotein 1 antigen	protein
GMP-EQU-Preg-ass-gly-1-Ab	Anti-Equine Pregnancy-associated glycoprotein 1 antibody	antibody(Fc/mFc)
GMP-EQU-proBNP-Ag	Equine pro B-type natriuretic peptide (proBNP) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-proBNP-Ab	Anti-Equine pro B-type natriuretic peptide (proBNP) antibody	antibody(Fc/mFc)
GMP-EQU-Prog-Ag	Equine Prog (Progesterone) antigen	protein
GMP-EQU-Prog-Ab	Anti-Equine Prog (Progesterone) antibody	antibody(Fc/mFc)
GMP-EQU-ProGRP-Ag	Equine Gastrin-releasing peptide (ProGRP) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-ProGRP-Ab	Anti-Equine Gastrin-releasing peptide (ProGRP) antibody	antibody(Fc/mFc)
GMP-EQU-Proinsulin-Ag	Equine Proinsulin (Proinsulin) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-Proinsulin-Ab	Anti-Equine Proinsulin (Proinsulin) antibody	antibody(Fc/mFc)
GMP-EQU-PSA-Ag	Equine Prostate-specific antigen (PSA) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-PSA-Ab	Anti-Equine Prostate-specific antigen (PSA) antibody	antibody(Fc/mFc)
GMP-EQU-p-tau181-Ag	Equine Phospho-tau181 (p-tau181) antigen	protein

Cat.No	Products Name	Products Type
GMP-EQU-p-tau181-Ab	Anti-Equine Phospho-tau181 (p-tau181) antibody	antibody(Fc/mFc)
GMP-EQU-p-tau202-205-Ag	Equine Phospho-tau202 and 205 (p-tau202 and 205) antigen	protein
GMP-EQU-p-tau202-205-Ab	Anti-Equine Phospho-tau202 and 205 (p-tau202 and 205) antibody	antibody(Fc/mFc)
GMP-EQU-p-tau212-214-Ag	Equine Phospho-tau212 and 214 (p-tau212 and 214) antigen	protein
GMP-EQU-p-tau212-214-Ab	Anti-Equine Phospho-tau212 and 214 (p-tau212 and 214) antibody	antibody(Fc/mFc)
GMP-EQU-p-tau217-Ag	Equine Phospho-tau217 (p-tau217) antigen	protein
GMP-EQU-p-tau217-Ab	Anti-Equine Phospho-tau217 (p-tau217) antibody	antibody(Fc/mFc)
GMP-EQU-p-tau231-Ag	Equine Phospho-tau231 (p-tau231) antigen	protein
GMP-EQU-p-tau231-Ab	Anti-Equine Phospho-tau231 (p-tau231) antibody	antibody(Fc/mFc)
GMP-EQU-p-tau413-Ag	Equine Phospho-tau413 (p-tau413) antigen	protein
GMP-EQU-p-tau413-Ab	Anti-Equine Phospho-tau413 (p-tau413) antibody	antibody(Fc/mFc)
GMP-EQU-p-tau422-Ag	Equine Phospho-tau422 (p-tau422) antigen	protein
GMP-EQU-p-tau422-Ab	Anti-Equine Phospho-tau422 (p-tau422) antibody	antibody(Fc/mFc)
GMP-EQU-PTH-Ag	Equine parathyroid hormone (PTH) antigen	protein
GMP-EQU-PTH-Ab	Anti-Equine parathyroid hormone (PTH) antibody	antibody(Fc/mFc)
GMP-EQU-relaxin-Ag	Equine RLN (Relaxin) antigen	protein
GMP-EQU-relaxin-Ab	Anti-Equine RLN (Relaxin) antibody	antibody(Fc/mFc)
GMP-EQU-S100-A8-Ag	Equine Calprotectin (S100-A8) antigen	protein
GMP-EQU-S100-A8-Ab	Anti-Equine Calprotectin (S100-A8) antibody	antibody(Fc/mFc)
GMP-EQU-S100-A9-Ag	Equine Calprotectin heavy chain (S100-A9) antigen	protein
GMP-EQU-S100-A9-Ab	Anti-Equine Calprotectin heavy chain (S100-A9) antibody	antibody(Fc/mFc)
GMP-EQU-S100B-Ag	Equine S100 proteins (S100B) antigen	protein

Cat.No	Products Name	Products Type
GMP-EQU-S100B-Ab	Anti-Equine S100 proteins (S100B) antibody	antibody(Fc/mFc)
GMP-EQU-SAHH-Ag	Equine S-adenosylhomocysteine (AdoHcy) hydrolase (SAHH) antigen	protein
GMP-EQU-SAHH-Ab	Anti-Equine S-adenosylhomocysteine (AdoHcy) hydrolase (SAHH) antibody	antibody(Fc/mFc)
GMP-EQU-SCC-Ag-Ag	Equine squamous cell carcinoma antigen (SCC Ag) antigen	protein
GMP-EQU-SCC-Ag-Ab	Anti-Equine squamous cell carcinoma antigen (SCC Ag) antibody	antibody(Fc/mFc)
GMP-EQU-Serum-amyloid-A-Ag	Equine SAA (Serum amyloid A) antigen	protein
GMP-EQU-Serum-amyloid-A-Ab	Anti-Equine SAA (Serum amyloid A) antibody	antibody(Fc/mFc)
GMP-EQU-sFlt-1-Ag	Equine soluble vascular endothelial growth factor receptor-2 (sFlt-1) antigen	protein
GMP-EQU-sFlt-1-Ab	Anti-Equine soluble vascular endothelial growth factor receptor-2 (sFlt-1) antibody	antibody(Fc/mFc)
GMP-EQU-SPA-Ag	Equine Staphylococcus protein A (SPA) antigen	protein
GMP-EQU-SPA-Ab	Anti-Equine Staphylococcus protein A (SPA) antibody	antibody(Fc/mFc)
GMP-EQU-Spirochaeta-Ag	Equine Spirochaeta antigen	protein
GMP-EQU-Spirochaeta-Ab	Anti-Equine Spirochaeta antibody	antibody(Fc/mFc)
GMP-EQU-sST2-Ag	Equine soluble ST2 (sST2) antigen	protein
GMP-EQU-sST2-Ab	Anti-Equine soluble ST2 (sST2) antibody	antibody(Fc/mFc)
GMP-EQU-sTREM2-Ag	Equine soluble fragment of triggering receptor expressed on myeloid cells 2 (sTREM2) antigen	protein
GMP-EQU-sTREM2-Ab	Anti-Equine soluble fragment of triggering receptor expressed on myeloid cells 2 (sTREM2) antibody	antibody(Fc/mFc)
GMP-EQU-Tau-Ag	Equine tau proteins (Tau) antigen	protein
GMP-EQU-Tau-Ab	Anti-Equine tau proteins (Tau) antibody	antibody(Fc/mFc)
GMP-EQU-TG-Ag	Equine Thyroglobulin (TG) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-TG-Ab	Anti-Equine Thyroglobulin (TG) antibody	antibody(Fc/mFc)
GMP-EQU-TIMP-1-Ag	Equine Matrix Metalloproteinase Inhibitor-1 (TIMP-1) antigen	protein

Cat.No	Products Name	Products Type
GMP-EQU-TIMP-1-Ab	Anti-Equine Matrix Metalloproteinase Inhibitor-1 (TIMP-1) antibody	antibody(Fc/mFc)
GMP-EQU-TK1-Ag	Equine Thymidine Kinase 1 (TK1) antigen	protein
GMP-EQU-TK1-Ab	Anti-Equine Thymidine Kinase 1 (TK1) antibody	antibody(Fc/mFc)
GMP-EQU-TNFa-Ag	Equine tumor necrosis factor alpha (TNFα) antigen	protein
GMP-EQU-TNFa-Ab	Anti-Equine tumor necrosis factor alpha (TNFα) antibody	antibody(Fc/mFc)
GMP-EQU-TPO-Ag	Equine Thyroid peroxidase (TPO) antigen	protein
GMP-EQU-TPO-Ab	Anti-Equine Thyroid peroxidase (TPO) antibody	antibody(Fc/mFc)
GMP-EQU-TrxR-Ag	Equine thioredoxin reductase (TrxR) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-TrxR-Ab	Anti-Equine thioredoxin reductase (TrxR) antibody	antibody(Fc/mFc)
GMP-EQU-TSH-Ag	Equine Thyroid-stimulating hormone (TSH) antigen	protein
GMP-EQU-TSH-Ab	Anti-Equine Thyroid-stimulating hormone (TSH) antibody	antibody(Fc/mFc)
GMP-EQU-VDBP-Ag	Equine vitamin D binding protein (VDBP) antigen	protein
GMP-EQU-VDBP-Ab	Anti-Equine vitamin D binding protein (VDBP) antibody	antibody(Fc/mFc)
GMP-EQU-VEGF-Ag	Equine Vascular Endothelial Growth Factor (VEGF) antigen	protein
GMP-EQU-VEGF-Ab	Anti-Equine Vascular Endothelial Growth Factor (VEGF) antibody	antibody(Fc/mFc)
GMP-EQU-VLP-1-Ag	Equine Visinin-like protein 1 (VLP-1) antigen	protein
GMP-EQU-VLP-1-Ab	Anti-Equine Visinin-like protein 1 (VLP-1) antibody	antibody(Fc/mFc)
GMP-EQU-YKL-40-Ag	Equine Chitinase-3-like protein 1 (YKL-40) antigen	protein
GMP-EQU-YKL-40-Ab	Anti-Equine Chitinase-3-like protein 1 (YKL-40) antibody	antibody(Fc/mFc)
GMP-EQU-a-defensin-1-Ag	Equine alpha defensin 1 (α-defensin 1) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-a-defensin-1-Ab	Anti-Equine alpha defensin 1 (α-defensin 1) antibody	antibody(Fc/mFc)
GMP-EQU-a-defensin-5-Ag	Equine alpha defensin 5 (α-defensin 5) antigen	peptide-(naked/BSA/OVA)

Cat.No	Products Name	Products Type
GMP-EQU-a-defensin-5-Ab	Anti-Equine alpha defensin 5 (α -defensin 5) antibody	antibody(Fc/mFc)
GMP-EQU-a-synuclein-Ag	Equine Alpha-synuclein (α -synuclein) antigen	protein
GMP-EQU-a-synuclein-Ab	Anti-Equine Alpha-synuclein (α -synuclein) antibody	antibody(Fc/mFc)
GMP-EQU-B-CTx-Ag	Equine β -isomerized C-terminal telopeptides (β -CTx) antigen	peptide-(naked/BSA/OVA)
GMP-EQU-B-CTx-Ab	Anti-Equine β -isomerized C-terminal telopeptides (β -CTx) antibody	antibody(Fc/mFc)

Avian

Cat.No	Products Name	Products Type
GMP-VT-P046-Tg003-Ag	Chlamydiaceae Outer Membrane Protein A antigen	protein
GMP-VT-P046-Tg003-Ab	Anti-Chlamydiaceae Outer Membrane Protein A antibody	antibody(Fc/mFc)
GMP-VT-P101-Tg002-Ag	Pasteurella multocida Outer Membrane Protein A antigen	protein
GMP-VT-P101-Tg002-Ab	Anti-Pasteurella multocida Outer Membrane Protein A antibody	antibody(Fc/mFc)
GMP-VT-P101-Tg003-Ag	Pasteurella multocida Outer Membrane Protein H antigen	protein
GMP-VT-P101-Tg003-Ab	Anti-Pasteurella multocida Outer Membrane Protein H antibody	antibody(Fc/mFc)
GMP-VT-P149-Tg001-Ag	Treponema antigen	protein
GMP-VT-P149-Tg001-Ab	Anti-Treponema antibody	antibody(Fc/mFc)
GMP-VT-P150-Tg001-Ag	IBV - Variant ItalyO2 Spike Protein S1 antigen	protein
GMP-VT-P150-Tg001-Ab	Anti-IBV - Variant ItalyO2 Spike Protein S1 antibody	antibody(Fc/mFc)
GMP-VT-P151-Tg001-Ag	IBV - Variant Arkansas Spike Protein S1 Gene antigen	protein
GMP-VT-P151-Tg001-Ab	Anti-IBV - Variant Arkansas Spike Protein S1 Gene antibody	antibody(Fc/mFc)
GMP-VT-P152-Tg001-Ag	IBV - Variant D274 Spike Protein S1 Gene antigen	protein
GMP-VT-P152-Tg001-Ab	Anti-IBV - Variant D274 Spike Protein S1 Gene antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P153-Tg001-Ag	IBV - Variant IB80 Spike Protein S Gene antigen	protein
GMP-VT-P153-Tg001-Ab	Anti-IBV - Variant IB80 Spike Protein S Gene antibody	antibody(Fc/mFc)
GMP-VT-P154-Tg001-Ag	IBV - Variant Massachusetts Spike Protein S1 antigen	protein
GMP-VT-P154-Tg001-Ab	Anti-IBV - Variant Massachusetts Spike Protein S1 antibody	antibody(Fc/mFc)
GMP-VT-P155-Tg001-Ag	Salmonella Gallinarum antigen	protein
GMP-VT-P155-Tg001-Ab	Anti-Salmonella Gallinarum antibody	antibody(Fc/mFc)
GMP-VT-P156-Tg001-Ag	Histomonas meleagridis antigen	protein
GMP-VT-P156-Tg001-Ab	Anti-Histomonas meleagridis antibody	antibody(Fc/mFc)
GMP-VT-P157-Tg001-Ag	Mycoplasma avium antigen	protein
GMP-VT-P157-Tg001-Ab	Anti-Mycoplasma avium antibody	antibody(Fc/mFc)
GMP-VT-P158-Tg001-Ag	Duck Baiyangdian Virus or Duck tembusu virus Envelope Protein antigen	protein
GMP-VT-P158-Tg001-Ab	Anti-Duck Baiyangdian Virus or Duck tembusu virus Envelope Protein antibody	antibody(Fc/mFc)
GMP-VT-P159-Tg001-Ag	Salmonella gallinarum , Salmonella pullorum Invasion Plasmid Antigen J antigen	protein
GMP-VT-P159-Tg001-Ab	Anti-Salmonella gallinarum , Salmonella pullorum Invasion Plasmid Antigen J antibody	antibody(Fc/mFc)
GMP-VT-P159-Tg002-Ag	Salmonella gallinarum , Salmonella pullorum SPIC Protein antigen	protein
GMP-VT-P159-Tg002-Ab	Anti-Salmonella gallinarum , Salmonella pullorum SPIC Protein antibody	antibody(Fc/mFc)
GMP-VT-P159-Tg003-Ag	Salmonella gallinarum , Salmonella pullorum Fimbrial Regulatory Protein antigen	protein
GMP-VT-P159-Tg003-Ab	Anti-Salmonella gallinarum , Salmonella pullorum Fimbrial Regulatory Protein antibody	antibody(Fc/mFc)
GMP-VT-P160-Tg001-Ag	ARV-A, ARV-D Viral Protein 6 antigen	protein
GMP-VT-P160-Tg001-Ab	Anti-ARV-A, ARV-D Viral Protein 6 antibody	antibody(Fc/mFc)
GMP-VT-P161-Tg001-Ag	IBV - Variant QX antigen	protein
GMP-VT-P161-Tg001-Ab	Anti-IBV - Variant QX antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P162-Tg001-Ag	IBV - Variant Q1 antigen	protein
GMP-VT-P162-Tg001-Ab	Anti-IBV - Variant Q1 antibody	antibody(Fc/mFc)
GMP-VT-P163-Tg001-Ag	GPV, MDPV Non-structural Protein 1 antigen	protein
GMP-VT-P163-Tg001-Ab	Anti-GPV, MDPV Non-structural Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P163-Tg002-Ag	GPV, MDPV Viral Protein 3 antigen	protein
GMP-VT-P163-Tg002-Ab	Anti-GPV, MDPV Viral Protein 3 antibody	antibody(Fc/mFc)
GMP-VT-P163-Tg003-Ag	GPV, MDPV Viral Protein 2 antigen	protein
GMP-VT-P163-Tg003-Ab	Anti-GPV, MDPV Viral Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P164-Tg001-Ag	APV Nucleoprotein antigen	protein
GMP-VT-P164-Tg001-Ab	Anti-APV Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P165-Tg001-Ag	Chlamydia avium Major Outer Membrane Protein antigen	protein
GMP-VT-P165-Tg001-Ab	Anti-Chlamydia avium Major Outer Membrane Protein antibody	antibody(Fc/mFc)
GMP-VT-P166-Tg001-Ag	Turkey Hepatitis Virus antigen	protein
GMP-VT-P166-Tg001-Ab	Anti-Turkey Hepatitis Virus antibody	antibody(Fc/mFc)
GMP-VT-P167-Tg001-Ag	AIV-A H5N1 subtype Nucleoprotein antigen	protein
GMP-VT-P167-Tg001-Ab	Anti-AIV-A H5N1 subtype Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P167-Tg002-Ag	AIV-A H5N1 subtype Hemagglutinin antigen	protein
GMP-VT-P167-Tg002-Ab	Anti-AIV-A H5N1 subtype Hemagglutinin antibody	antibody(Fc/mFc)
GMP-VT-P168-Tg001-Ag	NDV Nucleoprotein antigen	protein
GMP-VT-P168-Tg001-Ab	Anti-NDV Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P169-Tg001-Ag	IBV Nucleoprotein antigen	protein
GMP-VT-P169-Tg001-Ab	Anti-IBV Nucleoprotein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P170-Tg001-Ag	FAdV Fiber Protein antigen	protein
GMP-VT-P170-Tg001-Ab	Anti-FAdV Fiber Protein antibody	antibody(Fc/mFc)
GMP-VT-P171-Tg001-Ag	IBDV Viral Protein 3 antigen	protein
GMP-VT-P171-Tg001-Ab	Anti-IBDV Viral Protein 3 antibody	antibody(Fc/mFc)
GMP-VT-P171-Tg002-Ag	IBDV Viral Protein 4 antigen	protein
GMP-VT-P171-Tg002-Ab	Anti-IBDV Viral Protein 4 antibody	antibody(Fc/mFc)
GMP-VT-P171-Tg003-Ag	IBDV Viral Protein X antigen	protein
GMP-VT-P171-Tg003-Ab	Anti-IBDV Viral Protein X antibody	antibody(Fc/mFc)
GMP-VT-P171-Tg004-Ag	IBDV Viral Protein 2 antigen	protein
GMP-VT-P171-Tg004-Ab	Anti-IBDV Viral Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P172-Tg001-Ag	ALV Protein 27 antigen	protein
GMP-VT-P172-Tg001-Ab	Anti-ALV Protein 27 antibody	antibody(Fc/mFc)
GMP-VT-P173-Tg001-Ag	MDV Glycoprotein E antigen	protein
GMP-VT-P173-Tg001-Ab	Anti-MDV Glycoprotein E antibody	antibody(Fc/mFc)
GMP-VT-P173-Tg002-Ag	MDV Glycoprotein I antigen	protein
GMP-VT-P173-Tg002-Ab	Anti-MDV Glycoprotein I antibody	antibody(Fc/mFc)
GMP-VT-P174-Tg001-Ag	CAstV Capsid Protein 11672 antigen	protein
GMP-VT-P174-Tg001-Ab	Anti-CAstV Capsid Protein 11672 antibody	antibody(Fc/mFc)
GMP-VT-P175-Tg001-Ag	AIV Nucleoprotein antigen	protein
GMP-VT-P175-Tg001-Ab	Anti-AIV Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P175-Tg002-Ag	AIV Hemagglutinin antigen	protein
GMP-VT-P175-Tg002-Ab	Anti-AIV Hemagglutinin antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P176-Tg001-Ag	DHV Viral Protein 3 antigen	protein
GMP-VT-P176-Tg001-Ab	Anti-DHV Viral Protein 3 antibody	antibody(Fc/mFc)
GMP-VT-P176-Tg002-Ag	DHV Viral Protein 1 antigen	protein
GMP-VT-P176-Tg002-Ab	Anti-DHV Viral Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P177-Tg001-Ag	Mycoplasma gallisepticum MGC2 Protein antigen	protein
GMP-VT-P177-Tg001-Ab	Anti-Mycoplasma gallisepticum MGC2 Protein antibody	antibody(Fc/mFc)
GMP-VT-P177-Tg002-Ag	Mycoplasma gallisepticum MGC3 Protein antigen	protein
GMP-VT-P177-Tg002-Ab	Anti-Mycoplasma gallisepticum MGC3 Protein antibody	antibody(Fc/mFc)
GMP-VT-P177-Tg003-Ag	Mycoplasma gallisepticum Hemagglutinin Protein antigen	protein
GMP-VT-P177-Tg003-Ab	Anti-Mycoplasma gallisepticum Hemagglutinin Protein antibody	antibody(Fc/mFc)
GMP-VT-P177-Tg004-Ag	Mycoplasma gallisepticum MG Protein antigen	protein
GMP-VT-P177-Tg004-Ab	Anti-Mycoplasma gallisepticum MG Protein antibody	antibody(Fc/mFc)
GMP-VT-P178-Tg001-Ag	Mycoplasma synoviae MSPB Protein antigen	protein
GMP-VT-P178-Tg001-Ab	Anti-Mycoplasma synoviae MSPB Protein antibody	antibody(Fc/mFc)
GMP-VT-P178-Tg002-Ag	Mycoplasma synoviae MS Protein antigen	protein
GMP-VT-P178-Tg002-Ab	Anti-Mycoplasma synoviae MS Protein antibody	antibody(Fc/mFc)
GMP-VT-P179-Tg001-Ag	ILTV Glycoprotein B antigen	protein
GMP-VT-P179-Tg001-Ab	Anti-ILTV Glycoprotein B antibody	antibody(Fc/mFc)
GMP-VT-P180-Tg001-Ag	Mycoplasma meleagridis MM Protein antigen	protein
GMP-VT-P180-Tg001-Ab	Anti-Mycoplasma meleagridis MM Protein antibody	antibody(Fc/mFc)
GMP-VT-P181-Tg001-Ag	Ornithobacterium rhinotracheale OR77 Protein antigen	protein
GMP-VT-P181-Tg001-Ab	Anti-Ornithobacterium rhinotracheale OR77 Protein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P182-Tg001-Ag	ARV Structural Protein Sigma B antigen	protein
GMP-VT-P182-Tg001-Ab	Anti-ARV Structural Protein Sigma B antibody	antibody(Fc/mFc)
GMP-VT-P182-Tg002-Ag	ARV Structural Protein Sigma C antigen	protein
GMP-VT-P182-Tg002-Ab	Anti-ARV Structural Protein Sigma C antibody	antibody(Fc/mFc)
GMP-VT-P183-Tg001-Ag	CAV Viral Protein 1 antigen	protein
GMP-VT-P183-Tg001-Ab	Anti-CAV Viral Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P183-Tg002-Ag	CAV Viral Protein 2 antigen	protein
GMP-VT-P183-Tg002-Ab	Anti-CAV Viral Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P183-Tg003-Ag	CAV Viral Protein 3 antigen	protein
GMP-VT-P183-Tg003-Ab	Anti-CAV Viral Protein 3 antibody	antibody(Fc/mFc)
GMP-VT-P184-Tg001-Ag	Chlamydia psittaci Major Outer Membrane Protein antigen	protein
GMP-VT-P184-Tg001-Ab	Anti-Chlamydia psittaci Major Outer Membrane Protein antibody	antibody(Fc/mFc)
GMP-VT-P185-Tg001-Ag	AIV-A H1N1 subtype Nucleoprotein antigen	protein
GMP-VT-P185-Tg001-Ab	Anti-AIV-A H1N1 subtype Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P185-Tg002-Ag	AIV-A H1N1 subtype Hemagglutinin antigen	protein
GMP-VT-P185-Tg002-Ab	Anti-AIV-A H1N1 subtype Hemagglutinin antibody	antibody(Fc/mFc)
GMP-VT-P186-Tg001-Ag	AIV-A H3N2 subtype Nucleoprotein antigen	protein
GMP-VT-P186-Tg001-Ab	Anti-AIV-A H3N2 subtype Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P186-Tg002-Ag	AIV-A H3N2 subtype Hemagglutinin antigen	protein
GMP-VT-P186-Tg002-Ab	Anti-AIV-A H3N2 subtype Hemagglutinin antibody	antibody(Fc/mFc)
GMP-VT-P187-Tg001-Ag	AIV-A H7N2 subtype Nucleoprotein antigen	protein
GMP-VT-P187-Tg001-Ab	Anti-AIV-A H7N2 subtype Nucleoprotein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P187-Tg002-Ag	AIV-A H7N2 subtype Hemagglutinin antigen	protein
GMP-VT-P187-Tg002-Ab	Anti-AIV-A H7N2 subtype Hemagglutinin antibody	antibody(Fc/mFc)
GMP-VT-P188-Tg001-Ag	AIV-A H7N7 subtype Nucleoprotein antigen	protein
GMP-VT-P188-Tg001-Ab	Anti-AIV-A H7N7 subtype Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P188-Tg002-Ag	AIV-A H7N7 subtype Hemagglutinin antigen	protein
GMP-VT-P188-Tg002-Ab	Anti-AIV-A H7N7 subtype Hemagglutinin antibody	antibody(Fc/mFc)
GMP-VT-P189-Tg001-Ag	Mycoplasma iowae P45 Protein antigen	protein
GMP-VT-P189-Tg001-Ab	Anti-Mycoplasma iowae P45 Protein antibody	antibody(Fc/mFc)
GMP-VT-P190-Tg001-Ag	REV Glycoprotein 90 antigen	protein
GMP-VT-P190-Tg001-Ab	Anti-REV Glycoprotein 90 antibody	antibody(Fc/mFc)
GMP-VT-P191-Tg001-Ag	Salmonella Enteritidis Flagellin antigen	protein
GMP-VT-P191-Tg001-Ab	Anti-Salmonella Enteritidis Flagellin antibody	antibody(Fc/mFc)
GMP-VT-P192-Tg001-Ag	Avian Coronavirus Nucleoprotein antigen	protein
GMP-VT-P192-Tg001-Ab	Anti-Avian Coronavirus Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P192-Tg002-Ag	Avian Coronavirus Spike Glycoprotein S1 antigen	protein
GMP-VT-P192-Tg002-Ab	Anti-Avian Coronavirus Spike Glycoprotein S1 antibody	antibody(Fc/mFc)
GMP-VT-P193-Tg001-Ag	AEV Viral Protein 1 antigen	protein
GMP-VT-P193-Tg001-Ab	Anti-AEV Viral Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P193-Tg002-Ag	AEV Viral Protein 0 antigen	protein
GMP-VT-P193-Tg002-Ab	Anti-AEV Viral Protein 0 antibody	antibody(Fc/mFc)
GMP-VT-P193-Tg003-Ag	AEV Viral Protein 3 antigen	protein
GMP-VT-P193-Tg003-Ab	Anti-AEV Viral Protein 3 antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P194-Tg001-Ag	Avian Hepatitis E Virus Capsid Protein antigen	protein
GMP-VT-P194-Tg001-Ab	Anti-Avian Hepatitis E Virus Capsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P195-Tg001-Ag	aMPV Nucleoprotein antigen	protein
GMP-VT-P195-Tg001-Ab	Anti-aMPV Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P196-Tg001-Ag	aMPV-A & aMPV-B Nucleoprotein antigen	protein
GMP-VT-P196-Tg001-Ab	Anti-aMPV-A & aMPV-B Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P197-Tg001-Ag	aMPV-C Nucleoprotein antigen	protein
GMP-VT-P197-Tg001-Ab	Anti-aMPV-C Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P198-Tg001-Ag	aMPV-D Nucleoprotein antigen	protein
GMP-VT-P198-Tg001-Ab	Anti-aMPV-D Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P199-Tg001-Ag	ANV Capsid Protein antigen	protein
GMP-VT-P199-Tg001-Ab	Anti-ANV Capsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P200-Tg001-Ag	Avibacterium paragallinarum Haemagglutinin antigen	protein
GMP-VT-P200-Tg001-Ab	Anti-Avibacterium paragallinarum Haemagglutinin antibody	antibody(Fc/mFc)
GMP-VT-P201-Tg001-Ag	Bordetella avium Outer Membrane Protein A antigen	protein
GMP-VT-P201-Tg001-Ab	Anti-Bordetella avium Outer Membrane Protein A antibody	antibody(Fc/mFc)
GMP-VT-P202-Tg001-Ag	Campylobacter Protein 39 antigen	protein
GMP-VT-P202-Tg001-Ab	Anti-Campylobacter Protein 39 antibody	antibody(Fc/mFc)
GMP-VT-P202-Tg002-Ag	Campylobacter Protein 18 antigen	protein
GMP-VT-P202-Tg002-Ab	Anti-Campylobacter Protein 18 antibody	antibody(Fc/mFc)
GMP-VT-P203-Tg001-Ag	Duck Circovirus Capsid Protein antigen	protein
GMP-VT-P203-Tg001-Ab	Anti-Duck Circovirus Capsid Protein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P204-Tg001-Ag	DEV UL24 Antigen antigen	protein
GMP-VT-P204-Tg001-Ab	Anti-DEV UL24 Antigen antibody	antibody(Fc/mFc)
GMP-VT-P204-Tg002-Ag	DEV Glycoprotein B antigen	protein
GMP-VT-P204-Tg002-Ab	Anti-DEV Glycoprotein B antibody	antibody(Fc/mFc)
GMP-VT-P204-Tg003-Ag	DEV Glycoprotein C antigen	protein
GMP-VT-P204-Tg003-Ab	Anti-DEV Glycoprotein C antibody	antibody(Fc/mFc)
GMP-VT-P205-Tg001-Ag	DHV-1 Viral Protein 1 antigen	protein
GMP-VT-P205-Tg001-Ab	Anti-DHV-1 Viral Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P205-Tg002-Ag	DHV-1 Viral Protein 3 antigen	protein
GMP-VT-P205-Tg002-Ab	Anti-DHV-1 Viral Protein 3 antibody	antibody(Fc/mFc)
GMP-VT-P206-Tg001-Ag	EDSV Knob-S Protein antigen	protein
GMP-VT-P206-Tg001-Ab	Anti-EDSV Knob-S Protein antibody	antibody(Fc/mFc)
GMP-VT-P207-Tg001-Ag	Erysipelothrix spp. Surface Protective Antigen A antigen	protein
GMP-VT-P207-Tg001-Ab	Anti-Erysipelothrix spp. Surface Protective Antigen A antibody	antibody(Fc/mFc)
GMP-VT-P208-Tg001-Ag	FPV FPV140 Protein antigen	protein
GMP-VT-P208-Tg001-Ab	Anti-FPV FPV140 Protein antibody	antibody(Fc/mFc)
GMP-VT-P209-Tg001-Ag	Goose Circovirus Capsid Protein antigen	protein
GMP-VT-P209-Tg001-Ab	Anti-Goose Circovirus Capsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P210-Tg001-Ag	GHPV Viral Protein 1 antigen	protein
GMP-VT-P210-Tg001-Ab	Anti-GHPV Viral Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P211-Tg001-Ag	IBV - Variant 02 Spike Glycoprotein S2 Subunit antigen	protein
GMP-VT-P211-Tg001-Ab	Anti-IBV - Variant 02 Spike Glycoprotein S2 Subunit antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P212-Tg001-Ag	IBV - Variant 4/91 Spike Glycoprotein S3 Subunit antigen	protein
GMP-VT-P212-Tg001-Ab	Anti-IBV - Variant 4/91 Spike Glycoprotein S3 Subunit antibody	antibody(Fc/mFc)
GMP-VT-P213-Tg001-Ag	IBV - Variant D1466 Nucleoprotein antigen	protein
GMP-VT-P213-Tg001-Ab	Anti-IBV - Variant D1466 Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P214-Tg001-Ag	AIV-A H7N9 subtype Nucleoprotein antigen	protein
GMP-VT-P214-Tg001-Ab	Anti-AIV-A H7N9 subtype Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P214-Tg002-Ag	AIV-A H7N9 subtype Hemagglutinin antigen	protein
GMP-VT-P214-Tg002-Ab	Anti-AIV-A H7N9 subtype Hemagglutinin antibody	antibody(Fc/mFc)
GMP-VT-P215-Tg001-Ag	AIV-A H9N2 subtype Hemagglutinin antigen	protein
GMP-VT-P215-Tg001-Ab	Anti-AIV-A H9N2 subtype Hemagglutinin antibody	antibody(Fc/mFc)
GMP-VT-P215-Tg002-Ag	AIV-A H9N2 subtype H9 Hemagglutinin antigen	protein
GMP-VT-P215-Tg002-Ab	Anti-AIV-A H9N2 subtype H9 Hemagglutinin antibody	antibody(Fc/mFc)
GMP-VT-P215-Tg003-Ag	AIV-A H9N2 subtype Nucleoprotein antigen	protein
GMP-VT-P215-Tg003-Ab	Anti-AIV-A H9N2 subtype Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P216-Tg001-Ag	AIV-A N3 subtype Neuraminidase N3 antigen	protein
GMP-VT-P216-Tg001-Ab	Anti-AIV-A N3 subtype Neuraminidase N3 antibody	antibody(Fc/mFc)
GMP-VT-P217-Tg001-Ag	AIV-A N5 subtype Neuraminidase N5 antigen	protein
GMP-VT-P217-Tg001-Ab	Anti-AIV-A N5 subtype Neuraminidase N5 antibody	antibody(Fc/mFc)
GMP-VT-P218-Tg001-Ag	AIV-A N8 subtype Neuraminidase N8 antigen	protein
GMP-VT-P218-Tg001-Ab	Anti-AIV-A N8 subtype Neuraminidase N8 antibody	antibody(Fc/mFc)
GMP-VT-P219-Tg001-Ag	Salmonella enterica Pullorum Invasion Plasmid Antigen J antigen	protein
GMP-VT-P219-Tg001-Ab	Anti-Salmonella enterica Pullorum Invasion Plasmid Antigen J antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P220-Tg001-Ag	Riemerella anatipestifer Outer Membrane Protein H antigen	protein
GMP-VT-P220-Tg001-Ab	Anti-Riemerella anatipestifer Outer Membrane Protein H antibody	antibody(Fc/mFc)
GMP-VT-P221-Tg001-Ag	Salmonella Flagellin antigen	protein
GMP-VT-P221-Tg001-Ab	Anti-Salmonella Flagellin antibody	antibody(Fc/mFc)
GMP-VT-P221-Tg002-Ag	Salmonella Flagellin C antigen	protein
GMP-VT-P221-Tg002-Ab	Anti-Salmonella Flagellin C antibody	antibody(Fc/mFc)
GMP-VT-P221-Tg003-Ag	Salmonella Flagellin B antigen	protein
GMP-VT-P221-Tg003-Ab	Anti-Salmonella Flagellin B antibody	antibody(Fc/mFc)
GMP-VT-P222-Tg001-Ag	Turkey Astrovirus Capsid Protein antigen	protein
GMP-VT-P222-Tg001-Ab	Anti-Turkey Astrovirus Capsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P223-Tg001-Ag	THEV Hexon Protein antigen	protein
GMP-VT-P223-Tg001-Ab	Anti-THEV Hexon Protein antibody	antibody(Fc/mFc)
GMP-VT-P224-Tg001-Ag	Turkey Rhinotracheitis Fusion Protein antigen	protein
GMP-VT-P224-Tg001-Ab	Anti-Turkey Rhinotracheitis Fusion Protein antibody	antibody(Fc/mFc)
GMP-VT-P225-Tg001-Ag	Pigeon circovirus Capsid Protein antigen	protein
GMP-VT-P225-Tg001-Ab	Anti-Pigeon circovirus Capsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P226-Tg001-Ag	BFDV Capsid Protein antigen	protein
GMP-VT-P226-Tg001-Ab	Anti-BFDV Capsid Protein antibody	antibody(Fc/mFc)

Fish

Cat.No	Products Name	Products Type
GMP-VT-P227-Tg001-Ag	KHV KHV ORF 149 Protein antigen	protein

Cat.No	Products Name	Products Type
GMP-VT-P227-Tg001-Ab	Anti-KHV KHV ORF 149 Protein antibody	antibody(Fc/mFc)
GMP-VT-P228-Tg001-Ag	Salmon isavirus Haemagglutinin antigen	protein
GMP-VT-P228-Tg001-Ab	Anti-Salmon isavirus Haemagglutinin antibody	antibody(Fc/mFc)
GMP-VT-P228-Tg002-Ag	Salmon isavirus Nucleoprotein antigen	protein
GMP-VT-P228-Tg002-Ab	Anti-Salmon isavirus Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P229-Tg001-Ag	IPNV Viral Protein 2 antigen	protein
GMP-VT-P229-Tg001-Ab	Anti-IPNV Viral Protein 2 antibody	antibody(Fc/mFc)
GMP-VT-P230-Tg001-Ag	SVCV Glycoprotein antigen	protein
GMP-VT-P230-Tg001-Ab	Anti-SVCV Glycoprotein antibody	antibody(Fc/mFc)
GMP-VT-P231-Tg001-Ag	VHSV Nucleocapsid Protein antigen	protein
GMP-VT-P231-Tg001-Ab	Anti-VHSV Nucleocapsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P232-Tg001-Ag	IHNV Nucleoprotein antigen	protein
GMP-VT-P232-Tg001-Ab	Anti-IHNV Nucleoprotein antibody	antibody(Fc/mFc)

Multiple animals

Cat.No	Products Name	Products Type
GMP-VT-P002-Tg004-Ag	Anaplasma phagocytophilum antigen	protein
GMP-VT-P002-Tg004-Ab	Anti-Anaplasma phagocytophilum antibody	antibody(Fc/mFc)
GMP-VT-P026-Tg001-Ag	Leptospira Lipoprotein L32 antigen	protein
GMP-VT-P026-Tg001-Ab	Anti-Leptospira Lipoprotein L32 antibody	antibody(Fc/mFc)
GMP-VT-P031-Tg002-Ag	Anaplasma platys antigen	protein
GMP-VT-P031-Tg002-Ab	Anti-Anaplasma platys antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P100-Tg001-Ag	Brucella abortus/Brucella melitensis Outer Membrane Protein 2b antigen	protein
GMP-VT-P100-Tg001-Ab	Anti-Brucella abortus/Brucella melitensis Outer Membrane Protein 2b antibody	antibody(Fc/mFc)
GMP-VT-P100-Tg002-Ag	Brucella abortus/Brucella melitensis P17 Protein antigen	protein
GMP-VT-P100-Tg002-Ab	Anti-Brucella abortus/Brucella melitensis P17 Protein antibody	antibody(Fc/mFc)
GMP-VT-P100-Tg003-Ag	Brucella abortus/Brucella melitensis BP26 Protein antigen	protein
GMP-VT-P100-Tg003-Ab	Anti-Brucella abortus/Brucella melitensis BP26 Protein antibody	antibody(Fc/mFc)
GMP-VT-P100-Tg004-Ag	Brucella abortus/Brucella melitensis VirB5 Protein antigen	protein
GMP-VT-P100-Tg004-Ab	Anti-Brucella abortus/Brucella melitensis VirB5 Protein antibody	antibody(Fc/mFc)
GMP-VT-P100-Tg005-Ag	Brucella abortus/Brucella melitensis Outer Membrane Protein 28 antigen	protein
GMP-VT-P100-Tg005-Ab	Anti-Brucella abortus/Brucella melitensis Outer Membrane Protein 28 antibody	antibody(Fc/mFc)
GMP-VT-P100-Tg006-Ag	Brucella abortus/Brucella melitensis VirB12 Protein antigen	protein
GMP-VT-P100-Tg006-Ab	Anti-Brucella abortus/Brucella melitensis VirB12 Protein antibody	antibody(Fc/mFc)
GMP-VT-P101-Tg004-Ag	Pasteurella multocida antigen	protein
GMP-VT-P101-Tg004-Ab	Anti-Pasteurella multocida antibody	antibody(Fc/mFc)
GMP-VT-P143-Tg002-Ag	Babesia caballi P50 antigen	protein
GMP-VT-P143-Tg002-Ab	Anti-Babesia caballi P50 antibody	antibody(Fc/mFc)
GMP-VT-P174-Tg002-Ag	CAstV Capsid Protein antigen	protein
GMP-VT-P174-Tg002-Ab	Anti-CAstV Capsid Protein antibody	antibody(Fc/mFc)
GMP-VT-P233-Tg001-Ag	Staphylococcus saprophyticus antigen	protein
GMP-VT-P233-Tg001-Ab	Anti-Staphylococcus saprophyticus antibody	antibody(Fc/mFc)
GMP-VT-P234-Tg001-Ag	Trichinella Cystatin-Like Protein antigen	protein
GMP-VT-P234-Tg001-Ab	Anti-Trichinella Cystatin-Like Protein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P235-Tg001-Ag	Tritrichomonas suis Lipophosphoglycan antigen	protein
GMP-VT-P235-Tg001-Ab	Anti-Tritrichomonas suis Lipophosphoglycan antibody	antibody(Fc/mFc)
GMP-VT-P235-Tg002-Ag	Tritrichomonas suis TF1.17 Adhesin antigen	protein
GMP-VT-P235-Tg002-Ab	Anti-Tritrichomonas suis TF1.17 Adhesin antibody	antibody(Fc/mFc)
GMP-VT-P282-Tg001-Ag	E.coli Outer Membrane Protease T antigen	protein
GMP-VT-P282-Tg001-Ab	Anti-E.coli Outer Membrane Protease T antibody	antibody(Fc/mFc)
GMP-VT-P282-Tg002-Ag	E.coli antigen	protein
GMP-VT-P282-Tg002-Ab	Anti-E.coli antibody	antibody(Fc/mFc)
GMP-VT-P236-Tg001-Ag	Trichophyton mentagrophytes Keratinase antigen	protein
GMP-VT-P236-Tg001-Ab	Anti-Trichophyton mentagrophytes Keratinase antibody	antibody(Fc/mFc)
GMP-VT-P236-Tg002-Ag	Trichophyton mentagrophytes Secreted Protease Sub6 antigen	protein
GMP-VT-P236-Tg002-Ab	Anti-Trichophyton mentagrophytes Secreted Protease Sub6 antibody	antibody(Fc/mFc)
GMP-VT-P237-Tg001-Ag	Giardia Cyst Wall Protein 1 antigen	protein
GMP-VT-P237-Tg001-Ab	Anti-Giardia Cyst Wall Protein 1 antibody	antibody(Fc/mFc)
GMP-VT-P238-Tg001-Ag	Salmonella spp antigen	protein
GMP-VT-P238-Tg001-Ab	Anti-Salmonella spp antibody	antibody(Fc/mFc)
GMP-VT-P239-Tg001-Ag	Mycoplasma antigen	protein
GMP-VT-P239-Tg001-Ab	Anti-Mycoplasma antibody	antibody(Fc/mFc)
GMP-VT-P240-Tg001-Ag	Proteus antigen	protein
GMP-VT-P240-Tg001-Ab	Anti-Proteus antibody	antibody(Fc/mFc)
GMP-VT-P241-Tg001-Ag	Pseudomonas Elastase antigen	protein
GMP-VT-P241-Tg001-Ab	Anti-Pseudomonas Elastase antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P242-Tg001-Ag	Dermatophytes antigen	protein
GMP-VT-P242-Tg001-Ab	Anti-Dermatophytes antibody	antibody(Fc/mFc)
GMP-VT-P243-Tg001-Ag	Klebsiella pneumoniae antigen	protein
GMP-VT-P243-Tg001-Ab	Anti-Klebsiella pneumoniae antibody	antibody(Fc/mFc)
GMP-VT-P244-Tg001-Ag	WNV Envelope Protein antigen	protein
GMP-VT-P244-Tg001-Ab	Anti-WNV Envelope Protein antibody	antibody(Fc/mFc)
GMP-VT-P245-Tg001-Ag	campylobacter spp antigen	protein
GMP-VT-P245-Tg001-Ab	Anti-campylobacter spp antibody	antibody(Fc/mFc)
GMP-VT-P246-Tg001-Ag	HEV antigen	protein
GMP-VT-P246-Tg001-Ab	Anti-HEV antibody	antibody(Fc/mFc)
GMP-VT-P247-Tg001-Ag	Pestivirus antigen	protein
GMP-VT-P247-Tg001-Ab	Anti-Pestivirus antibody	antibody(Fc/mFc)
GMP-VT-P248-Tg001-Ag	Brucella antigen	protein
GMP-VT-P248-Tg001-Ab	Anti-Brucella antibody	antibody(Fc/mFc)
GMP-VT-P249-Tg001-Ag	Mycobacterium avium SahH antigen	protein
GMP-VT-P249-Tg001-Ab	Anti-Mycobacterium avium SahH antibody	antibody(Fc/mFc)
GMP-VT-P250-Tg001-Ag	Neospora caninum P38 antigen	protein
GMP-VT-P250-Tg001-Ab	Anti-Neospora caninum P38 antibody	antibody(Fc/mFc)
GMP-VT-P251-Tg001-Ag	Microsporium gypseum antigen	protein
GMP-VT-P251-Tg001-Ab	Anti-Microsporium gypseum antibody	antibody(Fc/mFc)
GMP-VT-P252-Tg001-Ag	Aspergillus brasiliensis antigen	protein
GMP-VT-P252-Tg001-Ab	Anti-Aspergillus brasiliensis antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P253-Tg001-Ag	Burkholderia mallei Hcp1 antigen	protein
GMP-VT-P253-Tg001-Ab	Anti-Burkholderia mallei Hcp1 antibody	antibody(Fc/mFc)
GMP-VT-P254-Tg001-Ag	Anaplasma phagocytophilum Major Surface Protein 5 antigen	protein
GMP-VT-P254-Tg001-Ab	Anti-Anaplasma phagocytophilum Major Surface Protein 5 antibody	antibody(Fc/mFc)
GMP-VT-P254-Tg002-Ag	Anaplasma phagocytophilum P44 antigen	protein
GMP-VT-P254-Tg002-Ab	Anti-Anaplasma phagocytophilum P44 antibody	antibody(Fc/mFc)
GMP-VT-P255-Tg001-Ag	Citrobacter freudii antigen	protein
GMP-VT-P255-Tg001-Ab	Anti-Citrobacter freudii antibody	antibody(Fc/mFc)
GMP-VT-P256-Tg001-Ag	RVFV Nucleoprotein antigen	protein
GMP-VT-P256-Tg001-Ab	Anti-RVFV Nucleoprotein antibody	antibody(Fc/mFc)
GMP-VT-P257-Tg001-Ag	FMDV Capsid Protein VP1 antigen	protein
GMP-VT-P257-Tg001-Ab	Anti-FMDV Capsid Protein VP1 antibody	antibody(Fc/mFc)
GMP-VT-P257-Tg002-Ag	FMDV Non-Structural Proteins 3ABC antigen	protein
GMP-VT-P257-Tg002-Ab	Anti-FMDV Non-Structural Proteins 3ABC antibody	antibody(Fc/mFc)
GMP-VT-P258-Tg001-Ag	B.b.s.l VlsE1 antigen	protein
GMP-VT-P258-Tg001-Ab	Anti-B.b.s.l VlsE1 antibody	antibody(Fc/mFc)
GMP-VT-P258-Tg002-Ag	B.b.s.l Outer Surface Protein C antigen	protein
GMP-VT-P258-Tg002-Ab	Anti-B.b.s.l Outer Surface Protein C antibody	antibody(Fc/mFc)
GMP-VT-P259-Tg001-Ag	Cryptosporidium parvum P23 Protein antigen	protein
GMP-VT-P259-Tg001-Ab	Anti-Cryptosporidium parvum P23 Protein antibody	antibody(Fc/mFc)
GMP-VT-P260-Tg001-Ag	Clostridium perfringens Epsilon Toxin antigen	protein
GMP-VT-P260-Tg001-Ab	Anti-Clostridium perfringens Epsilon Toxin antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P260-Tg002-Ag	Clostridium perfringens Beta2 Toxin antigen	protein
GMP-VT-P260-Tg002-Ab	Anti-Clostridium perfringens Beta2 Toxin antibody	antibody(Fc/mFc)
GMP-VT-P261-Tg001-Ag	Babesia canis Merozoite Surface Antigen 1 antigen	protein
GMP-VT-P261-Tg001-Ab	Anti-Babesia canis Merozoite Surface Antigen 1 antibody	antibody(Fc/mFc)
GMP-VT-P262-Tg001-Ag	Babesia gibsoni P50 antigen	protein
GMP-VT-P262-Tg001-Ab	Anti-Babesia gibsoni P50 antibody	antibody(Fc/mFc)
GMP-VT-P263-Tg001-Ag	Clostridium difficile Toxins A antigen	protein
GMP-VT-P263-Tg001-Ab	Anti-Clostridium difficile Toxins A antibody	antibody(Fc/mFc)
GMP-VT-P264-Tg001-Ag	Giardia lamblia Alpha-1 Giardin antigen	protein
GMP-VT-P264-Tg001-Ab	Anti-Giardia lamblia Alpha-1 Giardin antibody	antibody(Fc/mFc)
GMP-VT-P265-Tg001-Ag	Tritrichomonas gallinae Adhesion Protein antigen	protein
GMP-VT-P265-Tg001-Ab	Anti-Tritrichomonas gallinae Adhesion Protein antibody	antibody(Fc/mFc)
GMP-VT-P266-Tg001-Ag	C.botulinum type C and D Hemagglutinin C and D antigen	protein
GMP-VT-P266-Tg001-Ab	Anti-C.botulinum type C and D Hemagglutinin C and D antibody	antibody(Fc/mFc)
GMP-VT-P267-Tg001-Ag	Borrelia burgdorferi sensu lato Flagellin antigen	protein
GMP-VT-P267-Tg001-Ab	Anti-Borrelia burgdorferi sensu lato Flagellin antibody	antibody(Fc/mFc)
GMP-VT-P267-Tg002-Ag	Borrelia burgdorferi sensu lato P100 antigen	protein
GMP-VT-P267-Tg002-Ab	Anti-Borrelia burgdorferi sensu lato P100 antibody	antibody(Fc/mFc)
GMP-VT-P268-Tg001-Ag	Candida albicans antigen	protein
GMP-VT-P268-Tg001-Ab	Anti-Candida albicans antibody	antibody(Fc/mFc)
GMP-VT-P269-Tg001-Ag	Candida glabrata LOC1 Protein antigen	protein
GMP-VT-P269-Tg001-Ab	Anti-Candida glabrata LOC1 Protein antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P269-Tg002-Ag	Candida glabrata Acetyl-CoA Carboxylase antigen	protein
GMP-VT-P269-Tg002-Ab	Anti-Candida glabrata Acetyl-CoA Carboxylase antibody	antibody(Fc/mFc)
GMP-VT-P270-Tg001-Ag	Candida krusei Cytochrome C1 antigen	protein
GMP-VT-P270-Tg001-Ab	Anti-Candida krusei Cytochrome C1 antibody	antibody(Fc/mFc)
GMP-VT-P271-Tg001-Ag	Candida tropicalis Lanosterol 14-alpha Demethylase antigen	protein
GMP-VT-P271-Tg001-Ab	Anti-Candida tropicalis Lanosterol 14-alpha Demethylase antibody	antibody(Fc/mFc)
GMP-VT-P272-Tg001-Ag	E.faecalis Gelatinase antigen	protein
GMP-VT-P272-Tg001-Ab	Anti-E.faecalis Gelatinase antibody	antibody(Fc/mFc)
GMP-VT-P273-Tg001-Ag	EHDV Viral Protein 5 antigen	protein
GMP-VT-P273-Tg001-Ab	Anti-EHDV Viral Protein 5 antibody	antibody(Fc/mFc)
GMP-VT-P273-Tg002-Ag	EHDV Viral Protein 7 antigen	protein
GMP-VT-P273-Tg002-Ab	Anti-EHDV Viral Protein 7 antibody	antibody(Fc/mFc)
GMP-VT-P290-Tg001-Ag	E.coli O157:H7 Outer Membrane Protein C antigen	protein
GMP-VT-P290-Tg001-Ab	Anti-E.coli O157:H7 Outer Membrane Protein C antibody	antibody(Fc/mFc)
GMP-VT-P290-Tg002-Ag	E.coli O157:H7 Intimin Gamma 1 antigen	protein
GMP-VT-P290-Tg002-Ab	Anti-E.coli O157:H7 Intimin Gamma 1 antibody	antibody(Fc/mFc)
GMP-VT-P274-Tg001-Ag	S.cerevisiae Ubiquitin-Like-Specific Protease 1 antigen	protein
GMP-VT-P274-Tg001-Ab	Anti-S.cerevisiae Ubiquitin-Like-Specific Protease 1 antibody	antibody(Fc/mFc)
GMP-VT-P275-Tg001-Ag	S.enterica subsp. enterica Flagellin antigen	protein
GMP-VT-P275-Tg001-Ab	Anti-S.enterica subsp. enterica Flagellin antibody	antibody(Fc/mFc)
GMP-VT-P276-Tg001-Ag	S.enterica subsp. Typhimurium Flagellin antigen	protein
GMP-VT-P276-Tg001-Ab	Anti-S.enterica subsp. Typhimurium Flagellin antibody	antibody(Fc/mFc)

Cat.No	Products Name	Products Type
GMP-VT-P276-Tg002-Ag	S.enterica subsp. Typhimurium Outer Membrane Protein D antigen	protein
GMP-VT-P276-Tg002-Ab	Anti-S.enterica subsp. Typhimurium Outer Membrane Protein D antibody	antibody(Fc/mFc)
GMP-VT-P277-Tg001-Ag	S.aureus Fibronectin Binding Protein A antigen	protein
GMP-VT-P277-Tg001-Ab	Anti-S.aureus Fibronectin Binding Protein A antibody	antibody(Fc/mFc)
GMP-VT-P277-Tg002-Ag	S.aureus Secretory Antigen SsaA2 antigen	protein
GMP-VT-P277-Tg002-Ab	Anti-S.aureus Secretory Antigen SsaA2 antibody	antibody(Fc/mFc)
GMP-VT-P278-Tg001-Ag	T.gondii Surface Antigen 1 (SAG1) antigen	protein
GMP-VT-P278-Tg001-Ab	Anti-T.gondii Surface Antigen 1 (SAG1) antibody	antibody(Fc/mFc)
GMP-VT-P279-Tg001-Ag	Trichophyton rubrum Secreted Protease Sub6 antigen	protein
GMP-VT-P279-Tg001-Ab	Anti-Trichophyton rubrum Secreted Protease Sub6 antibody	antibody(Fc/mFc)
GMP-VT-P279-Tg002-Ag	Trichophyton rubrum Allergen Tri r 2 antigen	protein
GMP-VT-P279-Tg002-Ab	Anti-Trichophyton rubrum Allergen Tri r 2 antibody	antibody(Fc/mFc)
GMP-VT-P280-Tg001-Ag	Trichophyton tonsurans Allergen Tri t 4 antigen	protein
GMP-VT-P280-Tg001-Ab	Anti-Trichophyton tonsurans Allergen Tri t 4 antibody	antibody(Fc/mFc)
GMP-SMT-187-Ab	Anti-Growth hormone (hGH) antibody	antibody(Fc/mFc)
GMP-SMT-187-Ag	Growth hormone (hGH) antigen	chemicals-BSA/OVA
GMP-SMT-216-Ab	Anti-Vitamin B2 (Vb2) antibody	antibody(Fc/mFc)
GMP-SMT-216-Ag	Vitamin B2 (Vb2) antigen	chemicals-BSA/OVA
GMP-SMT-217-Ab	Anti-Vitamin B12 (Vb12) antibody	antibody(Fc/mFc)
GMP-SMT-217-Ag	Vitamin B12 (Vb12) antigen	chemicals-BSA/OVA
GMP-SMT-219-Ab	Anti-Biotin (BT) antibody	antibody(Fc/mFc)
GMP-SMT-219-Ag	Biotin (BT) antigen	chemicals-BSA/OVA

Cat.No	Products Name	Products Type
GMP-SMT-220-Ab	Anti-Lactoferrin (LF) antibody	antibody(Fc/mFc)
GMP-SMT-220-Ag	Lactoferrin (LF) antigen	chemicals-BSA/OVA
GMP-SMT-221-Ab	Anti-FITC antibody	antibody(Fc/mFc)
GMP-SMT-221-Ag	FITC antigen	chemicals-BSA/OVA
GMP-SMT-222-Ab	Anti-Histamine antibody	antibody(Fc/mFc)
GMP-SMT-222-Ag	Histamine antigen	chemicals-BSA/OVA
GMP-SMT-25-OH-(VD3+VD2)-Ab	Anti-25 hydroxyvitamin D (25-OH- (VD3+VD2)) antibody	antibody(Fc/mFc)
GMP-SMT-25-OH-(VD3+VD2)-Hapten	25 hydroxyvitamin D (25-OH- (VD3+VD2)) antigen	chemicals-BSA/OVA
GMP-SMT-25-OH-VD-3-Ab	Anti-25-hydroxy (OH) Vitamin 3 (25-OH-VD-3) antibody	antibody(Fc/mFc)
GMP-SMT-25-OH-VD-3-Ag	25-hydroxy (OH) Vitamin 3 (25-OH-VD-3) antigen	chemicals-BSA/OVA
GMP-SMT-25-OH-VD-Ab	Anti-25-hydroxy (OH) Vitamin (25-OH-VD) antibody	antibody(Fc/mFc)
GMP-SMT-25-OH-VD-Ag	25-hydroxy (OH) Vitamin (25-OH-VD) antigen	chemicals-BSA/OVA
GMP-SMT-Bisphenol-A-Ab	Anti-Bisphenol A antibody	antibody(Fc/mFc)
GMP-SMT-Bisphenol-A-Ag	Bisphenol A antigen	chemicals-BSA/OVA
GMP-SMT-Calcitonin-Ab	Anti-Calcitonin antibody	antibody(Fc/mFc)
GMP-SMT-Calcitonin-Ag	Calcitonin antigen	chemicals-BSA/OVA
GMP-SMT-CCP-Ab	Anti-cyclic citrullinated peptide (CCP) antibody	antibody(Fc/mFc)
GMP-SMT-CCP-Ag	cyclic citrullinated peptide (CCP) antigen	chemicals-BSA/OVA
GMP-SMT-CG-Ab	Anti-Cholyglycine (CG) antibody	antibody(Fc/mFc)
GMP-SMT-CG-Ag	Cholyglycine (CG) antigen	chemicals-BSA/OVA
GMP-SMT-cholanoic-acid-Ab	Anti-deoxycholic acid (cholanoic acid) antibody	antibody(Fc/mFc)
GMP-SMT-cholanoic-acid-Ag	deoxycholic acid (cholanoic acid) antigen	chemicals-BSA/OVA

Cat.No	Products Name	Products Type
GMP-SMT-Cortisol-Ab	Anti-Cortisol antibody	antibody(Fc/mFc)
GMP-SMT-Cortisol-Ag	Cortisol antigen	chemicals-BSA/OVA
GMP-SMT-cr-Ab	Anti-creatinine (cr) antibody	antibody(Fc/mFc)
GMP-SMT-cr-Ag	creatinine (cr) antigen	chemicals-BSA/OVA
GMP-SMT-DES-HES-Ab	Anti-Diethylstilbestrol/Hexestrol (DES/HES) antibody	antibody(Fc/mFc)
GMP-SMT-DES-HES-Ag	Diethylstilbestrol/Hexestrol (DES/HES) antigen	chemicals-BSA/OVA
GMP-SMT-DEX-Ab	Anti-Dexamethasone (DEX) antibody	antibody(Fc/mFc)
GMP-SMT-DEX-Ag	Dexamethasone (DEX) antigen	chemicals-BSA/OVA
GMP-SMT-DHEA-Ab	Anti-dehydroepiandrosterone (DHEA) antibody	antibody(Fc/mFc)
GMP-SMT-DHEA-Ag	dehydroepiandrosterone (DHEA) antigen	chemicals-BSA/OVA
GMP-SMT-Digoxin-Ab	Anti-digoxigenin (Digoxin) antibody	antibody(Fc/mFc)
GMP-SMT-Digoxin-Ag	digoxigenin (Digoxin) antigen	chemicals-BSA/OVA
GMP-SMT-E1-Ab	Anti-Estrone (E1) antibody	antibody(Fc/mFc)
GMP-SMT-E1-Ag	Estrone (E1) antigen	chemicals-BSA/OVA
GMP-SMT-E2-E3-Ab	Anti-Estradiol (E2/E3) antibody	antibody(Fc/mFc)
GMP-SMT-E2-E3-Ag	Estradiol (E2/E3) antigen	chemicals-BSA/OVA
GMP-SMT-EPO-Ab	Anti-Erythropoetin (EPO) antibody	antibody(Fc/mFc)
GMP-SMT-EPO-Ag	Erythropoetin (EPO) antigen	chemicals-BSA/OVA
GMP-SMT-FA-Ab	Anti-folic acid (FA) antibody	antibody(Fc/mFc)
GMP-SMT-FA-Ag	folic acid (FA) antigen	chemicals-BSA/OVA
GMP-SMT-Flugestone-Acetate-Ab	Anti-Flugestone Acetate antibody	antibody(Fc/mFc)
GMP-SMT-Flugestone-Acetate-Ag	Flugestone Acetate antigen	chemicals-BSA/OVA

Cat.No	Products Name	Products Type
GMP-SMT-Forchlorfenuron-Ab	Anti-Forchlorfenuron antibody	antibody(Fc/mFc)
GMP-SMT-Forchlorfenuron-Ag	Forchlorfenuron antigen	chemicals-BSA/OVA
GMP-SMT-G17-Ab	Anti-Gastrin 17 (G17) antibody	antibody(Fc/mFc)
GMP-SMT-G17-Ag	Gastrin 17 (G17) antigen	chemicals-BSA/OVA
GMP-SMT-HA-Ab	Anti-hyaluronic acid (HA) antibody	antibody(Fc/mFc)
GMP-SMT-HA-Ag	hyaluronic acid (HA) antigen	chemicals-BSA/OVA
GMP-SMT-HCY-Ab	Anti-Homocysteine (HCY) antibody	antibody(Fc/mFc)
GMP-SMT-HCY-Ag	Homocysteine (HCY) antigen	chemicals-BSA/OVA
GMP-SMT-HFABP-Ab	Anti-Human Fatty Acid Binding Protein (HFABP) antibody	antibody(Fc/mFc)
GMP-SMT-HFABP-Ag	Human Fatty Acid Binding Protein (HFABP) antigen	chemicals-BSA/OVA
GMP-SMT-IGFBP-1-Ab	Anti-Insulin-like growth factor binding protein 1 (IGFBP-1) antibody	antibody(Fc/mFc)
GMP-SMT-IGFBP-1-Ag	Insulin-like growth factor binding protein 1 (IGFBP-1) antigen	chemicals-BSA/OVA
GMP-SMT-LCA-Ab	Anti-Lithocholic Acid (LCA) antibody	antibody(Fc/mFc)
GMP-SMT-LCA-Ag	Lithocholic Acid (LCA) antigen	chemicals-BSA/OVA
GMP-SMT-LH-Ab	Anti-Luteinizing hormone (LH) antibody	antibody(Fc/mFc)
GMP-SMT-LH-Ag	Luteinizing hormone (LH) antigen	chemicals-BSA/OVA
GMP-SMT-LP-PLA2-Ab	Anti-Phospholipase A2 (LP-PLA2) antibody	antibody(Fc/mFc)
GMP-SMT-LP-PLA2-Ag	Phospholipase A2 (LP-PLA2) antigen	chemicals-BSA/OVA
GMP-SMT-MPA-Ab	Anti-mycophenolic acid (MPA) antibody	antibody(Fc/mFc)
GMP-SMT-MPA-Ag	mycophenolic acid (MPA) antigen	chemicals-BSA/OVA
GMP-SMT-Norethisterone-Ab	Anti-Norethisterone antibody	antibody(Fc/mFc)
GMP-SMT-Norethisterone-Ag	Norethisterone antigen	chemicals-BSA/OVA

Cat.No	Products Name	Products Type
GMP-SMT-Paracetamol-Ab	Anti-acetaminophen (Paracetamol) antibody	antibody(Fc/mFc)
GMP-SMT-Paracetamol-Ag	acetaminophen (Paracetamol) antigen	chemicals-BSA/OVA
GMP-SMT-PG-Ab	Anti-Prostaglandin (PG) antibody	antibody(Fc/mFc)
GMP-SMT-PG-Ag	Prostaglandin (PG) antigen	chemicals-BSA/OVA
GMP-SMT-Progesterone-Ab	Anti-Progesterone antibody	antibody(Fc/mFc)
GMP-SMT-Progesterone-Ag	Progesterone antigen	chemicals-BSA/OVA
GMP-SMT-Prolactin-Ab	Anti-Prolactin antibody	antibody(Fc/mFc)
GMP-SMT-Prolactin-Ag	Prolactin antigen	chemicals-BSA/OVA
GMP-SMT-RAA-Aldo-Ab	Anti-aldosterone (RAA/Aldo) antibody	antibody(Fc/mFc)
GMP-SMT-RAA-Aldo-Ag	aldosterone (RAA/Aldo) antigen	chemicals-BSA/OVA
GMP-SMT-RBP4-Ab	Anti-Retinol binding protein 4 (RBP4) antibody	antibody(Fc/mFc)
GMP-SMT-RBP4-Ag	Retinol binding protein 4 (RBP4) antigen	chemicals-BSA/OVA
GMP-SMT-rT3-Ab	Anti-Reverse triiodothyronine (rT3) antibody	antibody(Fc/mFc)
GMP-SMT-rT3-Ag	Reverse triiodothyronine (rT3) antigen	chemicals-BSA/OVA
GMP-SMT-SA-Ab	Anti-Streptavidin (SA) antibody	antibody(Fc/mFc)
GMP-SMT-SA-Ag	Streptavidin (SA) antigen	chemicals-BSA/OVA
GMP-SMT-SAH-Ab	Anti-S-Adenosylhomocysteine (SAH) antibody	antibody(Fc/mFc)
GMP-SMT-SAH-Ag	S-Adenosylhomocysteine (SAH) antigen	chemicals-BSA/OVA
GMP-SMT-TBG-Ab	Anti-Thyroid Binding Globulin (TBG) antibody	antibody(Fc/mFc)
GMP-SMT-TBG-Ag	Thyroid Binding Globulin (TBG) antigen	chemicals-BSA/OVA
GMP-SMT-Testosterone-Ab	Anti-Testosterone antibody	antibody(Fc/mFc)
GMP-SMT-Testosterone-Ag	Testosterone antigen	chemicals-BSA/OVA

Cat.No	Products Name	Products Type
GMP-SMT-TFF1-Ab	Anti-trefoil factor 1 (TFF1) antibody	antibody(Fc/mFc)
GMP-SMT-TFF1-Ag	trefoil factor 1 (TFF1) antigen	chemicals-BSA/OVA
GMP-SMT-TFF2-Ab	Anti-trefoil factor 2 (TFF2) antibody	antibody(Fc/mFc)
GMP-SMT-TFF2-Ag	trefoil factor 2 (TFF2) antigen	chemicals-BSA/OVA
GMP-SMT-TFF3-Ab	Anti-trefoil factor 3 (TFF3) antibody	antibody(Fc/mFc)
GMP-SMT-TFF3-Ag	trefoil factor 3 (TFF3) antigen	chemicals-BSA/OVA
GMP-SMT-Thyroxin-T3-Ab	Anti-Thyroxin T3 antibody	antibody(Fc/mFc)
GMP-SMT-Thyroxin-T3-Ag	Thyroxin T3 antigen	chemicals-BSA/OVA
GMP-SMT-Thyroxin-T4-Ab	Anti-Thyroxin T4 antibody	antibody(Fc/mFc)
GMP-SMT-Thyroxin-T4-Ag	Thyroxin T4 antigen	chemicals-BSA/OVA
GMP-SMT-VPA-Ab	Anti-Sodium valproate (VPA) antibody	antibody(Fc/mFc)
GMP-SMT-VPA-Ag	Sodium valproate (VPA) antigen	chemicals-BSA/OVA

Pseudotype | Human Host

Pseudotype Name	Gene Name	Host	Classification
Canine adenovirus 1	HARA E3	Canine	DNA
Felis catus (domestic cat)	GAPDH	Feline	DNA
African swine fever virus	ASFV	Swine	DNA
PRRSV and respiratory syndrome virus	NSP2	Swine	DNA
Canine parainfluenza virus	N	Canine	RNA
Canine morbillivirus	NP	Canine	RNA
Feline calicivirus	F65	Feline	RNA
Feline leukemia virus isolate	TP2R	Feline	RNA
Feline immunodeficiency virus	CHN17	Feline	RNA
Hog cholera virus	CAP	Swine	RNA
Bovine viral diarrhea virus	BVDV	Bovine	RNA
Feline calicivirus isolate	HRB-SS	Feline	RNA

Detection of SAA protein in cats, dogs and horses

Detection of SAA in feline plasma samples

We used the recommended paired recombinant feline SAA antigens (Cat. No. [GMP-FEL-Serum-amyloid-A-Ag02] and [GMP-FEL-Serum-amyloid-A-Ag01]) and feline plasma samples for simultaneous dilution testing. Figure 1 shows a representative dilution curve of recombinant antigens in mixed feline plasma samples using the paired [GMP-FEL-Serum-amyloid-A-Ab01] - [GMP-FEL-Serum-amyloid-A-Ab02] . The dilution curves of the paired [GMP-FEL-Serum-amyloid-A-Ab05] - [GMP-FEL-Serum-amyloid-A-Ab06] for healthy cats and inflammatory cats are shown in Figure 2 , A and B, respectively.

For antibodies that yield parallel dilution curves when detecting feline plasma samples and recombinant antigen, recombinant feline SAA can be used as a calibrator to detect SAA in natural samples. For the pairing [GMP-FEL-Serum-amyloid-A-Ab10] - [GMP-FEL-Serum-amyloid-A-Ab11] , we recommend using the recombinant antigen containing Tag tag (Cat. No. [GMP-FEL-Serum-amyloid- A-Ag02]) as a calibrator. For other recommended pairings, both recombinant antigens can be used.




Figure 1. Calibration curve of paired [GMP-FEL-Serum-amyloid-A-Ab01] - [GMP-FEL-Serum-amyloid-A-Ab02] prototype reagent (the antigens are label-free recombinant cat SAA - Catalog number [GMP-FEL-Serum-amyloid-A-Ag01] and mixed cat plasma), The SAA concentration of mixed cat plasma was 133 µg/ml and was diluted 219.5 times before testing.




Figure 2. Dilution curve of paired [GMP-FEL-Serum-amyloid-A-Ab05] - [GMP-FEL-Serum-amyloid-A-Ab06] prototype reagents for detection of recombinant antigen and plasma samples. For healthy individuals (N1, N8), the SAA concentrations were 0.55µg/ml and 0.35µg/ml respectively. Before testing, the plasma sample was diluted 35 times. The initial concentration of recombinant feline SAA (catalog number [GMP-FEL-Serum-amyloid-A-Ag01]) is 72 ng/ml. (B) Postoperative cats (A26, A30). SAA concentrations were 178 µg/ml and 217 µg/ml respectively. Before testing, the plasma sample was diluted 200 times. The initial concentration of recombinant feline SAA (catalog number [GMP-FEL-Serum-amyloid-A-Ag01]) is 200 ng/ml.

We used prototype reagents to test 40 samples from healthy cats and 33 samples from inflammatory cats. Compared with normal cats, the SAA levels in inflammatory cats were significantly higher (Figure 3).




Figure 3. SAA concentrations in healthy cats and cats with postoperative inflammation. Samples from healthy individuals were pre-diluted 100-2500 times, and samples from inflammatory cats were pre-diluted 8000 times.

Recombinant canine SAA sandwich immunoassay system

The dilution curve of the recombinant canine antigen tested using the paired [GMP-CAN-Serum-amyloid-A-Ab01]-[GMP-CAN-Serum-amyloid-A-Ab02] is shown in Figure 4




Figure 4. Dilution curve of paired [GMP-CAN-Serum-amyloid-A-Ab01] - [GMP-CAN-Serum-amyloid-A-Ab02] test recombinant canine antigen. [GMP-CAN-Serum-amyloid-A-Ab01] is coated in a microplate and blocked at 37°C for 1 hour (blocking solution contains 1% casein and 0.05% Tween 20). Recombinant canine SAA (Cat. No. [GMP-CAN-Serum-amyloid-A-Ag01]) was serially diluted using the same diluent and incubated at 37 degrees for 1 hour. Europium-labeled [GMP-CAN-Serum-amyloid-A-Ab02] was used as the detection antibody .

Figure 5 shows the test results of dog serum SAA paired with [GMP-CAN-Serum-amyloid-A-Ab02] - [GMP-CAN-Serum-amyloid-A-Ab05] . The picture shows that SAA levels in diseased individuals are significantly higher than in normal individuals.




Figure 5 . Use the pairing [GMP-CAN-Serum-amyloid-A-Ab02]- [GMP-CAN-Serum-amyloid-A-Ab05] to compare the SAA detection results of sick individuals and healthy individuals in dogs . See Figure 4 for the reagent preparation process. The serum of healthy individuals is pre-diluted 50 times, and the serum of sick individuals is pre-diluted 1000 (horse)/2000 times (dog). The diluent is a buffer containing 1% casein and 0.05% Tween 20.

Recombinant equine SAA sandwich immunoassay system

The dilution curve of the recombinant horse antigen tested using the paired [GMP-EQU-Serum-amyloid-A-Ab01]-[GMP-EQU-Serum-amyloid-A-Ab02] is shown in Figure 6.




Figure 6. Dilution curve of paired [GMP-EQU-Serum-amyloid-A-Ab01] - [GMP-EQU-Serum-amyloid-A-Ab02] testing recombinant horse SAA antigen. [GMP-EQU-Serum-amyloid-A-Ab01] was coated in a microplate and blocked for 1 hour at 37°C (blocking buffer contains 1% casein and 0.05% Tween 20). Recombinant horse SAA [GMP-EQU-Serum-amyloid-A-Ag01] was serially diluted using the same dilution buffer and incubated at 37°C for 1 hour. Europium-labeled [GMP-EQU-Serum-amyloid-A-Ab02] was used as the detection antibody.

Figure 7 shows the test results of horse serum SAA paired [GMP-EQU-Serum-amyloid-A-Ab01] - [GMP-EQU-Serum-amyloid-A-Ab02] . The picture shows that SAA levels in diseased individuals are significantly higher than in normal individuals.




Figure 7 . Use the pairing [GMP-EQU-Serum-amyloid-A-Ab02] -[GMP-EQU-Serum-amyloid-A-Ab05] to compare the SAA detection results of diseased equine individuals and healthy individuals. See Figure 6 for the reagent preparation process. The serum of healthy individuals was pre-diluted 50 times, and the serum of sick individuals was pre-diluted 1000 times. The diluent was a buffer containing 1% casein and 0.05% Tween 20.

A single immunoassay system that simultaneously detects feline, canine and equine SAA

Pairing [GMP-FEL-Serum-amyloid-A-Ab10]-[GMP-FEL-Serum-amyloid-A-Ab11] is not only suitable for the development of cat SAA detection reagents, but also suitable for the development of canine and equine SAA detection reagents. Figure 6 shows the detection results of SAA in serum samples of cats (A), dogs (B) and horses (C) by this pairing. Among them, the SAA results of diseased individuals were significantly higher than those of normal individuals




Figure 6 . Test results of paired [GMP-FEL-Serum-amyloid-A-Ab10]-[GMP-FEL-Serum-amyloid-A-Ab11] for cat (A), dog (B) and horse (C) serum samples. SAA19 is a coating antibody and is blocked for 1 hour (the blocking solution contains 2.5% sodium caseinate). Canine serum is pre-diluted 500 times, cat and horse serum is pre-diluted 800 times, the diluent is Tris buffer containing 0.01% CHAPS, and the washing solution and diluent are the same soluti

